**Strisce per il test delle urine per autodiagnosi – 2 parametri (Glucosio e Chetone)**  
**2-parameter (Glucose and Ketone) Urine Test**  
**Bandelettes d'analyse d'urine pour auto-diagnostic - 2 paramètres (glucose et corps cétoniques)**  
**Urinteststreifen zur Selbstdiagnose - 2 Parameter (Glucose und Keton)**  
**Tiras para la prueba de orina para autodiagnóstico – 2 parámetros (Glucosa y Cetona)**  
**Tiras para teste das urinas para auto-diagnoses – 2 parâmetros (Glucose e Cetona)**  
**Λωρίδες για το τεστ ούρων Για αυτοδιάγνωση– 2 παράμετροι (Γλυκόζη και Κετόνη)**  
**شرائح لفحص البول 2 – يتاذل اصي خشنتلل ميق (نوت يكوو زوكولج)**

MANUALE D'USO  
OPERATOR'S MANUAL  
MANUEL D'UTILIZATION  
BEDIENUNGSANLEITUNG  
MANUAL DE USO  
MANUAL DE USO  
Εγχειρίδιο χρήσης  
دليل للإرشادات

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## BIBLIOGRAFIA

A. H. Free and H. M. Free "Urinalysis, critical discipline of clinical science" CRC Critical Reviews in Clinical Laboratory Sciences, 481-531, 1972.

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GIMA Spa  
Via Marconi, 1  
Fabbricante 20060 Gessate (MI) - Italia

## INTENDED USE

The 2-parameter (Glucose and Ketone) Urine Test Strips contains solid phase reagent areas affixed to a plastic stick. They are provided as a dry reagent. The 2-parameter (Glucose and Ketone) Urine Test Strips provide tests for the semi-quantitative determinations of glucose and ketone in human urine.

## IMPORTANT NOTICE

YOU MUST READ THE INSERT BEFORE RUN THE TEST. IF YOU HAVE ANY POSITIVE RESULTS OR NEGATIVE RESULTS BUT THE SYMPTOMS CONTINUE, YOU SHOULD VISIT YOUR DOCTOR.

## TEST PRINCIPLES

The test principles are based on various dyes and reagent reactions with components of the urine that lead to colored components, which can be visually detected and/or measured by the instrument.

## MATERIALS PROVIDED

1. The 2-parameter (Glucose and Ketone) urine test strips (Each pouch contains three tests).
2. Color label chart.
3. Instructions for use.

## PRECAUTIONS

1. For in Vitro diagnostic use.
2. Do not touch test areas of strip.
3. Do not use after expiration date.
4. Do not refrigerate or freeze.
5. Store all test strips in sealed a pouch at room temperature, out of direct sun light. Do not remove the desiccant from pouch.

## SPECIMEN COLLECTION

1. Urine should be collected in a urine collection cup.
2. If testing cannot be done within an hour after voiding, refrigerate the specimen immediately. Return to room temperature before testing.

## STORAGE

1. Store at room temperature, out of direct sun light.
2. Do not use after expiration date. Shelf life is 24 months.
3. Do not refrigerate or freeze.
4. Store all test strips in sealed a pouch. Do not remove the desiccant from pouch.

## TEST PROCEDURE

1. Bring specimens to room temperature before use.
2. Remove a 2-parameter (Glucose and Ketone) test strip from the pouch.
3. Immerse test areas of the strip completely in urine and remove immediately to avoid dissolving of reagents.
4. To remove excess urine, run the edge of the strip against the rim of the urine container. Hold the strip in a horizontal position to prevent possible mixing of chemicals from adjacent reagent areas. Excess urine may also be removed by gently blotting the lengthwise edge on absorbent paper.

5. Compare the test results carefully with the color chart in good light.
6. Note: The optimal reading time of each test parameter varies from 30 seconds up to 2 minutes. Changes in color that appear only on the edges of the test areas or after more than 2 minutes are of no clinical significance.
7. If you have color vision issue, you need people to help you during comparison with color chat and interpretation on results.

## RESULTS

The results are obtained by direct comparison of the test strip with the color chart at the specified times.

## REAGENT COMPOSITION

**Glucose:** 10.54% w/w glucose oxidase (aspergillus, 250 IU), 0.2% w/w peroxidase (horseradish, 2,500 IU), 5.0% w/w potassium iodide and 84.3% non-reactive ingredients.  
**Ketone:** 4.5% w/w sodium nitroprusside and 95.5% w/w buffer.

## LIMITATIONS

**Glucose:** Large amounts of ketone bodies (50 mg/dl or greater) may decrease color development.

**Ketone:** Color reactions that could be interpreted as „positive“ may be obtained with urine specimens containing medium or large amounts of phenylketones.

## PERFORMANCE CHARACTERISTICS

### Accuracy:

The accuracy of Glucose and Ketone; 94.77% and 92.17% between Lay User and Reference Lab

We also found few participants show more mistakes in test results that may tells that he/ she may have color vision problems. If you have color sensation issue and need people help you to read the color chat.

### Sensitivity:

Glucose	100 - 250 mg/dl
Ketone	5 - 10 mg/dl

### Specificity:

Glucose	98.98%
Ketone	99.13%

## NORMAL VALUE REFERENCE

Glucose	Negative
Ketone	Negative

## QUESTIONS & ANSWERS

**Q:** If I see color change show in color chat in Glucose or Ketone reagent area, do I have Diabetes?

**A:** Yes. This should be interpreted as a positive result.

**Q:** Can the result change after standing for a certain length of time?

**A:** In order to prevent any incorrect results, the test result should not be interpreted after 2 minute.

**Q:** If the test result is positive, what do I do?

**A:** If the test result is positive, you should see your doctor. This is home-use test only for screening test done by non-professional. The earlier you see your doctor the better for further tests.

**Q:** If the test result is negative, what do I do?

**A:** A negative result means that no glucose and ketone has been detected. If your symptoms persist, repeat the test. If the second test is still negative and you still have not recovered from your symptoms, you should see your doctor.

**Q:** Can I get a false or misleading result?

**A:** This should not happen if you follow the instructions correctly. However, if you get the absorbent pad wet before using the test, then you may not get a reliable result. Certain rare medical conditions can also cause a false result. If you repeatedly get unexpected results, you should consult your doctor.

## REFERENCE

A. H. Free and H. M. Free “ Urinalysis, critical discipline of clinical science” CRC Critical Reviews in Clinical Laboratory Sciences, 481-531, 1972.

## Graphical Symbols used

	Attention, <sup>see</sup> instructions for use		Tests per kit		Manufacturer
	For <i>in vitro</i> diagnostic use only		Use by		Do not reuse
	Store between 4-30°C		Lot Number		Catalog <b>24065</b>
	Read instructions for use		Keep dry		Keep away from sunlight

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