Strisce per il test delle urine - Per autodiagnosi - 2 parametri (Leucociti e Nitriti)
2-parameter (Leucocytes and Nitrite) Urine Test Strips
Bandelettes d'analyse d'urine pour auto-diagnostic 2 paramètres (leucocytes et nitrites)
Urinteststreifen zur Selbstdiagnose - 2 Urinparameter (Leukozyten und Nitrite)
Tiras para la prueba de orina Para autodiagnóstico - 2 parámetros (Leucocitos y Nitritos)
Tiras para teste das urinas - Para auto-diagnoses - 2 parâmetros (Leucócitos e Nitritos)
Λωριδες για το τεστ ούρων Για αυτοδιάγνωση 2 παράμετροι (Λευκοκύτταρα και Νιτρώδη) ούρων
شرائح لفحص البول للتشخيص الذاتي 2 قيم (كريات بيضاء ونتريت)

MANUALE D'USO
OPERATOR'S MANUAL
MANUEL D'UTILIZATION
BEDIENUNGSANLEITUNG
MANUAL DE USO
MANUAL DE USO
Εγχειρίδιο χρήσης
دليل للإرشادات

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 **GIMA Spa**
Via Marconi, 1
Fabbricante 20060 Gessate (MI) - Italia

INTENDED USE

The 2-parameter (Leucocytes and Nitrite) Urine Test Strips contains solid phase reagent areas affixed to a plastic stick. They are provided as a dry reagent. The 2-parameter (Leucocytes and Nitrite) Urine Test Strips provide tests for the semi-quantitative determinations of Leucocytes and nitrite in human urine to diagnosis of urinary track infection.

IMPORTANT NOTICE

YOU MUST READ THE INSERT BEFORE RUN THE TEST. IF YOU HAVE ANY POSITIVE RESULTS OR NEGATIVE RESULTS BUT THE SYMPTOMS CONTINIUE, YOU SHOULD VISIT YOUR DOCTOR.

TEST PRINCIPLES

The test principles are based on various dyes and reagent reactions with components of the urine that lead to colored components, which can be visually detected and/or measured by the instrument.

MATERIALS PROVIDED

1. 2-parameter (Leucocytes and Nitrite) urine test strips (Each pouch contains three tests).
2. Color label chart.
3. Instructions for use.

MATERIALS REQUIRED BUT NOT PROVIDED

1. Urine collection cup.
2. Clock or timer.

PRECAUTIONS

1. For in Vitro diagnostic use.
2. Do not touch test areas of strip.
3. Do not use after expiration date. Shelf life is 24 months.
4. Do not refrigerate or freeze.
5. Store all test strips in sealed a pouch at room temperature, out of direct sun light. Do not remove the desiccant from pouch.

SPECIMEN COLLECTION

1. Urine should be collected in a urine collection cup.
2. If testing cannot be done within an hour after voiding, refrigerate the specimen immediately. Return to room temperature before testing.

STORAGE

1. Store at room temperature, out of direct sun light.
2. Do not use after expiration date. Shelf life is 24 months.
3. Do not refrigerate or freeze.
4. Store all test strips in sealed a pouch. Do not remove the desiccant from pouch.

TEST PROCEDURE

1. Bring specimens to room temperature before use.
2. Remove a 2-parameter (Leucocytes and Nitrite) test strip from the pouch.

3. Inspect the strip. (Discoloration or darkening of reagent test areas may indicate deterioration. Do not use the strip.)
4. Immerse test areas of the strip completely in urine and remove immediately to avoid dissolving of reagents.
5. To remove excess urine, run the edge of the strip against the rim of the urine container. Hold the strip in a horizontal position to prevent possible mixing of chemicals from adjacent reagent areas. Excess urine may also be removed by gently blotting the lengthwise edge on absorbent paper.
6. Compare the test results carefully with the color chart in good light.
7. **Note:** The optimal reading time of each test parameter varies from 30 seconds up to 2 minutes. Changes in color that appear only on the edges of the test areas or after more than 2 minutes are of no clinical significance.
8. If you have color vision issue, you need people to help you during comparison with color chart and interpretation on results.

RESULTS

The results are obtained by direct comparison of the test strip with the color chart at the specified times.

REAGENT COMPOSITION

Nitrite: 2.0% w/w p-arsanilic acid, 2.2% w/w a-naphthylamine and 95.8% w/w buffer.

Leukocytes: 0.1% w/w ester, 0.6% w/w diazonium salt, 40% w/w buffer and 59.3% w/w nonreactive ingredients.

LIMITATIONS

Nitrite: Any degree of uniform pink color development should be considered positive, however, pink spots or pink edges should not be interpreted as a positive result.

EXPECTED VALUES

Nitrite: Any degree of pink color after 30 seconds indicates a positive test for nitrite suggesting a clinically significant infection by bacteria. A negative test does not necessarily rule out bacterial infection.

Leukocytes: negative

PERFORMANCE CHARACTERISTICS

Accuracy:

The accuracy of Leukocytes and Nitrite; 89% and 93% between Lay User and Reference Lab. We also found few participants show more mistakes in test results that may tells that he/she may have color vision problems. If you have color sensation issue and need people help you to read the color chart.

Sensitivity:

Nitrite 0.06- 0.1 mg/dl

Leukocytes 5 - 15 cells/h

Specificity:

Nitrite 99.17%

Leukocytes 98.33%

NORMAL VALUE REFERENCE

Leukocytes Negative
Nitrite Negative

QUESTIONS & ANSWERS

Q: If I see color change show in color of Leukocytes or Nitrite reagent area, do I have urinary track infection?

A: Yes. This should be interpreted as a positive result that you may have urinary track infection.

Q: Can the result change after standing for a certain length of time?

A: In order to prevent any incorrect results, the test result should not be interpreted after 2 minute.

Q: If the test result is positive, what do I do?

A: If the test result is positive, you should see your doctor. This is home-use test only for screening test done by non-professional. The earlier you see your doctor the better for further tests.

Q: If the test result is negative, what do I do?

A: A negative result means you are ok. If your symptoms persist, repeat the test. If the second test is still negative and you still have not recovered from your symptoms, you should see your doctor.

Q: Can I get a false or misleading result?

A: This should not happen if you follow the instructions correctly. However, if you get the absorbent pad wet before using the test, then you may not get a reliable result. Certain rare medical conditions can also cause a false result. If you repeatedly get unexpected results, you should consult your doctor.

REFERENCE

A. H. Free and H. M. Free "Urinalysis, critical discipline of clinical science" Med. Aust., I, R47, 1960.

Graphical Symbols used

	Attention, ^{see} instructions for use		Tests per kit		Manufacturer
	For <i>in vitro</i> diagnostic use only		Use by		Do not reuse
	Store between 4-30°C		Lot Number		Catalog 24066
	Read instructions for use		Keep dry		Keep away from sunlight

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