
# GIMA

PROFESSIONAL MEDICAL PRODUCTS

Gima S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) Italy  
Italia: tel. 199 400 401 - fax 199 400 403  
Export: tel. +39 02 953854209/221/225 fax +39 02 95380056  
gima@gimaitaly.com - export@gimaitaly.com  
[www.gimaitaly.com](http://www.gimaitaly.com)

**BILANCIA GIMA PROFESSIONAL ANALOGICA**  
**GIMA PROFESSIONAL ANALOG SCALE**  
**BALANCE ANALOGIQUE PROFESSIONAL GIMA**  
**ANALOGWAAGE „PROFESSIONAL“ GIMA**  
**BÁSCULA PROFESIONAL ANALÓGICA GIMA**  
**BALANÇA PROFISSIONAL ANALÓGICA GIMA**  
**ΑΝΑΛΟΓΙΚΗ ΕΠΑΓΓΕΛΜΑΤΙΚΗ ΖΥΓΑΡΙΑ GIMA**  
**(GIMA) "ميزان مهني قياسي " جيما"**

**MANUALE D'USO E MANUTENZIONE**  
**USE AND MAINTENANCE BOOK**  
**INSTRUCTIONS DE FONCTIONNEMENT ET ENTRETIEN**  
**BETRIEBS UND WARTUNGS ANWEISUNGEN**  
**MANUAL DE USO Y MANTENIMIENTO**  
**MANUAL DE USO E MANUTENÇÃO**  
**ΕΓΧΕΙΡΙΔΙΟ ΧΡΗΣΗΣ ΚΑΙ ΣΥΝΤΗΡΗΣΗΣ**  
**دليل الإستعمال والرعاية**

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## FEATURES

The GIMA Professional analog scale is a device for measuring the weight of the human body. The wide measurement range makes it ideal for any kind of use.

### Technical Features

Maximum capacity: 160 kg

Precision: 1 kg

**Permitted tolerance: <100 kg ±1 kg >100 kg ±1,5 kg**

Dimensions: 320 x 70 x 430 mm

Weight of unit: 2,4 kg

Not certified product class III

## PRESCRIPTIONS

*Do not use the equipment in case it is damaged. Apply to your retailer.*

*Avoid precarious repairs. Repairs shall be carried out with original spare parts only, which shall be installed according to the intended use.*

Since the product is made of corrosion-proof materials suitable for the environmental conditions foreseen for its normal use, does not require special care, however it is necessary to store it in a closed place making sure that is protected from dust and dirt to assure its hygienic conditions. Moreover, it is recommended to store the product in a place which can be reached easily by the personnel in case of necessity.

## UNPACKING

*Always remember that packing elements (paper, cellophane, stitches, adhesive tape, etc.) can cut and/or hurt if they are not carefully handled. They shall be removed with adequate means and shall not be left at the mercy of irresponsible persons; the same is valid for tools used to remove packages (scissors, knives, etc.).*

After opening the packages, first of all it is necessary to check all pieces and parts composing the product. Check that they are all present and in perfect conditions.

### OPERATION

- Position the scale on the floor in a stable position (ex.: not on a rug)
- Verify that the pointer is set to zero
- Put and remove a weight from the scale a few times and verify that the pointer returns to its position exactly on the zero, if necessary reposition the pointer using the adjustment wheel.
- Step onto the scale and verify that both feet are positioned correctly.
- Try to stand as still as possible for a precise measurement
- To obtain valid results for statistical purposes, weigh yourself every day possibly at the same time. It is best to weigh yourself without clothes and in the morning as soon as you wake up.

*Caution: the scale can hold a maximum weight of 160 kg, do not overload the scale or it could get damaged!*

### MAINTENANCE

The device was designed to last a lifetime. It is resistant to chemical attack, except for M.E.K. ketones, diluting solutions and polish removers.

### GIMA WARRANTY CONDITIONS

Congratulations for purchasing a GIMA product.

This product meets high qualitative standards both as regards the material and the production.

The warranty is valid for 12 months from the date of supply of GIMA.

During the period of validity of the warranty, GIMA will repair and/or replace free of charge all the defected parts due to production reasons. Labor costs and personnel traveling expenses and packaging not included.

All components subject to wear are not included in the warranty.

The repair or replacement performed during the warranty period shall not extend the warranty.

The warranty is void in the following cases: repairs performed by unauthorized personnel or with non-original spare parts, defects caused by negligence or incorrect use.

GIMA cannot be held responsible for malfunctioning on electronic devices or software due to outside agents such as: voltage changes, electro-magnetic fields, radio interferences, etc.

The warranty is void if the above regulations are not observed and if the serial code (if available) has been removed, cancelled or changed.

The defected products must be returned only to the dealer the product was purchased from.

Products sent to GIMA will be rejected.

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