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Card *Syphilis Ultra Rapid Test* vs. TPHA

Metodo		TPHA		Risultati Totali
		Positivo	Negativo	
Card Syphilis Ultra Rapid Test	Risultati			
	Positivo	384	2	386
	Negativo	1	493	494
Risultati Totali		385	495	880

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BIBLIOGRAFIA

1. Claire M. Fraser. *Complete genome sequence of Treponema Pallidum, the syphilis spirochete*, Science 1998; 281 July: 375-381.
2. Center for Disease Control. *Recommendations for diagnosing and treating syphilis in HIV-infected patients*,

MMWR Morb. Mortal Wkly Rep. 1988; 37: 601

3. Aral R. Marx. *Crack, sex and STD*, Sexually Transmitted Diseases, 1991; 18:92-101

4. J.N. Wasserheit. *Epidemiological Synergy: Interrelationships between human immunodeficiency virus infection*

and other sexually transmitted diseases, Sexually Transmitted Diseases 1992; 19:61-77

5. Johnson Phillip C. *Testing for Syphilis*, Dermatologic Clinic 1994; 12 Jan: 9-17.

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 Fabbricante
 Innovacon, Inc.
 9975 Summers Ridge Road
 San Diego, CA 92121, USA
 Made in P.R.C.

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Rapid test for the diagnosis of syphilis for the qualitative determination of anti-*Treponema Pallidum* (TP) antibodies (IgG and IgM) in whole blood, serum or plasma.

Only for professional in vitro diagnostic use.

INTENDED USE

The card for the *Syphilis Ultra Rapid Test* (Blood/Serum/Plasma) is a rapid immunochromatographic test for the qualitative determination of antibodies (IgG and IgM) anti-*Treponema Pallidum* (TP) in whole blood, serum or plasma to be used as guide in the diagnosis of syphilis.

SUMMARY

Treponema Pallidum (TP) is the aetiological agent of the venereal disease syphilis. TP is a spirochete,

that is corkscrew-shaped bacteria with an external sheath and a cytoplasmic membrane.

1 Not much is known about this microorganism compared to other bacteria. According to the Center

for Disease Control (CDC), the incidence of syphilis infection has increased significantly since 1985.

2 Some of the determining factors that contributed to the spreading of this disease include “the epidemic” of cocaine-crack and the high incidence of prostitution among drug addicts.

3 A study highlighted a substantial epidemiologic link between catching the disease and spreading the HIV and syphilis virus.

4 Typical features of syphilis are the multiple clinical stages and the long periods of asymptomatic latent infection. Primary syphilis is marked by the presence of a syphiloma, or luetic ulceration, in the point of inoculation. The antibody response to TP can be detected within 4-7 days after the appearance of syphiloma. The infection can be detected until the patient is not properly treated.

5 The card for the *Syphilis Ultra Rapid Test* (Blood/Serum/Plasma) uses a combination with double antigen of a particle sensitized with a syphilis antigen immobilized on the membrane for the qualitative and selective determination of the anti-TP antibodies (IgG and IgM) in whole blood, serum or plasma.

PRINCIPLE

The card for the *Syphilis Ultra Rapid Test* (whole blood, serum or plasma) is a qualitative immunologic dose that uses the membrane strip for the determination of the anti-TP antibodies (IgG and IgM) in whole blood, serum or plasma. In this test, the recombinant antigen of syphilis is immobilized on the reaction section of the strip. After adding the specimen in the specific well of the card, the specimen reacts with the particles sensitized with the syphilis antigen. The mixture migrates on the entire strip and interacts with the immobilized syphilis antigen.

The structure of the double antigen test is able to determine both IgG antibodies and IgM antibodies in the sample. If a specimen contains anti-TP antibodies, a red strip will appear in the reactive area of the membrane. This indicates a positive result. If, on the other hand, a sample does not contain anti-TP antibodies, no red strip will appear and this indicates a negative result.

To control the procedure, a red strip will appear in the control line to indicate that an adequate volume was added and that the migration on the membrane took place properly.

REAGENT

The card contains particles sensitized with syphilis antigen and syphilis antigen adherent to the membrane.

PRECAUTIONS

- Only for professional *in vitro* diagnostic use. Do not use after date of expiration.
- Do not eat, drink or smoke in the area where the samples or kit are handled.
- Do not use the test if the bag is damaged.
- Handle all the specimen as potentially infective. Take all necessary measures against microbiologic risks in all analysis phases and follow the standard procedures for the proper disposal of the samples.
- During the analysis of the samples wear protective clothing: lab coat, disposable gloves and eye protection.
- After use, the card must be eliminated according to the local regulations in force.
- Humidity and temperature may negatively affect the results.

PRESERVATION AND STABILITY

Keep sealed bag at room temperature or in the refrigerator (2-30°C). The card is stable up to the date of expiration indicated on the sealed bag. Keep card in the sealed bag up until use. **DO NOT FREEZE.** Do not use after date of expiration.

TAKING AND PREPARING THE SAMPLES

- The card for the *Syphilis Ultra Rapid Test* (whole blood, serum or plasma) can be used with whole blood (as venous blood sample or lancing device), serum or plasma.
- To take samples of **whole blood with lancing device**:
 - Wash the patient's hands with warm water and soap or clean with gauze and alcohol. Let it dry.
 - Massage the hand, without touching the point of introduction of the lancet, from the wrist towards the tip of the middle finger or the ring finger.
 - Puncture the finger with a sterile lancet. Remove the first drops of blood.
 - Delicately rub the hand from the wrist towards the palm and the finger so that a sufficient quantity of blood gathers near the point of puncture.
 - How to place whole blood sample taken with lancing device on the card through the **test tube for capillary blood taking**:
 - Touch the tip of the test tube to take approx. 50 µl of blood. Do not form any air bubbles.
 - Position the bulb on the upper tip of the test tube, then press in order to push the whole blood on the sample well (S) of the test card.
 - How to put whole blood sample directly on card, **without using the test tube**:
 - Position the patient's finger so that the drop of blood is right above the card well (S).
 - Let 2 drops of whole blood fall in the center of the well (S) on the card, or move the patient's finger so that the drop of blood touches the center of the well (S). Do not touch well (S) with finger.
 - Separate the serum or the plasma from the blood to avoid hemolysis. Only use transparent samples, not hemolyzed.
 - Perform test right after taking samples. Do not leave samples at room temperature for long periods. The samples of serum and plasma can be kept at a temperature of 2-8°C for 3 days. For longer periods of preservation, freeze samples at temperature lower than -20°C. The venous whole blood samples collected must be kept at a temperature of 2-8°C if the test is performed within 2 days from when it was taken. Do not freeze blood samples. The whole blood samples taken with the lancing device must be analyzed immediately.
- Bring the samples at room temperature before the analysis. Let the frozen samples defrost completely and mix them accurately before the analysis. Do not freeze and defrost the samples several times.
- In case of shipping, pack the samples according to the local regulations in force for the shipment of aetiological agents.

COMPOSITION OF THE PACKAGE

Material Supplied

- Card
- Dropper
- Buffer (only for whole blood)
- Methodology

Material Needed But Not Supplied

- Containers for the collection of samples
- Lancets (only for taking blood sample with lancing device)
- Spinning unit
- Timer
- Test tubes for capillary blood taking and dispenser (only for taking blood sample with lancing device)

PROCEDURE

Bring card, buffer and serum or plasma samples and/or control at room temperature (15-30°C) before the analysis.

1. Bring the bag at room temperature before opening it. Remove the card from the sealed bag and use it as soon as possible.

2. Position the card on a clean flat surface.

For the **serum or plasma samples**: keep the dropper in vertical position and **transfer 3 drops of serum or plasma** (approx. 75 µl) onto the well (S) and start timer. See illustration below.


For **venous whole blood samples**: keep the dropper in vertical position and **transfer 2 drops of whole blood** (approx. 50 µl) onto the well (S), then **add 1 drop of buffer** (approx. 40 µl) and start timer. See illustration below.

For the **whole blood samples taken using the lancing device**:

• How to use a test tube for capillary blood taking: **fill test tube and transfer approx. 50 µl of whole blood sample** taken with lancing device on the well (S) of the card, **then add 1 drop of buffer** (approx. 40 µl) and start timer. See illustration below.

• How to use the drops directly: **let 2 drops of whole blood taken using the lancing device** (approx. 50 µl) fall on the well (S) of the card, **then add 1 drop of buffer** (approx. 40 µl) and start timer. See illustration below.

3. Wait until the red strip/s appear. **Read result after 10 minutes.** Do not read result after 30 minutes.

INTERPRETATION OF RESULTS

(See previous illustration)

POSITIVE:* **Two different red strips appear**, one in the control area (C), the other in the reactive area (T).

***NOTE:** The intensity of the red color in the reactive area (T) varies according to the concentration of the anti-*TP* antibodies in the sample. Every shade of red in the reactive area (T) must be considered a positive result.

NEGATIVE: **A red strip appears in the control area (C).** No red or pink strip appears in the reactive area (T).

NOT VALID: **No strip appears in the control area.** The reason why no strip appeared in the control area is that there is an insufficient volume of sample or the procedure was not performed correctly. Go over procedure and repeat the test with a new card. If the problem continues, do not use kit and contact local distributor.

QUALITY CONTROL

Every test includes a procedure control. The red strip that appears in the control area (C) is an internal control of the procedure and indicates that a sufficient sample volume was taken and the procedure was performed correctly.

No control standards are given with the kit; we recommend analyzing a positive and negative control as a correct laboratory practice to confirm the analytical procedure and verify the correct execution of the test.

LIMITS

1. The card *Syphilis Ultra Rapid Test* (Whole Blood/Serum/Plasma) is exclusively for *in vitro* diagnostic procedures. The test must be used to determine anti-*TP* antibodies only in whole blood, serum and plasma samples. It is not able to determine the quantitative value or the percentage of increase of anti-*TP* antibodies.

2. The card *Syphilis Ultra Rapid Test* (Whole Blood/Serum/Plasma) only indicates the presence of anti-*TP* antibodies in the sample and must not be used as the only diagnosis criteria for the *TP* infection.

3. As for all the diagnostic tests, the results must be interpreted in relation to the medical history of the patient.

4. If the test is negative and the clinical symptoms persist, we recommend performing other tests with different clinical methods. A negative result does not however exclude the possibility of *TP* infection.

EXPECTED RESULTS

The card *Syphilis Ultra Rapid Test* (Whole Blood/Serum/Plasma) was compared with one of the most popular TPHA tests available on the market and it offered an accuracy equal to or higher than 99.7%.

PERFORMANCE

Clinical sensitivity, specificity and accuracy

The card *Syphilis Ultra Rapid Test* (Whole Blood/Serum/Plasma) correctly identified samples of a serum-conversion chart and it was compared to one of the most popular TPHA tests available on the market using clinical samples. The results show that the relative sensitivity of the card for the *Syphilis Ultra Rapid Test* (Whole Blood/Serum/Plasma) is equal to 99.7%, whereas relative specificity is equal to 99.6%.

Card Syphilis Ultra Rapid Test vs. TPHA

Method		TPHA		Global Results
Card Syphilis Ultra Rapid Test	Results	Positive	Negative	
	Positive	384	2	386
	Negative	1	493	494
Global Results		385	495	880

Relative sensitivity: 99.7% (98.6-100.0%)*

Relative specificity: 99,6% (98,6-100,0%)*

Accuracy: 99,7% (99,0-99,9%)*

* Allowed intervals 95%

Accuracy
Intra-Session

The intra-session accuracy was determined analyzing 10 repetitions of 4 samples: one negative, one low positive, one medium positive and one high positive. The values of the negative sample, positive low, medium positive and high positive were identified correctly in over >99% of the cases.

Inter-Session

The inter-session accuracy was determined through 10 independent analysis on the same 4 sample: one negative, one low positive, one medium positive and one high positive. Three different batches of the Card for the Syphilis Ultra Rapid Test (Whole Blood/Serum/Plasma) were analyzed in 3 days with negative samples and low, medium and high positive samples. The samples were correctly identified in >99% of the cases.

REFERENCE

1. Claire M. Fraser. *Complete genome sequence of Treponema Pallidum, the syphilis spirochete*, Science 1998; 281 July: 375-381.
2. Center for Disease Control. *Recommendations for diagnosing and treating syphilis in HIV-infected patients*, MMWR Morb. Mortal Wkly Rep. 1988; 37: 601
3. Aral R. Marx. *Crack, sex and STD*, Sexually Transmitted Diseases, 1991; 18:92-101
4. J.N. Wasserheit. *Epidemiological Synergy: Interrelationships between human immunodeficiency virus infection and other sexually transmitted diseases*, Sexually Transmitted Diseases 1992; 19:61-77
5. Johnson Phillip C. *Testing for Syphilis*, Dermatologic Clinic 1994; 12 Jan: 9-17.

Graphical Symbols used

	Attention, see instructions for use		Tests per kit		Manufacturer
	For <i>in vitro</i> diagnostic use only		Use by		Do not reuse
	Store between 2-30°C		Lot Number		Catalog 29118

 Manufacturer
 Innovacon, Inc.
 9975 Summers Ridge Road
 San Diego, CA 92121, USA
 Made in P.R.C.

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MDSS GmbH