EN

iHealth®

**Wireless Activity and Sleep Tracker
Quick Start Guide**

Welcome

Welcome to your iHealth Wireless Activity and Sleep Tracker (the iHealth Tracker), a device that tracks your daily movement and rest. This Tracker, along with the free iHealth MyVitals app, can keep you motivated and help you stay on track to a more active and healthier lifestyle.

Box Contents

1 Tracker

2 Colored wristbands and matching waist clips

1 Quick Start Guide

1 USB cable for charging





1 Sleep wristband

Wearing the Tracker

Wear your iHealth tracker on your wrist or waist during the day to track your daily activity. Slip it in the sleep wristband at night to monitor your sleep.

Your iHealth Tracker is rain, splash and sweat proof. It is not waterproof and should not be worn while swimming or showering.

Mobile Device Compatibility

The iHealth Tracker works with the following iOS version 5.0 or higher devices:

iPod Touch®
5th Generation +

iPhone 4S +

iPad mini

iPad 3rd
Generation +

Companion Mobile App

Download the free “iHealth MyVitals” app from the App store. Follow the on-screen instruction to register and set up your iHealth ID.

Access iHealth Cloud

Upon setting up your iHealth ID, you will also have access to a free, secure iHealth Cloud account by using your ID and password. Go to www.ihealthlabs.com, then click on “Sign In”.

Battery Charge

The iHealth Tracker has a built-in rechargeable battery. Charge the iHealth Tracker for approximately two hours before first use. Plug one end of the charging cable into the USB port on the back of the Tracker and the other into your computer's USB port. A fully charged battery will typically last 5–7 days.

Note: If your Tracker's battery is drained, fully charge the battery and open your app to reestablish the *Bluetooth* connection before using it again.

First Time Set-up

Your iHealth Tracker uses *Bluetooth 4.0* Low Energy (BLE) technology. Connect your iOS device to an active internet, enable the *Bluetooth* on your iOS device (Setting->General->*Bluetooth*->On), launch the app to initiate the connection, the *Bluetooth* icon will light up and stop flashing when a successful connection is established, and follow the on-screen instructions for first-time connection.


Data Syncing

Press the “Sync” button on the “activity” and “sleep” card in the app to send your stats from the iHealth Tracker to the mobile app and the Cloud. Your Tracker can save up to 14 days’ worth of stats between syncing. Daily syncing is recommended.

Activity Tracking

Your iHealth Tracker is always tracking. View your steps, calories, and distance stats by pressing the button on the iHealth Tracker. Activity stats reset to zero at midnight.

Sleep Tracking

To start tracking your sleep, press and hold the button on the Tracker for 2 seconds until you see “sleep” on the Tracker display. When you are awake, press and hold the same button for 2 seconds to end sleep tracking and resume activity tracking. Sleep stats will appear on your app after syncing.

Silent Alarms

When you set a silent alarm from the iHealth MyVitals app, your Tracker will vibrate gently on your wrist to either wake you up or remind you to be more active.

Travel Mode

By Car

If you plan to drive a long distance, you have the option to put the Tracker on the “driving” mode. The iHealth Tracker will not count any steps when “driving” mode is turned on.

By Airplane

When traveling by air, you can put the tracker on the “airplane” mode to disable *Bluetooth* connection while in flight.

To set up either of these two options, go to the iHealth MyVitals app->Settings->Device Setting. Press and hold the button on the Tracker for 2 seconds to exit and resume activity tracking.

Visit www.ihealthlabs.com for additional product information. For Customer Service, please contact:
USA

Tel: 855.816.7705 (9:00 a.m. - 5:00 p.m. PST, Monday to Friday, except holidays)

Email: support@ihealthlabs.com

Europe

Tel: +33(0)1 44 94 04 81 (9:00 a.m. - 5:30 p.m. Monday to Friday, except holidays)

Email: support@ihealthlabs.eu

IMPORTANT INFORMATION REQUIRED BY THE FCC

This device complies with part 15 of the FCC Rules. Its operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by iHealth Lab Inc could void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

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