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48-HOUR CONTINUOUS CONTROL OF BODY TEMPERATURE (N-37-38-39-40°C)

Remove the Fever Monitor from its support and apply it onto the forehead. After 30 seconds, if several numbers appear together with the letter "N", the temperature corresponds to the highest number. If a fever is read, a clinical thermometer should be used. Carry out the test at room temperature, under the supervision of an adult. Do not expose to direct sunlight and keep away from heat, alcohol and solvents. The device is disposable. Do not keep it in continuous contact with the skin for more than 48 hours. (Issue of December 7, 2005).

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48-HOUR CONTINUOUS CONTROL OF BODY TEMPERATURE (N-37-38-HF)

Remove the Fever Monitor from its support and apply it onto the forehead. After 30 seconds, if several numbers appear together with the letter "N", the temperature corresponds to the highest number. If a fever is read, a clinical thermometer should be used. Carry out the test at room temperature, under the

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supervision of an adult. Do not expose to direct sunlight and keep away from heat, alcohol and solvents. The device is disposable. Do not keep it in continuous contact with the skin for more than 48 hours. (Issue of December 7, 2005).

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