MAXIASPEED

GUARANTEED

3

YEARS

HOSPITAL

MANUALE D'USO E MANUTENZIONE
GEBRAUCHS- UND INSTANDHALTUNGSANWEISUNGEN
USE AND MAINTENANCE MANUAL
MANUEL D'UTILISATION ET D'ENTRETIEN
MANUAL DE USO Y MANTENIMIENTO

CARATTERISTICHE TECNICHE/TECHNISCHE DATEN/TECHNICAL CHARACTERISTICS/CARACTÉRISTIQUES TECHNIQUES/ CARACTERÍSTICAS TÉCNICAS

Elettrocompressore a pistone con protettore termico
Elektrischer Hubkolbenkompressor mit Wärmeschutz
Piston electro-compressor with thermal protection
Électrocompresseur à piston avec protection thermique
Electrocompresor de pistón con protector térmico

Versione pompa singola: apparecchio alto flusso, alto vuoto
Ausführung mit Einzelpumpe: Hochflussgerät, Hochvakuum
Single pump version: high flow, high vacuum appliance
Version pompe unique: appareil à débit élevé, vide élevé
Versión con bomba única: aparato de alto flujo, alto vacío

Versione doppia pompa: apparecchio alto flusso, alto vuoto
Ausführung mit Doppelpumpe: Hochflussgerät, Hochvakuum
Double pump version: high flow, high vacuum appliance
Version double pompe: appareil à débit élevé, vide élevé
Versión con doble bomba: aparato de alto flujo, alto vacío

Fusibile/Sicherung/Fuse/Fusible/Fusible: T2,5A – 250V

- Tensione di alimentazione versione pompa singola/Anschlussspannung der Ausführung mit Einzelpumpe/Single pump version power supply voltage/Tension d'alimentation version pompe unique/Tensión de alimentación en versión con bomba única: **230V ~ 50Hz 300VA**
- Tensione di alimentazione versione doppia pompa/Anschlussspannung der Ausführung mit Doppelpumpe/Double pump version power supply voltage/Tension d'alimentation version double pompe/Tensión de alimentación en versión con doble bomba: **230V ~ 50Hz 330VA**
- Livello di vuoto regolabile/Regelbarer Vakuumpegel/Adjustable vacuum level/Niveau de vide réglable/Nivel de vacío regulable: **0 ± -0,90 bar (90 kPa)**
- Portata aria versione pompa singola/Saugleistung der Ausführung mit Einzelpumpe/Single pump version air flow/Débit d'air version pompe unique/Caudal de aire de versión con bomba única: **60 l/min.**
- Portata aria versione doppia pompa/Saugleistung der Ausführung mit Doppelpumpe/Double pump version air flow/Débit d'air version double pompe/Caudal de aire de versión con doble bomba: **90 l/min.**
- Dimensioni/Abmessungen/Dimensions/Dimensions/Dimensiones: **470 x 580 x 560 mm**
- Peso versione pompa singola/Gewicht der Ausführung mit Einzelpumpe/Single pump version weight/Poids version pompe unique/Peso de la versión con bomba única: **15 kg**
- Peso versione doppia pompa/Gewicht der Ausführung mit Doppelpumpe/Double pump version weight/Poids version double pompe/Peso de la versión con doble bomba: **17 kg**
- Rumorosità/Betriebsgeräusch/Noise level/Niveau sonore: **54 dBA**
- Classe di rischio secondo la Direttiva 93/42/EEC/Gefährdungsklasse nach Richtlinie 93/42/EEC/Hazard class according to Directive 93/42/EEC/Classe de risque selon la Directive 93/42/CEE/Clase de riesgo según la Directiva 93/42/EEC: **IIa**
- Condizioni di esercizio: Temperatura/Betriebsbedingungen: Temperatur/Operating conditions: Temperature/Conditions de service : Température/Condiciones de ejercicio: Temperatura **min. 0°C max 40°C** – Umidità aria/Luftfeuchtigkeit/Air humidity/Humidité de l'air/Humedad del aire: **min. 10% max 95%**
- Condizioni di conservazione: Temperatura/Lagerbedingungen: Temperatur/Storage conditions: Temperature/Conditions de conservation : Température/Condiciones de conservación: Temperatura **min. -10°C max 50°C** – Umidità aria/Luftfeuchtigkeit/Air humidity/Humidité de l'air/Humedad del aire: **min. 10% max 95%**
- Pressione atmosferica di esercizio/conservazione/Luftdruck Betrieb/Lagerung/Operating/storage atmospheric pressure/Pression atmosphérique de service/conservation/Presión atmosférica de ejercicio/conservación: **min. 690 hPa max 1060 hPa**

Apparecchio di tipo B/Gerät des Typs B/Type B appliance/Appareil de type B/Aparato de tipo B

Attenzione leggere attentamente le istruzioni d'uso!/Achtung! Die Gebrauchsanweisungen vor der Benutzung lesen!/Read the use instructions carefully!/Attention, lire attentivement les instructions d'utilisation!/Atención! leer atentamente las instrucciones de uso!

Apparecchio di classe II/Gerät der Klasse II/Class II appliance/Appareil de classe II/Aparato de clase II

Interruttore generale acceso/Hauptschalter Ein/Main switch on/Interrupteur général allumé/Interruptor general encendido

Interruttore generale spento/Hauptschalter Aus/Main switch off/Interrupteur général éteint/Interruptor general apagado

Interruttore pedale acceso/Pedalschalter Ein/Pedal switch on/Interrupteur pédale allumé/Interruptor de pedal encendido

Interruttore pedale spento/Pedalschalter Aus/Pedal switch off/Interrupteur pédale éteint/Interruptor de pedal apagado

Corrente alternata/Wechselstrom/Alternate current/Courant alternatif/Corriente alterna

Non utilizzare l'apparecchio mentre si fa il bagno o la doccia/Das Gerät nicht in der Badewanne oder beim Duschen benutzen/Do not use the appliance while taking a bath or shower/Ne pas utiliser l'appareil lorsqu'on prend un bain ou une douche/No utilizar el aparato al bañarse o ducharse

Monouso/Einweg/Disposable/Jetable/Desechable

Fusibile di protezione/Sicherung/Protection fuse/Fusible de protection/Fusible de protección

Ciclo di funzionamento: **Uso continuo**/Betriebsart: **Dauerbetrieb**/Operating cycle: **Continuous use**/Cycle de fonctionnement: **utilisation continue**/Ciclo de funcionamiento: **Uso continuo**

Grado di protezione/Schutzart/Protection level/Degré de protection/Grado de protección: **IP20**

0434 Conforme alla Direttiva 93/42/EEC/Entspricht Richtlinie 93/42/EEC/Conforms to Directive 93/42/EEC/Conforme à la Directive 93/42/CEE/Conforme a la Directiva 93/42/EEC

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The MAXIASPEED aspirator is a professional appliance suitable for hospitals, tracheotomy patients or for surgical applications. It is equipped with 4 anti-static wheels (8-9), two of which with a brake device (9), a telescopic pulling handle (2), vacuum regulator (3), vacuum gauge (1), two containers of 2 L. or 4 L. (7) protection device (6) against liquid entering the vacuum pump which interrupts the aspiration flow. It is lubrication free, easy to handle, simple to use, reliable, strong and quiet. The MAXIASPEED aspirator also has a switch (4) for switching the aspiration flow from one container to the other, very easily and without having to touch any possibly contaminated parts.

The MAXIASPEED aspirator has the following 3A accessories: 2 Containers – either 2 L. or 4 L. with protection device, power supply cable, 2 silicon connection tubes which can be sterilised Ø 8x14 - short, 2 silicon connection tubes which can be sterilised Ø 8x14 - long, 2 aspiration fittings, 2 anti-bacteria filters, pedal switch (optional).

N.B.: Only use original 3A accessories.

IMPORTANT WARNINGS

This is a medical device and it must be used by authorised personnel. It must be operated as instructed in this use manual. It is important for the operator to read and understand the use and maintenance information of the appliance. Contact your stockist for any questions you may have. MICROBIAL CONTAMINATION: where pathologies are present with microbial infection and contamination risks, perform the cleaning and disinfecting operations scrupulously after every use on the patient.

The manufacturer ensures that everything necessary is done to ensure that every appliance is of the highest quality and safety. Nevertheless, basic safety rules must always be followed, as for every electrical appliance.

- The appliance must only be used by specialised personnel.
- This appliance must be exclusively used for the purpose for which it was designed; in particular as a surgical aspirator. All other uses are improper and dangerous, and the manufacturer cannot be held responsible for any consequences deriving from improper use.
- Never use adaptors for power supply voltages different to those shown on the plate data label on the back of the appliance.
- Keep the cable far away from hot surfaces.
- The appliance is not suitable for use in the presence of anaesthetic mixes inflammable with air, or oxygen, or nitrogen oxide.
- Never handle the power supply plug with wet hands.
- Never leave the appliance close to water; do not submerge it in liquids, do not wet it. If this should occur, pull the plug out of the power supply socket immediately, before touching the appliance. Do not use the appliance if the plug or power supply cable are damaged, or wet (send it to your stockist immediately).
- The appliance casing is not protected against the entry of liquids.
- Always take out the power supply plug immediately after use.
- Maintenance and/or repair must only be performed by authorised personnel. Unauthorised repairs invalidate the guarantee.
- Ensure that the connections and container closing is performed carefully to prevent aspiration leakages.
- Do not upset the containers while they are connected to the appliance when running, because the liquid comes into direct contact with the hydrophobic anti-bacteria filter and immediately blocks aspiration; if this should happen, empty the container and replace the anti-bacteria filter.
- The appliance is equipped with a safety fuse, easily inspected in the event of a fault, located on the back of the appliance in the panel socket. Pull out the power supply plug before inspection.
- The anti-bacteria filter is a **disposable** element which must be replaced after every application.
- Do not lift the appliance using the transportation handle in any circumstance.
- Before each use, follow the cleaning and disinfection operations as shown in the "Cleaning and disinfection operations" paragraph of this instruction manual.

USAGE INSTRUCTIONS

Before each use, ensure that all the accessories are perfectly clean according to the instructions shown in the "CLEANING AND DISINFECTION OPERATIONS" paragraph.

1. Connect the appliance as illustrated in figures 2 and 2a: connect the ends of the short silicon tubes which can be sterilised respectively to the "INLET" connectors of the aspirator and to the "VACUUM" connector of the containers. Connect the long tube to the "PATIENT" connector of the selected plug.
2. Operation of the telescopic handle: for space reasons the aspirator is supplied with the telescopic handle lowered; to raise it, pull it up until the travel end hooks engage. To put the telescopic handle down again, press the travel end hooks on the external sides of the telescopic tubes with the tool provided (Fig. 2b) and push the handle downwards until it closes on the aspirator casing.
3. Using the vacuum regulator (3) it is possible to set the desired vacuum value (bar). By turning the

knob towards "+" a higher vacuum is obtained, while turning it towards "-" a lower vacuum is obtained: these values can be read on the vacuum gauge (1).

4. Start the appliance by switching the main switch to the "I" (ON) position (10). If the appliance is equipped with the pedal switch, to use the same, insert the loose pedal plug into its panel socket on the back of the appliance; switch the switch to the "I" (ON) (11) position located to the side of the main switch (10) which permits switching on/off using the pedal switch.
5. After the application, switch of the appliance, remove the power supply cable from the current socket and perform the cleaning operations as illustrated in the "CLEANING AND DISINFECTION OPERATIONS" paragraph.
6. N.B.: there is a very capacious compartment on the back of the appliance for putting the power supply cable back in after the application, and if the appliance is equipped with the pedal switch, also for this latter.

CLEANING AND DISINFECTION OPERATIONS


Detach the cannula, the manual regulator and the silicon tubes which can be sterilised $\varnothing 8 \times 14$ - long from the container. Disconnect the silicon tubes which can be sterilised $\varnothing 8 \times 14$ - short from the respective containers and from the respective anti-bacterial filters and empty the container/s. We recommend emptying the liquids collected into the biological material container.

In order to clean the accessories, rinsing them under hot running water is sufficient. Before performing the container cleaning operations empty their content into the biological material containers and clean them by sterilisation in an autoclave with a cycle at 121°C . The silicon aspiration tubes can be sterilised in an autoclave with a cycle at 121°C .

During the sterilisation operations the operator must use latex gloves and suitable clothing in order not to come into contact with any possibly contaminating substances.

The disposable anti-bacteria filter must be replaced after every application.

Never wash the appliance under water or by immersion; clean the outside of the appliance using a cloth moistened with detergent (non abrasive).

Socket for pedalswitch (where foreseen).
Cable and pedal switch compartment
Power supply socket with fuse (inspection type)



PERIODICAL CONTROL FOR APPLIANCE SAFETY

MAXIASPEED does not require maintenance and/or lubrication. Nevertheless, several simple checks must be performed before every use:

- Check the integrity of the casing, power supply cable and the pedal switch, if supplied.
- Block the aspiration connector with a finger and check that the vacuum level reaches -0.9 bar.
- Check that there are no disturbing noises suspicious of malfunctioning.
- Check that the cage is positioned correctly inside its housing. It must be in axis with the plug aspiration hole so that the floating valve can block it in the circumstance in which the aspirated liquid exceeds the maximum volume the container can receive. Do not connect the cage as shown in figure 3a; it must be connected as shown in figure 3b.



CAGE OUT OF AXIS

CAGE IN AXIS

- Check that the float is fitted in the correct position and that it is free to move inside the cage (dirt or encrustation may hinder its movement). Figure 4 shows the correct insertion.
- Check that the plug seal is inserted with the edging as illustrated in figure 5.

Electromagnetic compatibility – additional notes in conformity with the IEC 60601-1-2/A1:2004-09 Standard.

Emission		
Emission test	Conformità	Electromagnetic Environment - guidance
RF Emissions Cispr 11	Group 1	MAXIASPEED uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions Cispr 14	--	The suction equipment is suitable for use in all establishments and those directly connected to the public low-voltage power supply network that supplies buildings used to domestic purposes. It is possible to use the device in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, adopting provisions in the installation, such as a longer distance from potentially sensitive appliances.
Harmonic Emissions IEC 61000-3-2	Classe A Conforms	It is possible to use the device in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings user for domestic purposes.
Emissioni di fluttuazioni di tensione/flicker IEC 61000-3-3	Conforme	

Immunity aspects

MAXIASPEED is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity test	Test level EN 60601-1-2	Compliance Level	Electromagnetic environment - guide
Electrostatic discharge (ESD) EN 61000-4-2	± 6kV contact ± 8kV air	± 6kV contact ± 8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Burst/Fast Transient EN 61000-4-4	±2kV power supply lines	±2kV power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge EN 61000-4-5	±1kV differential mode ±2kV common mode	±1kV differential mode ±2kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines EN 61000-4-11	< 5% UT (>95% dip in UT) for 0,5 cycles 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles < 5% UT (>95% dip in UT) for 5 seconds	< 5% UT (>95% dip in UT) for 0,5 cycles 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles < 5% UT (>95% dip in UT) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions. It is recommended that the Laser/Pulsed Light BiosYag be powered from an uninterruptible power supply or a battery.
Power frequency magnetic field EN 61000-4-8	3 A/m	3 A/m	Magnetic power frequency fields should be that of a typical commercial or hospital environment.

Immunity aspects at r.f.

MAXIASPEED is intended for use in the electromagnetic environment specified below. The customer or the user of the navigator should assure that it is used in such an electromagnetic environment.

Immunity test	Test Level EN 60601-1-2	Compliance Level	Electromagnetic environment - guide
RF conducted EN 61000-4-6	3 Veff from 150kHz to 80MHz	3 Veff from 150kHz to 80MHz	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from that equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1,2 \sqrt[3]{P}$ da 150kHz a 80MHz $d = 1,2 \sqrt[3]{P}$ da 80 MHz a 800 MHz $d = 2,3 \sqrt[3]{P}$ da 800 MHz a 2,5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).
RF radiated EN 61000-4-3	3 V/m from 80MHz to 2,5GHz	3 V/m from 80MHz to 2,5GHz	

Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:

Recommended separation distances between portable and mobile RF communications equipment and the Laser/Pulsed Light Bios Yag

MAXIASPEED is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of the transmitter (m)		
	From 150kHz to 80MHz $d = 1,2 \cdot \sqrt{P}$	From 80MHz to 800MHz $d = 1,2 \cdot \sqrt{P}$	From 800MHz to 2,5GHz $d = 2,3 \cdot \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer

Notes:

(1) At 80 MHz and 800 MHz the separation distance for the higher frequency range applies.

(2) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

PROBLEMS, CAUSES AND SOLUTIONS

PROBLEMS	POSSIBLE CAUSES	SOLUTIONS
Excessive noise.	Damaged pump or blockages in the internal aspiration circuit	Send to the assistance circuit.
The unit switches on but does not aspirate	- Damaged pump - Vacuum regulator fully open. Connection tubes disconnected and/or badly connected, broken connection tubes. Container not in a vertical position, full, or defective overflow valve. Possible blockage of the hydraulic circuit inside the unit.	- Send to the assistance circuit. - Check the position of the vacuum regulator. Check the connections and the integrity of the tubes. Position the container in a vertical position, check the overflow valve (blocked) and/or replace the silicon tubes.
The vacuum rate cannot be regulated	Damage to the internal hydraulic system or blockage of the connection tubes to the aspiration unit.	Send to the assistance circuit.
When the appliance is switched on, the protection fuse always trips	Pump probably damaged or in short-circuit.	Send to the assistance circuit.
The vacuum gauge does not work	Liquid penetrating the pneumatic circuit.	Send to the assistance circuit.

Note: if you experience faults or malfunctioning problems different to those listed above, always and exclusively contact authorised assistance centres.

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GB

1. Vacuum gauge
2. Telescopic transport handle
3. Vacuum regulator
4. Container selector
5. Anti-bacterial filter
6. Protection device
7. Containers (2000cc or 4000cc)
8. Anti-static wheels without brake
9. Anti-static wheels with brake
10. Main switch
11. Operating switch with pedal

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INSTRUCTIONS FOR CORRECT PRODUCT DISPOSAL ACCORDING TO EUROPEAN DIRECTIVE 2002/96/EC.

At the end of its useful life the product must not be disposed of together with urban refuse.

It may be consigned to specialised differentiated collection sites provided by municipal authorities, or to stockists providing this service.

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This warranty certificate is valid only if returned to your dealer along with Receipt and Faulty Unit. Warranty does not cover damages caused by misuse, crashes or lack of attention.

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- 3A Health Care S.r.l. non risponde di ulteriori estensioni del periodo di garanzia assicurate da terzi.