GIMA

PROFESSIONAL MEDICAL PRODUCTS

Gima S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) Italy
Italia: tel. 199 400 401 - fax 199 400 403
Export: tel. +39 02 953854209/221/225 fax +39 02 95380056
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

MONITOR OSTETRICO

PC-8000Pro

OBSTETRIC MONITOR

PC-8000Pro

Manuale d'uso
User manual

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The Manual describes, in accordance with the PC-8000Pro Obstetric Monitor's features and requirements, main structure, functions, specifications, correct methods for transportation, installation, usage, operation, repair, maintenance and storage, etc. as well as the safety procedures to protect both the user and equipment. Refer to the respective chapters for details.

The Manual is published by Shenzhen Creative Industry Co., Ltd (Creative) and Creative has the ultimate right to explain the Manual. No part of this manual may be photocopied, reproduced or translated into another language without the prior written consent. Creative reserves the right to improve and amend it at any time without prior notice. Amendments will however be published in a new edition of this manual.

Version of This Manual: Ver1.0

May 5, 2009

All rights reserved.

Marks in the Manual:

-  Warnings: must be followed to avoid endangering the operator and the patient.
-  Notes: some important information and tips about operations and application.
-  Attentions: must be followed to avoid causing damage to the monitor.

Instructions to User

Dear Users,

Thank you very much for purchasing our product. Please read the following pages very carefully before using this equipment.

Read these instructions carefully before using this monitor. These instructions describe the operating procedures to be followed strictly. Failure to follow these instructions can cause monitoring abnormality, equipment damage and personal injury. The manufacturer is NOT responsible for the safety, reliability and performance issues and any monitoring abnormality, personal injury and equipment damage due to user's negligence of the operation instructions. The manufacturer's warranty service does not cover such faults.

- Do NOT operate this monitor in the circumstance where flammable gas or explosion factor exists.
- To ensure patient safety, do NOT place the monitor in any position that might cause it to fall on patient.
- The user must check that the equipment functions safely and ensure that it is in proper working condition before being used.
- The device connected to this monitor should be in accordance with IEC 60601-1.
- Please shut down the monitor and remove sensors while during MRI scanning, or else induced current could potentially cause burns. The monitor may affect the MRI image, and the MRI unit may affect the accuracy of monitor measurements.
- Although biocompatibility tests have been performed on all the applied parts, some exceptional allergic patients may still have anaphylaxis. Do NOT apply to those who suffer from anaphylaxis.
- All the cables and rubber tubes of the applied parts should be kept away from the patient's cervix to prevent any possible suffocation of the patient.
- All the parts of the monitor should NOT be replaced at will. If necessary, please use the components provided by the manufacturer or those of the same model and standards as the accessories along with the monitor which are provided by the same factory, otherwise negative effects concerning safety and biocompatibility, etc. may be caused.
- If the monitor falls off accidentally, please do NOT operate it until its safety and technical indexes have been carefully tested and positive testing results obtained.
- Before maintenance, please switch off power.
- This monitor is a professional medical device, can only be operated by trained personnel with qualification.
- Please peruse the relative content about the clinical restrictions and contraindication.
- Dispose of the expired device and its accessory according to applicable local regulations.

Table of Contents

CHAPTER 1	OVERVIEW.....	36
1.1	FEATURES.....	36
1.2	PRODUCT NAME AND MODEL	36
1.3	MAJOR APPLICATIONS AND SCOPE OF APPLICATION	36
1.4	OPERATING ENVIRONMENT	36
1.5	IMPACT ON THE ENVIRONMENT AND RESOURCE.....	36
1.6	SAFETY	36
CHAPTER 2	WORKING THEORIES	37
2.1	OVERALL STRUCTURE	37
2.2	COMPOSITION	37
2.3	WORKING THEORIES	37
CHAPTER 3	INSTALLATION AND CONNECTION.....	38
3.1	INSTALLATION.....	38
	3.1.1 Opening the Box and Check	38
	3.1.2 Connecting the AC Power Cable	38
	3.1.3 Starting the Monitor	38
3.2	APPEARANCE	39
	3.2.1 View of the Front Panel	39
	3.2.2 Left and Right Side View.....	40
	3.2.3 Rear View	40
	3.2.4 Bottom View	41
3.3	CONNECTION	41
	3.3.1 Ultrasonic Transducer Connection	41
	3.3.2 TOCO Transducer Connection	42
	3.3.3 FM Marker.....	42
	3.3.4 ECG Connection.....	42
	3.3.5 Blood Pressure Cuff Connection	44
	3.3.6 To connect the SpO ₂	45
	3.3.7 TEMP Transducer Connection.....	46
	3.3.8 Loading Printing Paper and Hanging Installation	47
CHAPTER 4	MONITORING SCREEN.....	47
4.1	PREPARATION PRIOR TO USE.....	47
4.2	INITIAL DISPLAY SCREEN	47
4.3	FREEZE AND REPLAY.....	49
4.4	SYSTEM MENU.....	49
	4.4.1 System Settings	50
	4.4.2 Color Settings	51
	4.4.3 Alarm Settings	51
	4.4.4 NIBP List.....	52
CHAPTER 5	TECHNICAL SPECIFICATIONS.....	52
5.1	FETAL MONITORING.....	52
5.2	ECG/RESP MONITORING	53
5.3	NIBP MONITORING.....	53
5.4	SpO ₂ MONITORING	53
5.5	OTHER TECHNICAL INFORMATION AND PARAMETERS	54
CHAPTER 6	DIMENSIONS AND WEIGHT	54
6.1	OVERALL DIMENSIONS AND WEIGHT	54
6.2	ACCESSORIES	54
CHAPTER 7	OPERATING PRINCIPLE OF THE MONITOR	55
7.1	THE PRINCIPLE FOR MEASURING THE FETAL HEART RATE	55
7.2	THE PRINCIPLE FOR MEASURING THE UTERINE PRESSURE AND FETAL MOVEMENT COUNT	55
7.3	HOW TO OBTAIN HIGH QUALITY ECG AND ACCURATE HEART RATE VALUE	55

7.4	NIBP MONITORING.....	56
	7.4.1 The Principle for Measurement of Blood Pressure	56
	7.4.2 Points to be Noted in Measurement of Blood Pressure	56
	7.4.3 Contraindication	56
7.5	SpO ₂ /PULSE MONITORING	57
	7.5.1 SpO ₂ Monitoring	57
	7.5.2 Points to be Noted in SpO ₂ and Pulse Measuring	57
CHAPTER 8	TROUBLESHOOTING	58
8.1	NO DISPLAY ON THE SCREEN	58
8.2	EXCESSIVE ECG SIGNAL INTERFERENCE OR TOO THICK BASELINE	58
8.3	NO BLOOD PRESSURE AND PULSE OXYGEN MEASURES.....	58
8.4	ABNORMAL FHR DETECTION	58
8.5	ABNORMAL UTERINE CONTRACTION PRESSURE DETECTION	59
CHAPTER 9	MAINTENANCE AND SERVICE.....	59
9.1	TECHNICAL MAINTENANCES	59
9.2	USER MAINTENANCES	60
9.3	CLEANING, STERILIZATION AND DISINFECTION	60
9.4	CLEANING, STERILIZATION AND DISINFECTION OF ACCESSORIES.....	60
9.5	STORAGE.....	60
9.6	TRANSPORTATION	60

Chapter 1

OVERVIEW

1.1 Features

This monitoring system can be used to monitor fetal heart rate, uterine contraction pressure, fetal movement count, mother's ECG, respiratory value, body temperature value, non-invasive blood pressure value, pulse oxygen saturation value, pulse rate and other physiological parameters.

- 2.1" high resolution color LCD display of fetal heart rate curve, uterine contraction pressure curve, and marks of fetal movement;
- Manual fetal movement marking function;
- Display of mother's ECG wave, respiratory wave and pulse oxygen saturation wave;
- Complete ECG monitoring through selection by leads;
- Waveform frozen function;
- Unique pulse oxygen saturation measuring device is used to ensure accuracy of the pulse oxygen and pulse rate measures;
- Flexible upper and lower alarm limit setting functions;
- Blood pressure may be measured in the mode of "adult/infant/Neonate", which may be selected via the menu, to better suit the adult, infant and neonate patient;
- Internal printer is optional;
- The modular equipment design allows the user to choose different monitoring parameters and form a monitor according to his own need;

1.2 Product Name and Model

Product Name: Obstetric Monitor
Model: PC-8000Pro

1.3 Major Applications and Scope of Application

This PC-8000Pro Obstetric Monitor is a multi-functional apparatus designed for monitoring the comprehensive physiological indices of the mother/fetus. With the functions of real-time recording and display of fetal heart rate, uterine contraction pressure and fetal movement change, it allows comprehensive analysis of fetal physiological and pathological conditions, thus providing the basis for judging the labor mode and timing. In the meantime, it may be used to observe a pregnant woman's health status by testing her physiological indices. This apparatus is applicable for use in major, medium and small-size hospitals, clinics, practitioner's office and patient's home. The operation should be performed by qualified professionals only.

1.4 Operating environment

1. Ambient temperature range: 5°C~40°C
Relative humidity: 30%~80%
Atmospheric pressure: 70 kPa ~106kPa
2. This apparatus should be situated in a place protected against direct sunlight, so as to prevent too high temperature inside it.
3. Do not use this apparatus in an environment with toxic or inflammable gas.
4. This apparatus should be fixed on a stand, so as to prevent possible shock.
5. Do not use this apparatus in combination with any equipment other than those expressly permitted in these instructions.

1.5 Impact on the Environment and Resource

Low

1.6 Safety

Class I Type BF, CF

Chapter 2

WORKING THEORIES

2.1 Overall Structure

The overall structure of this monitor is shown as Fig.2.1.


Figure 2.1 Overall Structure

2.2 Composition

The monitor consists of the main unit and the corresponding functional components (ECG leads, non-invasive blood pressure cuff, SpO₂ probe, temperature transducer, ultrasonic transducer, uterine contraction pressure probe and FM marker).

2.3 Working Theories

PC-8000Pro Obstetric Monitor, which performs physiological parameter measurement through different modules, is a product of module design. It consists of ECG module, NIBP module, SpO₂ module, fetal monitoring module and the main unit etc.

1. ECG module collects heart rate, respiration waveforms through ECG leads and collects temperature data through the temperature probes as well.
2. SpO₂ module collects the data of pulse rate, pulse oxygen saturation (SpO₂) and SpO₂ volume waveform via SpO₂ probe.
3. NIBP module collects blood pressure data, including diastolic, systolic and mean arterial pressure through NIBP cuff. The cuffs are designed for adult, infant and neonate respectively, and NIBP measurement has three modes: adult, infant and neonate.
4. Fetal monitoring module acquires fetal heart rate, uterine contraction value and the number of fetal movement through ultrasonic transducer, uterine contraction pressure probe and FM marker respectively.
5. The main unit consists of main board, multi-function board, and keyboard. Multi-functional board performs data communication among main board, ECG module, SpO₂ module, NIBP module and fetal monitoring module.

Chapter 3

INSTALLATION AND CONNECTION

3.1 Installation

3.1.1 Opening the Box and Check

1. Open the packaging, take out the monitor accessories from the box carefully and place it in a safe stable and easy to watch position.
2. Open the users' manual to sort the accessories according to the packing list.
 - Inspect the accessories for any mechanical damages
 - Check all the exposed leads and inserted accessories

Please contact the local dealer or our company in case of any problems. We are to offer you the best solution for your satisfaction.

3.1.2 Connecting the AC Power Cable

Connecting procedures:

- Make sure that the AC power supply is (100-250)VAC, (50/60)Hz.
- Use the power cable prepared by the manufacturer. Insert one end of it to the power port of the monitor and the other end to the grounded three-phase power jack.
- Connect the monitor to the grounding port with the provided ground cable.

Caution: ensure that the monitor is grounded correctly.

The provided battery of the monitor must be recharged after transportation or storage. So if the monitor is switched on without being connected to the AC power socket, it may not work properly due to insufficient power supply.

3.1.3 Starting the Monitor

The system performs self-test and enters initial display after switching on the monitor, and the orange alarm indicator blinks to inform that the user can begin operating it.

- Check all the applicable functions to make sure that the monitor works normally.
- If the built-in battery is applied please recharge it after using the monitor to ensure sufficient power storage.
-  Do not use the monitor to monitor the patient if there are indications of damage or reminders of error. Please contact the local dealer or our company.
After the monitor is switched off, it's recommended to delay 1 minute to start it again.

3.2 Appearance

3.2.1 View of the Front Panel


Figure 3.1 Front View

1. **TFT screen**
2. **Transducer holder hole**
3. **“ON/OFF”**: power switch. Press it to turn the monitor on/off.
4. **“+”**: press it to increase the volume;
5. **“-”**: press it to decrease the volume;
6. **“FHR1/FHR2”**: shift between the FHR sound of FHR1 and the FHR sound of FHR2 when monitoring twins. (Optional configuration)
7. **“PRINT”**: print data;
8. **“ZERO”**: press this key, the TOCO will resume to initial pressure. During monitoring, press “ZERO” key to zero the TOCO.
9. **Indicators**: “ON” indicates power indicator; “CHARGE” indicates battery charging indicator; “ALARM” indicates alarm indicator; “PRINT” indicates printing indicator
10. **“SILENCE”**: press it for suspend/active the alarm sound.
11. **“FREEZE”**: press it for freeze the waveforms; press it again for the real-time waveforms.
12. **“EVENT MARKER”**: doctor or nurse can use this key for marking important event or special event. The icon “↓” will be shown in the upper area of the TOCO curve.
13. **“NIBP”**: start/stop NIBP measurement.
14. **“INFO”**: press this button for 6 times and then rotate the navigation knob to change the network No..
15. **Navigation Knob**: it is the main operating key; for select system parameter and change its setting.

3.2.2 Left and Right Side View


Figure 3.2 Left Side View


Figure 3.3 Right Side View

The mains power connector is on the left panel, as shown in Figure 3.2; Functional connectors are on the right panel, as shown in Figure 3.3. The following are the details:

1. **Loudspeaker**
2. **“FM”**: fetal movement marker connector;
3. **“TOCO”**: TOCO connector;
4. **FHR sensor connector**
5. **ECG cable connector**
6. **Temperature sensor connector**
7. **SpO₂ sensor connector**
8. **NIBP cable connector**

3.2.3 Rear View


Figure 3.4 Rear View

1. **VGA connector**
2. **Network and USB connector**
3. **Handle**
4.  : Equal potential grounding terminal
5. **Fan**

3.2.4 Bottom View


Figure 3.5 Bottom View

1. **Fixing hole for hanging**
2. **Battery cover**
3. **Fuse holder: 2×T2L250V Φ 5×20mm**
4. **Icon and Symbol Instructions**

with type BF applied part

Caution! Refer to this manual.

3.3 Connection

3.3.1 Ultrasonic Transducer Connection

Connect the ultrasonic transducer to the yellow connector labeled “US” on the right panel, make sure it is ready and then follow the following procedure to fix the transducer.

1. Let the pregnant woman lie low or semi-lie.
2. Tie the abdominal belt on the patient, as shown in Figure 3.6.
3. Use stethoscope or feel out the fetal contour with hand to find the approximate position of fetal heart.
4. Spread some gel on the detecting area on the abdomen of the patient and on the front surface of the transducer.
5. Move the transducer slowly on the surface of the abdomen of the patient to trace FHB (fetal heart beat), and then adjust the transducer site, angle and the volume.
6. Fix the transducer on the belt where a clear sound can be sensed.


Figure 3.6 Ultrasonic and TOCO Transducer Location

Notice for FHR Monitoring:

- 🔊 Maybe it is hard to find fetal heart site for too obese pregnant women or other pregnant women with special disease, in this occasion please choose other measuring methods.
- 🔊 During long time monitoring lying in the supine position should be avoided owing to the possibility of supine hypotension. Sitting up or lateral position are preferable and may be more comfortable to the mother.
- 🔊 This monitor's ultrasonic emission power is lower than the regulation, but in order to ensure the pregnant women and fetus' health, it is suggested that the monitoring time is as short as possible when meeting the measuring requirement.
- 🔊 The ultrasonic transducer is nonwatertight, so it can not be used while the pregnant woman delivers in the water.
- 🔊 The ultrasonic transducer has ceramic wafer, please avoid vibrating the transducer and make the transducer fall down.
Do not start monitoring unless an audible fetal heart signal is present.
- 🔊 Fetal movement may affect the measuring result, please according to the actual condition check and adjust the transducer site during measurement.
- 🔊 Positions with strong placental sounds (swishing) or umbilical pulse sounds (indistinct pulse sounds at umbilical cord) should be avoided.
- 🔊 If the fetus is in the cephalic presentation and the mother is supine, the clearest heart sound will normally be found on the middle below umbilicus.
- 🔊 When monitoring the woman who has baby less than 28 weeks, the condition hard to measure may occur due to weak signal.
- 🔊 It is not recommended to monitor the pregnant women in the third pregnant phase.
Because at that time TOCO period is long and TOCO intensity is strong, continual fetus movement may bring the difficulty of monitoring FHR.

3.3.2 TOCO Transducer Connection

1. Connect the TOCO transducer to the connector labeled "TOCO" on the right panel.
2. Tie the abdominal belt on the patient and place the TOCO transducer in the middle position below navel, as shown in Figure 3.6.
3. Fix the TOCO transducer on the belt;
4. Press "ZERO" key to set the current pressure as the initial pressure when there is no uterine constriction.

Notice for TOCO Monitoring:

- 🔊 Press the central part of the TOCO transducer softly when necessary and observe the change of monitored TOCO values to see if it works normally.
Ultrasonic coupling gel should not be used on TOCO transducer or its detecting area.
The TOCO transducer is nonwatertight, so it can not be used while the pregnant woman delivers in the water.

3.3.3 FM Marker

1. Connect the FM marker to the connector labeled "FM" on the right panel.
2. When pregnant woman senses fetal movement, she can press FM marker to mark the movement, meanwhile the monitor will beep for prompting mark is made successfully (If the pregnant woman press the FM marker several times in 5 minutes, only one will be recorded.)

Note: the FM marker is a patient operated hand-held push key. The mother is usually instructed to press this key and hold it for 1 second when she senses fetal movement. The system will record it as one fetal movement and a mark of solid triangle mark "▲" shall be shown in the lower area of the TOCO waveform channel.

3.3.4 ECG Connection

ECG measurement is to collect the ECG signal via the ECG electrodes. Electrode connects the patient and the lead. The lead connects the monitor. The locations of the electrodes are very important for obtaining accurate ECG signals.

1. Connect the cable to the right-panel connector marked with the ECG icon.
2. Select electrodes to be used. Use only one type of electrode on the same patient to avoid variations in electrical resistance. For ECG monitoring, it is strongly recommended to use silver/silver chloride electrodes. When dissimilar metals are used for different electrodes, the electrodes may be subject

to large offset potentials due to polarization. Using dissimilar metals may also increase recovery time after defibrillation.

3. Prepare the electrode sites according to the electrode manufacturer's instructions.

The locations of the electrode are in the following figure:


Figure 3.7 Electrode Location

Note: If skin rash or other unusual symptoms develop, remove electrodes from patient.

4. After starting the monitor, if the electrodes become loose or disconnected during monitoring, the system will display "LEAD OFF" on the screen to alarm the operator.
5. The ECG leads and their corresponding locations are as follows:

Symbol		Position
RA		The intersection between the centerline of the right clavicle and Rib 2
LA		The intersection between the centerline of the left clavicle and Rib 2
LL		Left part of the lower abdomen
RL		Right part of the lower abdomen
C(V) Chest electrode	C1(V1)	Between the right edge of the breast bone and Rib 4
	C2(V2)	Between the left edge of the breast bone and Rib 4
	C3(V3)	Halfway between Points V2 and V4
	C4(V4)	Intersection between the centerline of the clavicle and Rib 5
	C5(V5)	Horizontal intersection between the front line of the left armpit and Point V4
	C6(V6)	Horizontal intersection between the center line of the left armpit and Point V4

Note: Lead V or C (see V(C) chest electrodes above) may be located at any of the six positions as shown in the table above and the figure below (Points V1, V2, V3, V4, V5 and V6), in other words, the chest electrodes (Lead V or C) should be located at only one of the six positions to be measured. Only measure of the position where such leads are located will be taken.

 It might not display ECG wave with 3 leads. The 5 leads should be used to have ECG wave.

Safety Instructions for ECG Monitoring

-  Obstetric Monitor can only be equipped with ECG leads provided by our company; using ECG leads supplied by other companies may cause improper performance or poor protection while using defibrillator.
-  Electric parts of electrodes, leads and cable are forbidden to contact any other electric parts (including ground).

- 🔔 Obstetric Monitor can resist against defibrillator and electrosurgical unit. Readings may be inaccurate for a short time after or during using defibrillator or electrosurgical unit.
- 🔔 Transient caused by cable circuitry blocks while monitoring may be similar to the real heartbeat waveform, as a result resistance heart rate alarm rings. If you put the electrodes and cable in proper places according to this manual's instructions and the instructions for using electrode, the chance of this transient occurring will be decreased.
- 🔔 To the patient with pacemaker, due to that this device has been designed to provide resistance to pacemaker signal interference, generally the pacemaker pulse is not counted in heart rate measurement and calculation, but when the cycle time of pacemaker pulse is over 2ms, it may be counted. In order to reduce this possibility, observe the ECG waveforms on the screen carefully and do NOT rely entirely on the heart rate display and alarm system of this monitor when monitoring this kind of patients.
- 🔔 Besides the improper connection with electrosurgical unit may cause burns, the monitor may be damaged or arouse deviations of measurement. You can take some steps to avoid this situation, such as do NOT use small ECG electrodes, choosing the position which is far away from the estimated Hertzian waves route, using larger electrosurgical return electrodes and connecting with the patient properly.
- 🔔 ECG leads may be damaged while using defibrillator. If the leads are used again, please do the functional check first.

3.3.5 Blood Pressure Cuff Connection

1. Connect the cable to the right-panel connector marked with the NIBP icon.
2. Select a cuff of appropriate size according to the age of the subject. Its width should be 2/3 of the length of the upper arm. The cuff inflation part should be long enough to permit wrapping 50-80% of the limb concerned. See the table below for the dimensions:

Note: The size of the cuff selected should suit the subjects while measuring.

Cuff Model	Arm Circumference	Cuff Width
Neonate Cuff	6.0cm~9.5cm	3cm
Small-sized Infant Cuff	6cm~11cm	4.5cm
Middle-sized Infant Cuff	10cm~19cm	8cm
Large-sized Infant Cuff	18cm~26cm	10.6cm
Adult Cuff	25cm~35cm	14cm

- 🔔 Prior to use of the cuff, empty the cuff until there is no residual air inside it to ensure accurate measurement.
- 🔔 When putting on the cuff, unveil and wrap it around the upper arm evenly to appropriate tightness as shown below.


Figure 3.8 Cuff Position

Safety Instructions for NIBP Monitoring

- 🔔 When taking the measure of an infant or neonate's (less than 10 years old) blood pressure, do NOT operate in the adult mode. The high inflation pressure may cause lesion or even body putrescence.

- 🔊 It is recommended to take the blood pressure measurement manually. Automatic or continuous measurement should be used at the presence of a doctor/nurse.
- 🔊 NIBP monitoring is prohibited to those who have severe hemorrhagic tendency or with sickle cell disease, or partial bleeding will appear.
- 🔊 Do NOT bind NIBP cuff on limbs with transfusion tube or intubations or skin lesion area, otherwise, damages may be caused to the limbs.
- 🔊 Pay attention to the color and sensitivity of the limb when measuring NIBP; make sure the blood circulation is not blocked. If blocked, the limb will discolor, please stop measuring or remove the cuff to other positions. Doctor should examine this timely.
- 🔊 When unplugging the cuff, hold the head of the connector and pull it out.
- 🔊 Do NOT twist the cuff tube or put heavy things on it.

3.3.6 To connect the SpO₂

SpO₂ probe is very delicate equipment. Please follow the steps and procedures in operating it. Failure to operate it correctly can cause damage to the SpO₂ probe.

Operation procedure:

1. Connect the SpO₂ probe to the right panel's jack labeled "SpO₂". When unplugging the probe, be sure to hold the head of the connector and pull it out.
2. Insert one finger into the probe (index finger, middle finger or ring finger with proper nail length) according to the finger mark on the probe, shown as below.


Figure 3.9 Demonstration for SpO₂ probe

3. If the neonate SpO₂ sensor is used, please follow Figure 3.10 to connect.


Figure 3.10 Neonate SpO₂ sensor placement

Safety Introductions for SpO₂ Monitoring

- 🚫 Continuous use of fingertip SpO₂ sensor may result in discomfort or pain, especially for those patients with microcirculatory problem. It is recommended that the sensor should NOT be applied to the same finger for over two hours.
- 🚫 SpO₂ measuring position must be examined more carefully for some special patient. Do NOT install the SpO₂ sensor on the finger with edema or fragile tissue.
- 🚫 Do NOT put the SpO₂ sensor and pressure cuff on the same limb, otherwise the NIBP measuring will affect SpO₂ measuring and cause the alarm error.

- 🔔 Do NOT use the damaged SpO₂ sensor.
- 🔔 Please do not allow the cable to be twisted or bended.
- 🔔 Please do not use nail polisher or other cosmetic product on the nail.
- 🔔 The fingernail should be of normal length.
- 🔔 The SpO₂ sensor can not be immersed into water, liquor or cleanser completely, because the sensor has no capability of waterproofness.

3.3.7 TEMP Transducer Connection

Obstetric Monitor has two TEMP transducers to measure different body temperature.

Connecting methods:

1. Attach the transducer to the patient firmly;
2. Connect it to the connector labeled "TEMP" on the left panel.

Note: When unplugging the probe, be sure to hold the head of the connector and pull it out.

3.3.8 Loading Printing Paper and Hanging Installation

Operation Procedure for Loading Printing Paper:

1. Push both tabs of the printer shield toward middle direction with two hands, meanwhile, and pull it to open the shield.
2. Put the prepared folded printing paper into the printing paper box flatly towards right side, pull a part of the printing paper out of the shield and make the back of printing paper faced to the rubber roller.
3. Push the shield back in position and secure it until the sound "click" is heard.

Note: When placing the folded printing paper, please let the page with grid upward and the green strip near the left side.

Procedure for Hanging Monitor

Hanging function is optional. Please use the provided screw to fix the hanging panel in an appropriate position. Then hang the monitor by hooking the hanging hole on the bottom panel of the monitor. When hanging the monitor stably, you can loose your hands.

Please make sure that the hanging panel is fixed firmly and the monitor is hung stably, or the monitor may fall off causing personal injury or monitor damage.

Chapter 4 MONITORING SCREEN

4.1 Preparation Prior to Use

Before turning on the power, connect the probe to the corresponding outlet, and load the recording paper. Finally, connect the power cable of the apparatus to 100~250V mains outlet. Then turn on the monitor, the power indicator will light.

4.2 Initial Display Screen

Once the system is powered on, the initial monitoring screen will appear as shown in the figure below. In normal condition, this screen is the most commonly used one in monitoring of patients.


Figure 4.1 Initial Screen

On-screen display

- 1. **FHR monitoring:** the current FHR curve displays in this area.
- 2. **UC monitoring:** UC curve displays in this area and its value is shown on the right.
- 3. **ECG monitoring:** the real-time ECG of the patient being monitored and the information of ECG lead, ECG wave gain and ECG scanning speed display in this area. As shown in Figure 4.1, the current lead selection is "II", ECG wave gain is "X1" and ECG scanning speed is "25mm/s".

" X1 25mm/s "

X1- ECG wave gain (amplification fold):

"Auto" for automatic wave pattern gain "X1/2" for half reduction of the wave pattern.

"X1" for no change of wave pattern "X2" for one time amplification of the wave pattern.

- ECG lead marks: I, II, III, aVR, aVL, aVF, and V (V1~V6) available. When in the mother's screen, press the "Lead" key to switch among them.

ECG scanning speed: three options: 12.5, 25 and 50mm/second.

- 4. **SpO₂ Plethysmograph**
- 5. **Respiration waveform**
- 6. **System Information**

"  ": FHR volume icon. Press "+" / "-" to increase/decrease the volume.

: alarm silence icon. Pressing “SILENCE” key can make the monitor stay in alarm status; “ ” alarm icon, indicate when alarm event occurs, the monitor will alarm and the alarm indicator will be on.

“BED 03”: Network No.

- **7. Fetal Parameters**

“FHR1 124” - Fetal heart rate.

“TOCO 35” - uterine contraction pressure.

- **8. NIBP Parameters**

“mmHg”: the current blood pressure unit

“115” - systolic pressure value.

“66” - diastolic pressure value.

“86” - mean arterial pressure value.

“Adult” - indicates the current status information of the non-invasive blood pressure module. It varies with the mode of the non-invasive blood pressure measurement.

“Manu”-blood pressure measuring mode is Manual. Two options: “Manual” or “Auto”;

“_ _min” - indicates the time cycle of the auto NIBP measurement. When blood pressure measuring mode is “Auto”, it is valid.

- **9. ECG Parameters**

“HR 61” - the current heart rate value is 61bpm

“♥” -cardiac impulse sync marker

- **10. Temperature Parameters**

“TEMP 36.8” -body temperature value

“□”-body temperature unit

- **11. Respiration Parameters**

“Resp 14”-the current respiration rate

“BPM”- respiration rate UNIT

- **12. SpO₂ Parameters**

“SpO₂% 91”- pulse oxygen saturation value

“ PR 62”- pulse rate value measured by SpO₂ probe;

 -strength graph of SpO₂

- **13. System Menu**: after pressing or rotating the Navigation Knob, the system menu will be activated. Rotate the knob to shift the blue selecting box and then press the knob for entering the corresponding setup window. Three options: “System”, “Colors” and “Alarms”.

- **14. “2006-05-25 16:31:26”**-the current date and time. Date format: yyyy-mm-dd; time format: hh:mm:ss.

- **FM mark “y”**: during FHR monitoring, press FM maker to mark the FM of this time.

- **Event mark “z”**: during FHR monitoring, if important event or special event occurs, doctor or nurse can press “Event Marker” to mark the event time.

Operating instructions

- **“+” key**: press it to increase the volume;

- **“-” key**: press it to decrease the volume;

- **“PRINT” key**: print data;

- **“ZERO” key**: press this key, the TOCO will resume to initial pressure. During monitoring, press “ZERO” key to zero the TOCO.

- **“SILENCE” key**: press it for suspend/active the alarm sound.

- **“FREEZE” key**: press it for freeze the waveforms; press it again for the real-time waveforms.

- **“NIBP” key**: start/stop NIBP measurement. (this function is optional)

4.3 Freeze and Replay

During FHR monitoring, press “Freeze” key for freezing the monitoring screen, meanwhile “Freeze & Replay” will be displays in the lower right area above the system menu. The display screen is as shown in Figure 4.2.


Figure 4.2 Freeze and Replay Screen

Differences between Freeze and Replay Screen and Initial Screen

- The waveforms on the current screen are frozen;
- On Freeze and Replay Screen, you can review the historical data. Leftward/Rightward rotate the knob to scroll the page left/right for reviewing the trend graph and event mark in monitoring period;
- On Freeze and Replay Screen, only the current waveforms are frozen, the real-time waveform and data are still undergoing. The latest monitoring data can be reviewed by rightward rotating the knob to last page;

Press “FREEZE” key again to exit the Freeze and Replay Screen.


4.4 System Menu

On System Setup, Color Setup or Alarm Setup screen, please follow the operating instructions below to set the parameter values:

- Step1: Move the triangle cursor by rotating the knob to choose the item.
- Step2: Press the knob to confirm the selection, at this time the triangle cursor turns red;
- Step3: Rotate the knob to adjust the value or change the setting;
- Step4: Press the knob again to finish setting.
- Step5: Move the triangle cursor to “Exit”, press the knob for exiting from the setup window.

4.4.1 System Settings

On the Initial Screen, rotate the knob and let the selecting box stay in system option, then press the knob for entering into System Setup window, as shown in Figure 4.3.

System Setup Exit

No. [06] Date year [2008] Mon [03] Day [27] Time hour [10] Min [36] ◀

Language [English] Work mode [Demo] Display [Wave]

FHR Volume [5] FHR Speed [X1]

ECG Speed [25] Mm/s Lead [II] Gain [X1] Filter [Diag.]

Temp Unit [C]

NIBP Type [Adult] Node [Manual] Cycle [0] Min Unit [mmHg]

Resp Ctrl [OFF] Print Mode [1m/min]

Figure 4.3 System Setup

On-screen Display:

- **No.:** Network number;
- **Date year Month Day:** set the current date;
- **Time hour min:** set the current time;
- **Language:** operating system language.
- **Work mode:** set the monitoring mode. Two options: “Real” and “Demo”. If the selection is “Demo”, all the monitoring signals will be generated by the monitor.
- **Display:** for shifting display screen. “Wave”: ECG, SpO₂ and respiration waveforms will be displayed in the lower area of the main screen. If this option is changed to be “Data”, the waveforms displayed in the lower area of the main screen will be replayed by the NIBP list.
- **FHR Volume:** adjust the volume of FHR, 0~7 adjustable;
- **FHR Speed:** scan speed: “X1”, “X2”, “X4”; “X1” means FHR speed is 10mm/min; “X2” means FHR speed is 20mm/min; “X4” means FHR speed is 30mm/min;
- **ECG Speed:** the display speed of ECG curve. The options of 12.5, 25 and 50mm/second are provided. The default setting when powered on is 25 mm/second.
- **Lead:** seven options: I, II, III, aVR, aVL, aVF and V.
- **Gain:** ECG curve gain. The options of x0.5, x1, x2 and auto are provided. The default setting when powered on is x1.
 “x0.5”: half size of the basic waveform
 “x1”: the basic waveform
 “x2”: twice size of the basic waveform
- **Filter:** ECG filter modes. The options of monitoring, surgery and diagnosis are provided.
- **TEMP Unit:** body temperature unit. The default setting when powered on is “°C”. It may be changed to “°F” (1°F = 1°C*9/5+32).
- **NIBP Type:** When “Adult” is selected, the non-invasive blood pressure will be measured in the adult mode. When “Child” is selected, the NIBP measurement will be taken in the child mode. When “Neonate” is selected, it means the NIBP measurement will be taken in the neonate mode.
- **NIBP Mode:** “manual” and “auto” optional. The default setting when powered on is “Manual”. When “Manual” is selected, the blood pressure will be measured in a manual mode, in which case, whenever the “NIBP” key is pressed, the blood pressure measure will be taken. To stop or suspend the ongoing measurement, press this key another time. When “Auto” is selected, non-invasive blood pressure will be automatically measured.
- **NIBP Cycle:** When the “Auto” mode is selected, the operator may set the cycle value.
- **NIBP Unit:** blood pressure unit. The default setting when powered on is “mmHg”. Alternatively, it may be set to “kPa” (1kPa=7.5mmHg).
- **Resp Control:** turn “on” or turn “off” the respiration monitoring function;
- **Print Mode:** select the print speed: “1cm/min, 2 cm/min, 3 cm/min” for printing FHR curve and TOCO curve, or select “ECG” for printing ECG, plethysmogram and respiration waveform.

4.4.2 Color Settings

On the Initial Screen, rotate the knob and let the selecting box stay in color option, then press the knob for entering into System Setup window, as shown in Figure 4.4. Adjust the color of each parameter on this screen.


Figure 4.4 Color Setup

4.4.3 Alarm Settings

On the Initial Screen, rotate the knob and let the selecting box stay in alarms option, then press the knob for entering into System Setup window, as shown in Figure 4.5.


Figure 4.5 Alarm Setup

On-screen Display:

- "ON"/ "OFF": turn the Alarm on/off.
- Alarm limits of parameters

Items	High Limit	Low limit
FHR	160~190 bpm	90~120 bpm
HR	32~250 bpm	30~248 bpm
RR	1~99bpm	0~98bpm
Temperature	34~45°C	26~40°C
SpO ₂	71~99%	70~98%
Systolic Pressure	240 mmHg	60mmHg
Mean Arterial Pressure	200 mmHg	40 mmHg
Diastolic Pressure	180mmHg	30mmHg

4.4.4 NIBP List

On System Setup screen, set the item “Display” as “Data”, and then exit the System Setup screen. The NIBP list will be displayed in the lower area of the main screen, as shown in Figure 4.6. The latest 24 groups of NIBP values can be displayed here, if the record is more than 24 groups, rotate the knob to view the earlier records. Press the knob to exit from the NIBP list screen.

Time	Sys	Dias	Mean	PR	FHR 1/2	TOCO	HR	Temp	RR	SpO2
15:06	122	64	92	77	156/130	35	60	36.4	17	96
15:16	112	60	94	77	156/130	35	60	36.4	17	96
15:26	122	64	92	77	156/130	35	60	36.4	17	96
15:36	122	64	92	77	156/130	35	60	36.4	17	96
15:46	122	64	92	77	156/130	35	60	36.4	17	96

Figure 4.6 NIBP List

Chapter 5 TECHNICAL SPECIFICATIONS

5.1 Fetal Monitoring

A: Ultrasonic probe

Ultrasonic working frequency: 1MHz tolerance: $\pm 5\%$ ISATA: $\leq 10\text{mW/cm}^2$
 Fetal heart rate measuring range: 65~210bpm

Accuracy of fetal heart rate: $\pm 2\text{bpm}$

FHR alarm: when the FHR value exceeds the preset value, the monitor will alarm with 30s and pressing Alarm key will stop the current alarm.

B: TOCO Transducer

Zero setting mode: Manual

TOCO measuring range: 0%~100% tolerance: $\pm 10\%$

C: Marker Mode

Mode: Manual and AUTO

5.2 ECG/RESP Monitoring

ECG amplification bandwidth:	Monitoring mode: 1 Hz~25Hz (+0.4dB -3.0 dB)	Diagnostic mode: 1 Hz~75Hz (+0.4dB -3.0 dB)
ECG measuring range:	20~250bpm (Accuracy: $\pm 2\%$ or ± 2 bpm, whichever is greater.)	
Common-mode rejection ratio:	CMRR ≥ 80 dB	
ECG noise level:	$\leq 30\mu$ VP-P	
ECG input loop current:	$\leq 0.1\mu$ A	
Input impedance:	≥ 5 M Ω	
Time constant	Monitoring mode: ≥ 0.3 s	Diagnostic mode: ≥ 3.2 s
ECG gain:	X1/2	5mm/mV $\pm 5\%$
	X1	10mm/mV $\pm 5\%$
	X2	20mm/mV $\pm 5\%$
RESP monitoring range	0~100rpm (Accuracy: $\pm 5\%$ or ± 2 rpm, whichever is greater.)	

5.3 NIBP Monitoring

Measuring technology:	Automatic measurement by smart oscillating method	
Initial cuff inflation pressure:	Adult mode	180mmHg
	Infant mode	120mmHg
	Neonate mode	90mmHg
Max. inflation pressure allowed:	Adult: 300mmHg Infant/Neonate: 140mmHg Accuracy: $\pm 10\%$ or ± 10 mmHg, whichever is greater	
Factory-set measuring mode:	Manual, auto mode	
Measuring intervals:	1-240 minutes adjustable, interval tolerance: < 10 s	
Monitoring range:	Systolic pressure: 50~240mmHg	
	Mean arterial pressure: 25~200mmHg	
	Diastolic pressure: 15~180mmHg	
Measurement accuracy	± 8 mmHg	

5.4 SpO₂/Temperature Monitoring

SpO ₂ sensor:	Dual-wavelength LED	
Pulse oxygen saturation range:	0~99%	
Accuracy of pulse oxygen saturation:	75~99%	$\pm 2\%$
	50~74%	$\pm 3\%$
	Lower than 50% will not be defined.	
Pulse range:	30 ~ 240 times/minute	
Accuracy:	$\pm 2\%$ or ± 2 bpm, whichever is greater.	
Temperature measuring range:	28.0~45.0°C	
Accuracy	± 0.2 °C	
Temperature response time:	≤ 150 s	

5.5 Other Technical Information and Parameters

Mains power supply:	100~250VAC, 50Hz/60Hz
Power consumption:	<100VA
Internal power supply:	Li-ion, rechargeable 4 hours to full recharge during monitoring 6 hours to full recharge when monitor is off
Display mode:	TFT display
Alarm mode:	Audible & visible alarm
Networking interface:	Net port
Printing paper	Z-Fold, 112mm×90mm
The type of protection against electric shock:	Class I
The degree of protection against electric shock:	With Type BF, CF applied part
The degree of protection against harmful ingress of liquids:	Ordinary equipment without protection against ingress of water

Chapter 6

DIMENSIONS AND WEIGHT

6.1 Overall dimensions and weight

The product is packed in high quality corrugated cartons with foam inside to protect the apparatus against damage in the handling process.

Gross weight: Details see the indication on the outer package

Dimension: 360L×230W×120 H (mm)

6.2 Accessories

1	Ultrasonic probe	One set
2	Uterine contraction probe	One set
3	Fetal movement probe	One piece
4	Abdominal belt	Two pieces
5	Ultrasonic coupling gel	One bottle
6	ECG cable	One set
7	Blood pressure cuff	One piece
8	Blood oxygen probe	One piece
9	Body temperature probe	One piece
10	Power cable	One piece
11	Grounding wire	One piece
12	Disposable electrode	Ten pieces
13	User Manual	One copy
14	Warranty	One copy
15	Quality Certificate	One copy

Note: The accessories provided may vary with the model of the main equipment.

Chapter 7

OPERATING PRINCIPLE OF THE MONITOR

7.1 The principle for measuring the fetal heart rate

This monitor tests the fetal heart rate through non-invasive ultrasonic Doppler Effect. As is known, ultrasonic wave propagating at a given frequency will be reflected when encountering an obstacle. If it is a still obstacle, the back wave will share the same frequency with the transmitted wave. Once the obstacle moves, the frequency of the back wave will be changed. The higher rate the object moves at, the bigger frequency change will take place. This is the so-called Doppler Effect. With the apparatus, the ultrasonic probe is placed on the abdominal surface of the pregnant woman. The fetal heartbeat is motion with respect to the ultrasonic wave probe. When the transmitted wave encounters the fetal heart, the back wave will develop offset frequency. With the offset frequency, we can work out the fetal heart rate and frequency.


7.2 The principle for measuring the uterine pressure and fetal movement count

Measurement of the uterine pressure is, in practice, external measurement of pressure caused by uterine contraction by means of a pressure-sensing device (pressure probe), which is to be placed on the abdomen of the pregnant woman. It is used to assess the count, duration and relative strength of the uterine contraction. The index measures are relative values with respect to the indices while in the equipoise (the uterine pressure value is shown in percentage).

By fetal movement counting, while in the real-time monitoring process, whenever the pregnant woman feels fetal movement, a fetal movement mark may be left on the uterine pressure wave pattern automatically or manually by pressing the fetal movement probe.

7.3 How to Obtain High Quality ECG and Accurate Heart Rate Value

The electrocardiogram (ECG or EKG) is primarily a tool for evaluating the electrical events within the heart. The action potentials of cardiac-muscle cells can be viewed as batteries that cause charge to move throughout the body fluids. These currents represent the sum of the action potentials occurring simultaneously in many individual cells and can be detected by recording electrodes at the surface of the skin. The figure below shows the system of the heart.

First of all, the hospital should be equipped with a 100~250V power supply system with a typical grounding wire. If big interference in ECG continues, connect one end of the grounding wire provided with this equipment to the grounding wire on the back panel of this monitor, and the other end to the special grounding wire, water pipe or radiator.

A common ECG plate electrode used together with this monitor has short shelf life. Generally, the shelf life is only one month after the package is opened. When outdated plate electrode is used, due to skin's contact impedance and big electrode potential, the chance of interference will be increased, and the ECG baseline will have an unstable inclination. Therefore, always use valid plate electrodes.

7.4 NIBP Monitoring

7.4.1 The Principle for Measurement of Blood Pressure

Blood pressure may be measured in an invasive way (whereby the sensor will be inserted into blood vessel directly) or a non-invasive way. The non-invasive way includes several methodologies, such as the Korotkoff Sound Method and oscillating method. The Korotkoff Sound Method is used as a conventional way, whereby stethoscope is used to measure the blood pressure. By the oscillating method, an inflation pump will fill the air, and release it slowly. A computer will record change of the cuff pressure when the air is released. With this record, the blood pressure value will be determined. First of all, make sure the signal quality judgment by computer meets the requirements of accurate calculation (such as sudden limb movement or cuff being hit during the measurement). If the answer is negative, give up the calculation. If the answer is positive, proceed with calculation of the blood pressure value.

As change of the blood pressure is recorded by electric sensor, which sensitivity is much higher than that of human ears, the oscillating method uses different definitions for measurement of diastolic pressure, mean arterial pressure and systolic pressure from the Korotkoff Sound Method. When the oscillating method is used, the circuit in the measuring apparatus will separate the amplitude of the cuff pressure from its change with pulsation. With the oscillating method, the blood pressure at the maximum amplitude of cuff pressure is defined as the mean arterial pressure. The blood pressure at amplitude of cuff pressure forward reduced according to proper proportion is defined as systolic pressure, while the blood pressure at amplitude of cuff pressure backward reduced according to proper proportion is defined as diastolic pressure. The maximum change of pulse pressure occurs at these two points. They are equivalent to the point with pulse sound and the point without pulse sound respectively in the Korotkoff Sound Method.

When the risk of invasive monitoring method outweighs its advantage of accuracy, non-invasive monitoring method shall be used.

7.4.2 Points to be Noted in Measurement of Blood Pressure

Like common non-invasive blood pressure measurement, improper operation may cause inaccurate or blank result or misunderstanding of the measuring information when the oscillating method is used to take the measure of blood pressure. This point needs particular attention of the operators.

1. Requirements of the cuff:
 - 1) Appropriate cuff should be selected according to the age of the subject.
 - 2) Remember to empty the residual air in the cuff before the measurement is commenced.
 - 3) Locate the cuff in such a way that the "☒" mark is at a location where the clearest pulsation of brachial artery is observed.
 - 4) The cuff should be tightened to a degree where insertion of one finger is allowed.
 - 5) The lower end of the cuff should be 2cm above the elbow joint.
2. The subject should lie on the back so that the cuff and the heart are in a horizontal position and the most accurate measure is taken. Other postures may lead to inaccurate measurement.
3. Do not speak or move before or during the measurement. Care should be taken so that the cuff will not be hit or touched by other objects.
4. The measures should be taken at appropriate intervals. Continuous measurement at too short intervals may lead to pressed arm, reduced blood flow and lower blood pressure, and resulting inaccurate measure of blood pressure. It is recommended the measure be taken at intervals of more than two minutes.
5. With the oscillating method, when blood pressure is measured, the inflation pressure of the cuff will be automatically adjusted according to the previous measure. Generally, the initial inflation pressure is 180mmHg (for the adult mode) or 100mmHg (for the infant mode) or 80 mmHg (for the neonate mode) when it is powered on. Following that, 50mmHg (for the adult mode) or 30mmHg (for infant mode) or

10mmHg (for the neonate mode) will be added on the basis of the last measurement of systolic pressure. In this way, when the blood pressure rises or the subject is changed, the blood pressure meter may fail in giving the result after the first-time inflation. This monitor will automatically adjust the inflation pressure until the measure is taken, after that, up to four measures will be allowed.

6. When an adult subject is monitored, the machine may fail in giving the blood pressure measure if the infant or neonate mode is selected.

7.4.3 Contraindication

1. Serious angiospasm, vasoconstriction, or too weak pulse.
2. When extremely low or high heart rate or serious arrhythmia of the subject occurs. Especially auricular fibrillation will lead to unreliable or impossible measurement.
3. Do not take the measurement when the subject is connected with an artificial heart-lung machine.
4. Do not take the measurement when the subject uses diuresis or vasodilator.
5. When the subject is suffering from major hemorrhage, hypovolemic shock and other conditions with rapid blood pressure change or when the subject has too low body temperature, the reading will not be reliable, for reduced peripheral blood flow will lead to reduced arterial pulsation.
6. Subject with hyperadiposis;

In addition, statistics show that 37% people report blood pressure difference of no less than 0.80kPa(6mmHg) between the left and right arms, and 13% people report difference of no less than 1.47kPa (11mmHg).

Note: Some practitioners may report big discreteness or abnormal value of the blood pressure measures when the oscillating method is used. As a matter of fact, the so-called “big discreteness” must be a term in the sense of statistical significance of mass data. Abnormal data may be observed in some individual cases. It is normal in the scientific experiments. It may be caused by an apparent reason, or by an unknown factor in some cases. Such individual doubtful experimental data may be identified and eliminated using the special statistical technique. It is not a part of this manual. The practitioner may eliminate the apparently unreasonable data according to the experience.

7.5 SpO₂/Pulse Monitoring

7.5.1 SpO₂ Monitoring

This monitor measures the pulse oxygen saturation (SpO₂) and pulse by means of the radiograph of infrared light and the red light emitted by LED through body's peripheral areas (such as fingers), whereby the photoelectric detecting circuits will analyze the absorptivity of the oxyhemoglobin and reduced hemoglobin respectively, and give the photoabsorption rates before and after pulsation. Using the measure of photoabsorption change due to pulsatory arterial blood flow caused by PLETH waveform, the SpO₂ can be obtained.

7.5.2 Points to be Noted in SpO₂ and Pulse Measuring

1. The fingers should be properly placed (see the attached illustration of this instruction manual), or else it may cause inaccurate measurement result.
2. Make sure that capillary arterial vessel beneath the finger is penetrated through by red and infrared lights.
3. The SpO₂ sensor should not be used at a location or limb tied with arterial or blood pressure cuff or receiving intravenous injection.
4. Do not fix the SpO₂ sensor with adhesive tape, or else it may result in venous pulsation and consequential inaccurate measurement result of SpO₂.
5. Make sure the optical path is free from any optical obstacles like adhesive tape.
6. Excessive ambient light may affect the measuring result. It includes fluorescent lamp, dual ruby light, infrared heater, and direct sunlight etc.
7. Strenuous action of the subject or extreme electrosurgical interference may also affect the accuracy.
8. Please do not use the SpO₂ sensor when having the MRI, or burn may be caused by faradism.

Clinical Limit

1. As the measure is taken on the basis of arteriole pulse, substantial pulsating blood stream of subject is required. For a subject with weak pulse due to shock, low ambient/body temperature, major bleeding, or use of vascular contracting drug, the SpO₂ waveform (PLETH) will decrease. In this case, the measurement will be more sensitive to interference.

2. For those with a substantial amount of staining dilution drug (such as methylene blue, indigo green and acid indigo blue), or carbon monoxide hemoglobin (COHb), or methionine (Me+Hb) or thiosalicylic hemoglobin, and some with icterus problem, the SpO₂ determination by this monitor may be inaccurate.
3. The drugs such as dopamine, procaine, prilocaine, lidocaine and butacaine may also be a major factor blamed for serious error of SpO₂ measurements.
4. As the SpO₂ value serves as a reference value for judgement of anemic anoxia and toxic anoxia, the measurement result of some patients with serious anemia may also present as good SpO₂ value.

Chapter 8

TROUBLESHOOTING

Note: In case of trouble of this machine in service, follow the instructions below to eliminate the problem first. If the attempt fails, contact the dealer in your local area or the manufacturer.

Do NOT open the monitor without permission

8.1 No Display on the Screen

Shut down the machine and unplug the power cable. Use a universal meter to check if the outlet has proper voltage, if the power cable is in good condition, and if the power cable is properly connected with this apparatus or outlet. Remove the fuse from the back cover of this machine, and make sure it is in good condition.

8.2 Excessive ECG Signal Interference or too Thick Baseline

1. Check if the plate electrodes are properly located, and if valid plate electrodes are used.
2. Check whether the lead wires are properly inserted. If no ECG curve displayed, check if the ECG lead wires are broken.
3. Make sure the mains outlet has standard grounding wire.
4. Check if the grounding wire of the apparatus properly grounded.

8.3 No Blood Pressure and Pulse Oxygen Measures

1. Check if the blood pressure cuff is properly wrapped around the arm according to the operating instructions, if the cuff leaks, and if the inlet is closely connected with the NIBP jack on the side panel. Check if the indicator of the pulse oxygen sensor flashes and if the pulse oxygen probe is properly connected to the SpO₂ jack on the side panel.
2. If the problems still exist, please contact the manufacturer.

8.4 Abnormal FHR Detection

1. When neither fetal heart beat sound nor FHR curve can be detected, please check the transducer first to see whether it is placed correctly.
2. When the fetal heart beat sound can be detected but the FHR curve is disordered or in-and-out, please adjust the transducer for it is likely that it is placed at the abdominal area of the fetus.
3. The FM (Fetal Movement) or FHR curve may become abnormal after the mother changes her body position due to the movement of the fetal heart that makes the transducer deviate from the heart.
4. When the fetus moves downwards with the constriction of the uterus, the TOCO record may get disordered or interrupted, which can be avoided by fixing the transducer 1cm-2cm below the normal position.
5. The coupling gel lessens if the examination lasts for a long time so that the transducer may not work normally. Please supplement the gel in time.
6. Difficulties may be caused by the fetus occipital-posterior position that makes the fetal back moves towards that of the mother's, for the fetus sometimes may descend and turn. It is impossible to move the ultrasonic transducer to the mother's back, so it is advisable to move it to the middle of the abdomen below the umbilicus.
7. The FHR curve detected may be desultory during monitoring or it may occur frequently, which means that the transducer is not placed in the optimum position.
8. The following two reasons may be responsible for the abnormally low or inaccurate FHR value after a period's monitoring. First, the ultrasonic transducer may deviate from the optimum position because

the mother moves her body; second, the inaccurate FHR value may result from the fetal position movement.

9. If the transducer is not placed in the proper position during monitoring, what detected by the monitor may be fetal pulse or umbilical cord blood flow. This would account for why there is only FHR but no audible rhythmical FHB (Fetal Heart Beat) sound. If the best position still can't be found out after close examination, it is necessary to perform other test to determine whether the fetus is all right.

8.5 Abnormal Uterine Contraction Pressure Detection

1. When Uterine Contraction Pressure is normal but the monitor fails to display the Uterine Contraction Pressure curve, please check whether the transducer is placed correctly.
2. The Uterine Contraction Pressure transducer that is placed at the uterine fundus area may not able to record the changes of the Uterine Contraction Pressure because it descends with the progresses of labor. So it is advisable to move the transducer downwards, fix it and perform zero-setting.
3. The abdominal belt may get loose easily after the pregnant woman turns over, thus the Uterine Contraction Pressure curve can't be recorded clearly. Please fasten the belt again when necessary.
4. The ultrasonic coupling gel that may cause damage to the Uterine Contraction Pressure transducer and detection failure can't be applied to the Uterine Contraction Pressure transducer. Please clear it thoroughly when necessary.

Chapter 9

MAINTENANCE AND SERVICE

The monitor should be properly maintained to ensure its maximum performance and long service life. In addition to the warranty period, the company also offers long-term service for each customer. It is important that the users read and follow the operating instructions, important information and maintenance measures.

9.1 Technical Maintenances

Before using the monitor, the checks below should be carried out:

- ☞ Check the monitor for any mechanical damage;
- ☞ Inspect the exposed parts and the inserted parts of all the leads and the accessories;
- ☞ Examine all the functions of the monitor that are likely to be used for patient monitoring, and ensure that it is in good working condition;
- ☞ Make sure that the monitor is grounded properly.
- ☞ Pay close attention to the fluctuation of the local power supply voltage. A manostat is recommended when necessary.

In case any indication of damage about the function of the monitor is detected and proven, it is not allowed to apply it to the patient for any monitoring. Please contact the local dealer or our company, and we are to offer the best solution as soon as possible for your satisfaction.

After each maintenance or the yearly maintenance, the monitor must be thoroughly inspected by qualified personnel, including function and safety examinations.

- ⚠ If the hospital fails to carry out a satisfactory maintenance program about the monitor, it may get disabled and endangered the patient's safety and health.
- ⚠ If there is any indication of cable damage or deterioration, they are prohibited from any further use.
- 🔔 The adjustable units in the monitor such as potentiometer are not allowed to adjust without permission to avoid unnecessary failures that affect normal application.

9.2 User Maintenances

- 1 Never use this machine in an environment with inflammable gas.
- 2 Avoid electric shock. The power cable should be plugged into 100-250VAC outlet with grounding wire. Do not use an outlet in poor condition. If the condition allows, use a power supply system with regulators.
- 3 Do not unplug the power cable with the wet hand. Keep the hand clean and dry whenever touch the cable.
- 4 Keep it away from corrosive substances, explosive substances, high temperature and moisture.
- 5 If the monitor gets wet, please stop operating it. When it is transferred from cold environment to warm and humid environment, please not use it immediately.

- 6 Do not let the transducer immerse into the liquid. When perform disinfection to the transducer, please wipe its surface with disinfectant impregnated wipe. Do not spray any liquid on it directly.
- 7 When clean the transducer belt, the water temperature should not exceed 60℃.
- 8 High temperature steam disinfection of the accessories is not permissible.

9.3 Cleaning, Sterilization and Disinfection

- 🔌 Switch off monitor and disconnect the power cable before cleaning.
- 🗑️ The monitor should be kept from dust.
- 🗑️ It is recommended to clean the outer shell and screen of the monitor to keep it clean. Only non-corrosive cleanser such as clear water is permitted.
- 🗑️ Use soft cloth with alcohol to wipe the surface of the monitor and transducers, and dry it with dry and clean cloth or simply air-dry.
- 🗑️ The monitor can be sterilized or disinfected, please clean the monitor first.
- 🔌 Do not let the liquid cleanser flow into the connector jack of the monitor to avoid damage.
- 🔌 Clean the exterior of the connector only.
- 🔔 Dilute the cleanser.
- 🔔 Do not use scrub materials.
- 🔔 Do not let any liquid flow into the shell or any parts of the monitor.
- 🔔 Do not let the cleanser and disinfectant stay on its surface.
- 🔔 Do not let any parts of the monitor immerse into the liquid.
- 🔔 Do not pour the disinfectant on its surface while sterilization.

9.4 Cleaning, Sterilization and Disinfection of Accessories

It is recommended to use a cloth dampened with isopropyl alcohol 70%, a 10% aqueous solution of sodium hypochlorite (bleach), a 2% glutaraldehyde solution, ammonia, mild soap or disinfectant spray cleaner to clean the accessories (including sensor, cable and plugs).

- 🔌 Do not use damaged accessories.
- 🔌 Accessories can not be entirely immersed into water, liquor or cleanser.
- 🔌 Do not use radial, steam or epoxyethane to disinfect accessories.

9.5 Storage

If the apparatus will not be used for a long period of time, wipe it clean and keep it packed in a dry place, free from dust.

Storage environment:	ambient temperature:	-20~60℃
	relative humidity:	10%~95%
	atmospheric pressure:	50 kPa~107.4kPa

9.6 Transportation

This monitor should be transported by land (vehicle or railway) or air in accordance with the contractual terms. Do not hit or drop it with force.

Disposal: *The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment. For further information on recycling points contact the local authorities, the local recycling center or the shop where the product was purchased. If the equipment is not disposed of correctly, fines or penalties may be applied in accordance with the national legislation and regulations.*

GIMA WARRANTY CONDITIONS

Congratulations for purchasing a GIMA product. This product meets high qualitative standards both as regards the material and the production. The warranty is valid for 12 months from the date of supply of GIMA.

During the period of validity of the warranty, GIMA will repair and/or replace free of charge all the defected parts due to production reasons. Labor costs and personnel traveling expenses and packaging not included. All components subject to wear are not included in the warranty. The repair or replacement performed during the warranty period shall not extend the warranty. The warranty is void in the following cases: repairs performed by unauthorized personnel or with non-original spare parts, defects caused by negligence or incorrect use. GIMA cannot be held responsible for malfunctioning on electronic devices or software due to outside agents such as: voltage changes, electro-magnetic fields, radio interferences, etc. The warranty is void if the above regulations are not observed and if the serial code (if available) has been removed, cancelled or changed.

The defected products must be returned only to the dealer the product was purchased from. Products sent to GIMA will be rejected.