MASCHERA CON CUSCINETTO IN SILICONE
MASQUE AVEC COUSSINET DE SILICONE
MASKE MIT SILIKONKISSENDICHTUNG
MASCARILLA CON COJIN DE SILICONA
SILICONE CUSHION MASK

REF 34250 - 34252 - 34253 - 34254 - 34255

MANUALE D'USO E MANUTENZIONE
USE AND MAINTENANCE BOOK
INSTRUCTIONS DE FONCTIONNEMENT ET ENTRETIEN
BETRIEBS UND WARTUNGS ANWEISUNGEN
MANUAL DE USO Y MANTENIMIENTO

ATTENZIONE: Gli operatori devono leggere e capire completamente questo manuale prima di utilizzare il prodotto.
ATTENTION: The operators must carefully read and completely understand the present manual before using the product.
AVIS: Les opérateurs doivent lire et bien comprendre ce manuel avant d'utiliser le produit.
ACHTUNG: Die Bediener müssen vorher dieses Handbuch gelesen und Verstanden haben, bevor sie das Produkt benutzen.
ATENCIÓN: Los operadores tienen que leer y entender completamente este Manual antes de utilizar el producto.

M34250-M-Rev.2.1-09.14

CE 0476

Last cover page

Simbologia / Index of symbols

REF	Codice prodotto Product code Numero de reference Katalognummer Número de catálogo		Teme l'umidità Keep dry Garder au sec Trocken halten Conservar en lugar seco
LOT	Numero di lotto Lot number Numero de lot Chargenkennung Código del lote		Riparare da luce diretta Keep away from sunlight Eviter exposition au soleil Vor sonnenlicht schützen Mantener lejos de la luz solar
	Leggere le istruzioni Read instructions Voir les instructions Siehe Gebrauchsanweisung Consulte las instrucciones		Data di fabbricazione Manufacturing date Date de fabrication Datum der Herstellung Fecha de fabricación
	Fabbricante Manufacturer Fabricant Hersteller Fabricante		Non riutilizzare Do not reuse Usage unique Nicht zur wiederverwendung No reutilizar

 GIMA Spa - Via Marconi, 1
 20060 Gessate (MI) - Italia
 Made in Taiwan

Instruction for use

Silicone Cushion Mask

With Inflated cuff

1. Warning

U.S. federal law restricts this device to sale by or on the order of a Physician (USA only).

For use by CPR trained personnel only. The proper application of a face mask to obtain a tight seal should be trained in particular. Make sure that such personnel are made familiar with the contents of this manual.


2. Description

Silicone Cushion Mask are designed for use with manual resuscitators, automatic resuscitators and ventilators.

Size 2,3,4,5 and 6 with a 22mm ISO connector. They can be autoclavable at 134 degree centigrade repeatedly.

Transparent dome can observe patient's condition. The soft cuff inflation can be adjusted.

■ SILICONE-CUSHION MASK

3. Operating instructions:

Ventilation with a cushion mask. Apply the mask firmly to the patient's face to achieve a tight seal. Hold the mask tight against the face while maintaining the correct head tilt to keep the airway open.

4. Cleaning disinfecting sterilizing

Method		Applicable	
Cleaning	Washing	Manual washing	√
		Washing machine	√
	Disinfecting	W.M.heat disinfecting	√
		Boiling	√
		Chemical	√
	Autoclaving	121°C	√
134°C		√	

The table shows applicable methods for cleaning, disinfecting and sterilizing for mask after each use.

Detergents and chemical disinfectants


Numerous brands of detergents and Chemical disinfectants are available. Use only brands that are compatible with the mask materials to avoid reduction in the lifetime of the materials.

Note: If residuals of detergents or disinfectants are allowed to dry on the mask materials, due to surface may be affected and the product lifetime may become reduced.

4.1 Disassembly of the mask (recommended)

Unbutton or remove mask cuff from dome.

Caution! Before autoclaving remove stopper from the inflation tube.

4.2 Cleaning of parts

Manual washing or automatic washing machine

Wash the parts in warm water using a detergent suitable for the mask materials. Rinse all parts thoroughly in clean water to remove all detergent residues.

If surface cleaning and/or disinfecting of the outer surfaces of the mask are carried out, make sure that the detergent and/or disinfectant is compatible with the materials of the mask and make sure to remove the detergent/disinfectant by rinsing with water.

An automatic washing machine with a programme designed for washing anesthesia breathing equipment may be used.

4.3. Disinfecting and/or sterilizing

Select heated or chemical disinfection according to local standards for disinfection

and the table of applicable methods.

P4

Washing machine Heat disinfecting program

Automatic washing machines designed for medical equipment will normally include program cycles for heat disinfection.

Boiling

Use clean water, heat and boil parts for 10 minutes to disinfect.

Chemical disinfecting

Follow the instructions of the manufacturer of the chemical disinfectant for dilution and exposure time. After exposure of the mask parts to the chemical disinfectant, rinse thoroughly in clean water to remove all residues.

Autoclaving

Use standard autoclaving equipment adjusted for max. 134°C.

Caution! Before autoclaving remove stopper from the inflation tube.

4.4 Drying and cooling

Leave the parts to dry and/or cool completely before reassembling the mask.

4.5 Inspection of parts

After cleaning, disinfecting and/or sterilizing carefully inspect all parts for damage or excessive wear and replace, if necessary. Some methods may cause discoloration of rubber parts without impact on their lifetime. In case of materials deterioration e.g. cracking, the parts should be replaced.

4.6 Assembling

Reassemble the mask.

5. Specifications

Weight

Size 2: app. 40 g.

Size 3: app. 62 g.

Size 4: app. 67 g.

Size 5: app. 78 g.

Size 6: app. 92 g

P5

Operating temperature range: 40°C to 50°C (40°F to 122°F)

Storage temperature range: -40°C to 70°C (40°F to 158°F)

Parts Materials

Mask dome: Polysulphone or Polycarbonate

Mask cuff: Silicone rubber

Mask connector: Silicone rubber

Hook: Stainless

Order Information:

34252 Durable Silicone Cushion Mask #2: 22mm ISO connector

34253 Durable Silicone Cushion Mask #3: 22mm ISO connector

34254 Durable Silicone Cushion Mask #4: 22mm ISO connector

34255 Durable Silicone Cushion Mask #5: 22mm ISO connector

34250 Durable Silicone Cushion Mask #6: 22mm ISO connector

NOTE: Polycarbonate is not a perfect durable material

The product life will vary due to different sterilization method.

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