CE
0123

AEROSOL AIR THERAPY
AIR THERAPY NEBULIZER SYSTEM
AIR THERAPY COMPRESOR-NEBULIZADOR
AIR THERAPY INHALATOR
AIR THERAPY AÉROSOL

MANUALE D'USO
USER MANUAL
HANDBUCH
MANUAL DE INSTRUCCIONES
MODE D'EMPLOI

Moretti S.p.A.
Via Bruxelles 3 - Melegnano
52022 Cavriglia (Arezzo)
Tel. +39 055 96 21 11

Fax. +39 055 96 21 200
www.morettispa.com
info@morettispa.com

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AIR THERAPY is a piston-type compressor nebulizer system working at 230V/50Hz. High performance with any type of drug, ideal for intensive hospital and clinic use.

Manufactured with high thermal and electric insulation plastic chassis in compliance with the latest European Safety regulations.

The oil-free piston compressor has long durability and is equipped with the highly efficient HI-FLO jet nebulizer to guarantee quick and accurate drug delivery. The device is designed for easy transport and holding and is recommended for atomising antibiotics and bronchodilator drugs. The medical device is designed for continuous use.

GENERAL WARNING

READ INSTRUCTION MANUAL CAREFULLY BEFORE USE

DRUG ADMINISTRATION MUST BE UNDER MEDICAL CONTROL

THE INSTRUMENT MUST NOT BE DISASSEMBLED. FOR A TECHNICAL SERVICE ALWAYS CONTACT MORETTI SpA

IMPORTANT SAFETY RULES

- On opening the packaging, check the integrity of the appliance, paying particular attention to the presence of damage to the plastic parts, which may make access possible to internal live parts and also to breakage and / or peeling of the power supply cable. **In these cases don't connect the plug to the electric socket. Carry out these controls before each use;**
- before connecting the appliance always check that the electric data indicated on the data label and the type of plug used, correspond to those of the mains electricity to which it's to be connected;
- If the plug supplied with the appliance is incompatible with the mains electricity socket, contact qualified staff for replacement of the plug with a suitable type. The use of simple or multiple and / or extension adapters is not generally recommended. Whenever their use is indispensable, use those in compliance with safety regulations, however paying attention not to exceed the maximum power supply limits, which are indicated on the adapters and extensions;
- Never leave the appliance inserted if not necessary disconnect the plug from the mains power supply when it is not being used;
- Respect the safety regulations indicated for electrical appliances and particularly:
 - Only use original accessories and components provided by the manufacturer to guarantee the highest efficiency and safety of the device;
 - Never immerge the appliance into water;
 - Position the appliance on flat stable surfaces;
 - Position the device in a way that the air inlets on the back aren't obstructed;
 - Never use the device in environments which have anaesthetic mixtures inflammable with air, oxygen or nitric oxide;
 - Don't touch the device with wet hands and always prevent the appliance coming into contact with liquids;
 - The use of this device by children and / or incompetent person always requires the careful surveillance of an adult in possession of their full mental faculties;
 - The medical device, and most of all the nebulae, must be kept out of children's reach as it contains small parts that could be swallowed;
 - Don't leave the appliance connected to the power supply socket when not in use;
 - Don't pull the power supply cable to disconnect the plug remove the plug from the mains socket correctly;
 - Preserve and use the medical device in environments protected from atmospheric factors and at a distance from heat sources;
- For repairs, exclusively contact GIMA SpA technical service and request the use of original spare parts. Failure to comply with the above can jeopardise the safety of the device;
- This medical device must be destined exclusively for the use for which it has been designed and described in this manual. It must therefore be used as an aerosol therapy system.** Any different use must be considered incorrect and therefore dangerous; the manufacturer cannot be considered liable for damage caused by improper, incorrect and / or unreasonable use or if the appliance is used in electrical plants that are not in compliance with the regulations in force;
- Particular precautions must be made concerning electromagnetic compatibility. The medical device must be installed and used according to information supplied with the accompanying documents;
- Store the accessories out of reach of children. Children and people with learning difficulties must only use the medical device under the strict supervision of an adult with full mental faculties. Keep the ampoule out of reach of children under 36 months as it contains small parts that may be swallowed accidentally. **Never leave the device unattended in places accessible to minors and / or the disabled.**
- None of electric or mechanical parts have been designed to be repaired by customers or end-users. Don't open the device, do not mishandle the electric / mechanical parts. Always contact MORETTI SpA technical assistance.
- Using the device in environmental conditions different than those indicated in this manual may harm seriously the safety and the technical characteristics of the same;
- The materials used to contain the drugs are made with highly stable thermoplastic polymers that are resistant against chemicals. Such materials were tested with commonly used drugs (Salbutamol, Beclometasone dipropionate, Acetylcysteine, Budesonide, Ambr oxol) and no interaction phenomenon was observed. Interactions cannot however be excluded given the variety and the continuous evolution of the drugs that are used. Remember to:
 - To consume the drugs as quickly as possible after opening its package;
 - To avoid keeping the drug in the tray-like container for too long and to clean it immediately after every application;
 - If the tray-like container presents any abnormal situation (such as softening or cracks), do not introduce any solution and do not proceed with the inhalation. Contact the technical service and describe the methods and type of drugs used.
- Remember to:
 - Only use this device with medicines prescribed by your doctor;
 - Carry out the treatment only using the accessory indicated by the doctor according to the pathology.


The manufacturer cannot be held liable for accidental or indirect damages should the device be modified, repaired without authorization or should any of its component be damaged due to accident or misuse.
Any minimal modification / repair on the device voids the warranty and does not guarantee the compliance with the technical requirements provided by the MDD 93/42/EEC (and subsequent changes) and its normatives.

TECHNICAL CHARACTERISTICS

Model	AIR THERAPY (LT130)
Typology (MDD 93/42/EEC)	Class IIa Medical device
Power Feeding	230V - / 50 Hz
Power Consumption	184 VA
Fuse	F 1 x 1.6A L 250V
Max Pressure	250 kPa (2.5 Bar)
Max Air Flow	16 l/min
Operating Pressure	130 kPa (1.3 Bar)
Operating Air Flow	5.2 l/min a 130 kPa
Neb-Rate (with 4ml of 0.9% NaCl solution)	0.50 ml/min with 4ml of 0.9% NaCl solution
MMAD	3.25 μ m
GSD	3.45
Weight	2.20 Kg
Size	255 x 190 x 165 (h) mm
Noise Level (measured as specifications of EN 13544-1)	Approx. 57dB (A)
Duty Cycle (to 40°C and 110% operating voltage)	Non-Stop Operated
Min Capacity Nebulizer	2ml
Max Capacity Nebulizer	6ml
Working Condition	Room Temperature: 10 \div 40° C Room humidity percentage 30 \div 75% RH Altitude: 0 \div 2000m s.l.m
Conservation condition and Transport	Room Temperature: - 40 \div 70° C Room humidity Percentage: 10 \div 100% RH

SIMBOLOGY

	Class II isolation equipment
CE 0123	CE marking in conformity with EC directive 93/42/EEC and subsequent changes
	Warning, consult the instruction manual
	Type B equipment
	Conservation temperature: -25 \div 70°C
	Keep in a cool, dry place
~	Alternate Current
Hz	Mains Frequency
I	ON
0	OFF
LOT	Lot Number
SN	Serial Number
REF	Identification device
	Manufacturer by: CA-MI Srl – Via Ugo La Mal nr.13, 43010 Pilastro (PR) Italia

MMAD = Mass Median Aerodynamic Diameter
 GSD = Geometric Standard Deviation

NB: The measures and curves are not valid for the high viscosity suspension drug.

STANDARD ACCESSORIES

ACCESSORIES	
HI-FLO KIT	Nebulizer HI-FLO (LTR237), Adult Mask (LTR232), Pediatric Mask (LTR233), Air Tube (LTR236), Mouth-piece (LTR231) and Nosepiece (LTR234)

For each individual patient it's recommended to use the nebulizer for 6 months or for a maximum of 120 treatments. The nebulizer must be replaced after a long period of inactivity, if it is deformed or broken, or if the nebulizer nozzle is blocked by dry medicine, dust, ecc. **Only use the original nebulizer supplied by MORETTI SpA with the device**
 Use the "nose piece" accessory only if expressly indicated by your doctor and paying attention **NEVER** to introduce inside the nose the nasal bifurcation, but only bring it as close as possible.


In the presence of infection or microbial contamination prone pathologies, we recommend using your personal accessories and nebulizer (always consult your doctor).

The device is equipped with a filter that removes any impurities from the air that was sucked in by the compressor. The air filter must be replaced every 25 hours of functioning or when it result particularly worn. Only use original filter.

The mask and tube must be replaced as soon as the materials they are made of show signs of deterioration.

DON'T USE THE DEVICE WITHOUT AIR FILTER

- 1- Air Tube
- 2- Nebulizer Tank
- 3- Nebulization Nozzle
- 4- Nebulizer Top
- 5- Mouthpiece
- 6- Adult Mask
- 7- Pediatric Mask
- 8- Nosepiece

CLEANING OF ACCESSORIES

Proceed as follows to disinfect the accessories:

1. Turn the upper part of the nebulizer in an anti-clockwise direction;
2. Disconnect the internal pisper at the base of the nebulizer using the fingers;
3. After using the appliance, disassemble the nebulizer and clean all parts in warm water, rinse carefully and remove excess water using a soft cloth and leave to dry in a clean place.
4. To disinfect the accessories (except the tube) use denatured alcohol or hypochlorite solutions easily obtainable from chemists. Before using after cleaning, rinse the accessories in warm water until all traces of the disinfectant have been eliminated, dry and store in a dry and dust free place.

DO NOT BOIL OR PUT IN AUTOCLAVE MASK AND AIR TUBE

After each treatment clean thoroughly each component of the nebulizer (except air tube) removing medication residual and possible impurities. Clean all parts in warm water. Rinse thoroughly making sure that all deposits are washed away and let dry

The nebulizer HI-FLO, the mouthpiece and the nosepiece must be disinfected by boiling in water (for max. 10 minutes).

Attention: Never use a cleaning brush or put sharp objects into the jet holes as this will damage the nebulizer

CLEANING DEVICE

Use a soft dry cloth with not – abrasive and not – solvent detergents. The device's plug must be removed from the wall socket before proceeding with any cleaning procedures.

Guidance and manufacturer's declaration – Electromagnetic Emissions

The aerosol AIR THERAPY is intended for use in the electromagnetic environment specified below. The customers or the user of the aerosol AIR THERAPY should assure that it's used in such an environment.

Emissions Test	Compliance	Electromagnetic environment - guidance
Irradiated / Conducted emissions CISPR11	Group 1	The aerosol AIR THERAPY only used RF energy only for its internal functioning. Therefore its RF emissions are very low and are not cause interference in proximity of any Electronic appliances.
Irradiated / Conducted emissions CISPR11	Class [B]	The aerosol AIR THERAPY can be used in all environments, including domestic and those connected directly to the public mains distribution that supplies power to environments used for domestic scopes or environments feeds to you from batteries.
Harmonic emissions IEC/EN 61000-3-2	Class [A]	
Voltage fluctuations / flicker emissions IEC/EN 61000-3-3	Complies	

Guidance and manufacturer's declaration – Immunity Emissions

The aerosol AIR THERAPY is intended for use in the electromagnetic environment specified below. The customers or the user of the aerosol AIR THERAPY should assure that it's used in such an environment.

Immunity Test	Level indicated by the IEC 60601-1-2	Compliance Level	Electromagnetic environments - guidance
Electrostatic discharge (ESD) IEC/EN 61000-4-2	$\pm 6kV$ on contact $\pm 8kV$ in air	The device doesn't change its state	Floors should be wood, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC/EN 61000-4-4	$\pm 2kV$ power supply lines $\pm 1kV$ for input / output lines	The device doesn't change its state	Mains power quality should be that of a typical commercial environment or hospital
Surge IEC/EN 61000-4-5	$\pm 1kV$ differential mode	The device doesn't change its state	Mains power quality should be that of a typical commercial environment or hospital
Loss of voltage, brief voltage interruptions and variations IEC/EN 61000-4-11	5% U_T (>95% dip U_T) for 0.5 cycle 40% U_T (>60% dip U_T) for 5 cycle 70% U_T (>30% dip U_T) for 25 cycle <5% U_T (>95% dip U_T) for 5 sec	-	Mains power quality should be that of a typical commercial environment or hospital if the user of the aerosol AIR THERAPY request that the appliance operates continuously, the use of a continuity unit is recommended.
Magnetic field IEC/EN 61000-4-8	3A/m	The device doesn't change its state	The power frequency magnetic field should be measured in the intended installation location to assure that it's sufficiently low.

Note U_T is the value of the power supply voltage

This section contains information regarding the conformity of the compliance with the IEC 60601-1-2 Standard.

The AIR THERAPY aerosol is an electro-medical device that requires particular precautions regarding electro-magnetic compatibility and which must be installed and commissioned according to the electro-magnetic compatibility information supplied.

Mobile and portable RF communication appliances (mobile phones, transceivers, etc.) can affect the medical system.

The use of accessories, transducers and cables different to those specified, with the exception of transducers and cables sold by the appliance and system manufacturer as spare parts, can lead to an increase in emissions or in a decrease of the immunity of the device or system.

Guidance and manufacturer's declaration – Immunity Emissions			
The aerosol AIR THERAPY is intended for use in the electromagnetic environment specified below. The customers or the user of the aerosol AIR THERAPY should assure that it's used in such an environment.			
Immunity Test	Level indicated by the IEC 60601-1-2	Compliance level	Electromagnetic environments - guidance
Conducted Immunity IEC / EN 61000-4-6	3Vrms 150kHz to 80MHz (for non life-supporting devices)	$V_1 = 3 \text{ V rms}$	<p>The portable and mobile RF communication devices, including cables, must not be used closer to the AIR THERAPY device, than the separation distance calculated by the equation applicable to the transmitter frequency. Recommended separation distance</p> $d = [3.5 / V_i] \sqrt{P}$ $d = [12 / E_i] \sqrt{P} \text{ from 80 MHz to 800MHz}$ $d = [23 / E_i] \sqrt{P} \text{ from 800 MHz to 2.5 GHz}$ <p>Where P is the maximum nominal output voltage of the transmitter in Watt (W) depending on the manufacturer of the transmitter and the recommended separation distance in metres (m). The intensity of the field from the fixed RF transmitters, as determined by an electro-magnetic study of the site^{a)}, could be lower than the level of conformity of each frequency interval^{b)}. It is possible to check for interference in proximity to devices identified by the following symbol:</p> 
Radiated Immunity IEC / EN 61000-4-3	3V/m 80MHz to 2.5GHz (for non life-supporting devices)	$E_1 = 3 \text{ V / m}$	
<p>Note 1: At 80 MHz and 800 MHz the interval with the highest frequency is applied Note 2: These guide lines may not be applicable in all situations. The electro-magnetic propagation is influenced by the absorption and by reflection from buildings, objects and people. a) The field intensity for fixed transmitters such as the base stations for radiotelephones (mobile and cordless) and terrestrial mobile radio, amateur radio devices, radio AM and FM transmitters and TV transmitters can not be theoretically and accurately foreseen. To establish an electro-magnetic environment generated by fixed RF transmitters, an electro-magnetic study of the site should be considered. If the field intensity measured in the place where the device will be used surpasses the above mentioned applicable level of conformity, the normal functioning of the device should be monitored. If abnormal performance arises, additional measures such as changing the device's direction or positioning may be necessary. b) The field intensity on an interval frequency of 150 kHz to 80 MHz should be less than 3 V/m.</p>			

Recommended separation distance between portable and mobile radio-communication devices and the monitor			
The aerosol AIR THERAPY is intended to operate in an electro-magnetic environment where RF irradiated interferences are under control. The client or operator of the AIR THERAPY device can help prevent electro-magnetic interference by keeping a minimum distance between the portable and mobile RF communication devices (transmitters) and the AIR THERAPY device, as recommended below, in relation to the radio-communication maximum output power.			
Maximum nominal output power of the Transmitter W	Separation distance from the frequency transmitter (m)		
	150 kHz to 80 MHz $d = [3.5 / V_i] \sqrt{P}$	80 MHz to 800 MHz $d = [12/E_i] \sqrt{P}$	800 MHz to 2.5 GHz $d = [23/E_i] \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters with a maximum nominal output power not shown above, the recommended separation distance in metres (m) can be calculated using the equation applicable to the transmitter frequency, where P is the maximum nominal output power of the transmitter in Watt (W) depending on the transmitter's manufacturer.			
Note 1: At 80 MHz and 800 MHz the interval with the highest frequency is applied			
Note 2: These guide lines may not be applicable in all situations. The electro-magnetic propagation is influenced by the absorption and by the reflection from buildings, objects and people.			

INSTRUCTION FOR USE

- Place the device on a flat, stable and clean surface and plug it in the wall socket. Make sure the power cord is thoroughly unrolled to avoid dangerous overheating. In case the power cord is damaged, contact MORETTI SpA technical assistance for replacement.
- Prepare the HI-FLO nebulizer opening the upper part and pouring the drug prescribed by your doctor into the lower tank. Close the nebulizer.
- Connect the air tube into the air outlet placed above the knob of nebulization.
- Connect the other end of air tube into the bottom of the nebulizer.
- Connect the selected accessory to the nebulizer: child mask, adult mask, mouthpiece or nosepiece.
- Make sure the air filter is placed properly in its seat in the bottom of the device.
Air Filter replacement: Open the air filter cover, remove the filter and insert the new one. Place back the cover.
- Press the ON/OFF switch to position **I** to start nebulization.
- To interrupt or stop the treatment press again the ON/OFF switch.
- Adjust the nebulization speed by turning the knob towards MIN for longer treatments, or towards MAX for quicker treatments.
- After treatment has been completed press the ON/OFF switch to position **0** and pull out the plug from the wall socket.
- Wash the nebulizers and the accessories as explained in the Cleaning section.
- Place back power cord and accessories into the compartments.

Always use the nebulizer facing upwards so that substances and / or medicines cannot escape from the nebulizer during the normal use.

**NEVER INHALE IN HORIZONTAL POSITION
NEVER BEND THE NEBULIZER OVER 60°**

MAKE SURE THAT CHILDREN AND/OR MENTALLY ILL PEOPLE DO NOT USE THE DEVICE WITHOUT ADULT SURVEILLANCE

INSTRUCTIONS FOR THE CORRECT DISPOSAL OF THE PRODUCT FOR THE PURPOSE AND EFFECT OF EUROPEAN COMMUNITY DIRECTIVE 2002/96/EC:

At the end of its useful life the product must not be disposed of with domestic waste.

It may be delivered to the relevant separate waste collection centres prepared by municipal administrations, or at retailers that supply this service. Disposing of the product separately prevents possible negative consequences for the environment and for health, deriving from inadequate disposal.

It also allows the recovery of materials of which it's made up in order to obtain an important saving of energy and resources. In order to point out the obligation to dispose of the electro-medical appliances separately, the product has the sign of the crossed mobile waste container applied to it.

MAINTENANCE

The **AIR THERAPY** atomiser does not need maintenance or lubrication.

Before use always check correct functioning and safety of the device. Carry out disinfection as described in the "CLEANING ACCESSORIES" section. Unpack the instrument and **always check** integrity of plastic parts and feeding cable, they might have been damaged during previous use. Connect cable to electrical network and turn switch on.

Close the air outlet with one finger to make sure that noise produced is regular and there is no malfunctioning.

With the air outlet always closed check the correct functioning of the nebulization regulator by turning the knob from MIN to MAX. Make sure the indicator of the pressure meter is working correctly.

Verify that the atomiser is not damaged by previous use (it was badly put away or badly knocked).

A protection fuse (**F 1x1.6A L 250V**) reachable from exterior and it situated in the plug protects the instrument. For use replacing, always check the type and the range indicated.

The manufacturer will provide upon request electric diagrams, components list, description, setting instructions and any other information that can help the technical assistance staff for product repair.

Fault type	Cause	Solution
1. The device doesn't work	a) The plug may be misplaced in the wall socket b) Thermal protector may be on (the device has been working beyond its limits and / or near heat sources)	a) Make sure the plug is properly placed in the wall socket. Make sure the ON/OFF switch is in position I . b) Switch off the device by pressing the switch to position 0 and left the motor cool down for at least 30 minutes.
2. Low Nebulization	Clogged Nebulizer Tank	Clean and disinfect the nebulizer tank as explained in the instruction manual
3. Low Nebulization	Clogged Nebulizer Tank	If cleaning was not successful change cruet
4. Absence of Nebulization	Clogged Nebulizer Tank	Check that the nebulizer contains medication; Make sure that the nebulizer is not clogged; Check the connection between the compressor air outlet port and the accessories
5. Slow Nebulization	Highly dense drug	Dilute drug in physiological liquid
6. Noisy Device	Extended use	Call retainer or Technical service MORETTI SpA
Fault 1 - 2 - 3 - 4 - 5 - 6	No solution with previous items	Call retainer or Technical service MORETTI SpA

If the unit doesn't nebulizer once the above conditions have been checked, we suggest to contact your dealer or technical service MORETTI SpA.

BEFORE EVERY CHECKING OPERATION, IN CASE OF ANOMALIES OR BAD FUNCTIONING, PLEASE CONTACT MORETTI TECHNICAL SERVICE.

MORETTI SpA DOES NOT GIVE GUARANTEE IF INSTRUMENT, AFTER THE TECHNICAL SERVICE CHECKING, APPEARS TO BE TAMPERED.

ANY MINIMAL MODIFICATION / REPAIR ON THE DEVICE VOIDS THE WARRANTY AND DOES NOT GUARANTEE THE COMPLIANCE WITH THE TECHNICAL REQUIREMENTS PROVIDED BY THE MDD 93/42/EEC DIRECTIVE (AND SUBSEQUENT CHANGES) AND ITS NORMATIVES.

RULES FOR RETURNING AND REPAIRING

COMPLYING WITH THE NEW EUROPEAN RULES, MORETTI INDICATES THE IMPORTANT POINTS TO PROTECT INSTRUMENT AND OPERATORS HYGIENE. THESE RULES MUST BE RESPECTED IN ORDER TO GUARANTEE HYGIENE AND SAFETY TO ALL THE PEOPLE OPERATING WITH THE INSTRUMENT TO OBTAIN QUALITY AND WELL BEING.

Every returned instrument will be hygienically checked before repairing. If MORETTI finds instrument not suitable for repairing due to clear signs of internal or external contamination, the same will be returned to customer with specification of **NOT REPAIRED INSTRUMENT**, accompanied by an explanation letter. MORETTI will decide if contamination is due to bad functioning or misuse. If contamination is due to bad functioning, MORETTI will substitute the instrument, only if a **SALE RECEIPT** and **STAMPED GUARANTEE** accompany the same.

MORETTI is not responsible for contaminated accessories, they will be substitute at customer's expenses.

For this reason it is **COMPULSORY** to carefully disinfect the external part of the instrument and accessories with a cloth soaked in methylated spirits or hypochlorite-based solutions. Put the instrument and accessories in a bag with indication of disinfecting. We also request to specify the kind of fault, in order to speed up repairing procedures. To this end, please read the instructions carefully in order to avoid damaging the equipment through improper use. Always specify the fault encountered so that MORETTI can establish whether it falls into the category of the faults covered by the guarantee.

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