CE  
0123

LEM

**Aspiratore Chirurgico Aspimed 2.2**  
**Surgical Aspirator Aspimed 2.2**  
**Aspirateur Aspimed 2.2**  
**Absauger Aspimed 2.2**  
**Aspirador Aspimed 2.2**

**MANUALE D'USO**  
**INSTRUCTION MANUAL**  
**MODE D'EMPLOI**  
**HANDBUCH**  
**MANUAL DE ISTRUCCIONES**

Moretti S.p.A.  
Via Bruxelles 3 - Melegnano  
52022 Cavriglia (Arezzo)  
Tel. +39 055 96 21 11

Fax. +39 055 96 21 200  
[www.morettispa.com](http://www.morettispa.com)  
[info@morettispa.com](mailto:info@morettispa.com)

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**ASPIMED 2.2** Surgical aspirator is a portable unit, working with 230V ~ / 50 Hz network electricity, designed for the aspiration of bodily fluids in adult and children. It's particularly suitable for nasal, oral or tracheal aspiration of mucus, catarrh or blood after minor surgical procedures and can be used in post-operative therapy at home or conveniently transported from one hospital ward to another. Easily portable equipment designed for continuous use.

Made of highly heat resistant, electrically insulated plastic material in conformity with the latest European safety standard, the product is supplied with a complete polycarbonate autoclavable jar with overflow valve and it is equipped with aspiration regulator and vacuum indicator located on the front panel.

## GENERAL WARNING

**READ INSTRUCTION MANUAL CAREFULLY BEFORE USE**

**ONLY HIGHLY QUALIFIED STAFF USE RESERVED**

**THE INSTRUMENT MUST NOT DISASSEMBLED  
FOR TECHNICAL SERVICE ALWAYS CONTACT MORETTI S.p.A.**

## IMPORTANT SAFETY RULES

1. Check the condition of the unit before each use. The surface of the unit should carefully inspected for visual damage. Check the mains cable and **do not connect to power** if damage is apparent;
2. Before connecting the appliance always check that the electric data indicated on the data label and the type of plug used, correspond to those of the mains electricity to which it's to be connected;
3. If the plug supplied with the appliance is incompatible with the mains electricity socket, contact qualified staff for replacement of the plug with a suitable type. The use of simple or multiple and / or extension adapters is not generally recommended. Whenever their use is indispensable, use those in compliance with safety regulations, however paying attention not to exceed the maximum power supply limits, which are indicated on the adapters and extensions;
4. Respect the safety regulations indicated for electrical appliances and particularly:
  - Use original components and accessories provided by the manufacturer to guarantee the highest efficiency and safety of the device;
  - The device can be used only with the bacteriological filter;
  - Never immerge the appliance into water;
  - Place instrument on stable and flat surfaces;
  - Position the device in a way that the air inlets on the back aren't obstructed;
  - Don't use in the presence of inflammable substances such as anaesthetic, oxygen or nitrous oxide;
  - Don't touch the device with wet hands and always prevent the appliance coming into contact with liquids;
  - Keep off the reach of children or not capable people without supervision;
  - Don't leave the appliance connected to the power supply socket when not in use;
  - Don't pull the power supply cable to disconnect the plug remove the plug from the mains socket correctly;
  - Preserve and use the medical device in environments protected from atmospheric factors and at a distance from heat sources;
  - Don't use the device thoracic drainage.
5. For repairs, exclusively contact technical service and request the use of original spare parts.  
Failure to comply with the above can jeopardise the safety of the device;
6. **Use only for the purpose intended.** Don't use for anything other than the use defined by the manufacturer.  
The manufacturer will not be responsible for damage due to improper use or connection to an electrical system not complying with current regulation.
7. Particular precautions must be made concerning electromagnetic compatibility. The medical device must be installed and used according to information supplied with the accompanying documents;
8. Instrument and accessory discharging must be done according to current regulations in the country of use.
9. None of electric or mechanical parts have been designed to be repaired by customers or end-users. Don't open the device, do not mishandle the electric / mechanical parts. Always contact technical assistance
10. Using the device in environmental conditions different than those indicated in this manual may harm seriously the safety and the technical characteristics of the same.

**The manufacturer cannot be held liable for accidental or indirect damages should the device be modified, repaired without authorization or should any of its component be damaged due to accident or misuse.**

**Any minimal modification / repair on the device voids the warranty and does not guarantee the compliance with the technical requirements provided by the MDD 93/42/EEC (and subsequent changes) and its normatives.**

## TECHNICAL CHARACTERISTICS

Model	<b>ASPIMED 2.2 (LTA260)</b>
Typology (MDD 93/42/EEC)	Medical Device Class Ila
Classification UNI EN ISO 10079-1	HIGH VACUUM / LOW FLOW
Main Voltage	230 V ~ / 50 Hz
Power consumption	184 VA
Fuse	F 1 x 1.6A L 250 V
Maximum suction aspiration (without jar)	-75kPa (- 0.75 bar)
Minimum suction aspiration (without jar)	Less -25kPa (- 0.25 bar)
Maximum flow (without jar)	16 l/min
Weight	2.5 Kg
Dimension	357 x 193 x 215 mm
Functioning (to 35°C and 110 % operating voltage)	<b>NON-STOP OPERATED</b>
Working condition	Room temperature: 5 ÷ 35 °C Room humidity percentage: 30 ÷ 75 % RH Atmospheric pressure: 700 ÷ 1060 hPa
Conservation condition	Room temperature: - 40 ÷ 70 °C Room humidity percentage: 10 ÷ 100% RH

## SYMBOLS

	Class 2 isolation equipment
<b>CE 0123</b>	CE marking in conformity with EC directive 93/42/CEE and subsequent changes
	Warning, consult the instruction manual
	Type B equipment
	Fuse
	Keep in a cool, dry place
	Conservation temperature: - 40 ÷ 70° C
~	Alternate current
Hz	Mains frequency
I	ON
0	OFF
	Manufacturer: CA-MI Srl - Via Ugo La Malfa n°13, 43010 Pilastrò (PR) Italy
<b>LOT</b>	Batch production
<b>SN</b>	Serial number
<b>REF</b>	Model / Ref number
<b>IP21</b>	Degree of protection an electrical device provides in the case of accidental or intentional contact with the human body or with objects, and protection in the case of contact with water.

This section contains information regarding the conformity of the compliance with the IEC 60601-1-2 Standard.

The ASPIMED 2.2 surgical aspirator is an electro-medical device that requires particular precautions regarding electro-magnetic compatibility and which must be installed and commissioned according to the electro-magnetic compatibility information supplied.

Mobile and portable RF communication appliances (mobile phones, transceivers, etc.) can affect the medical system.

The use of accessories, transducers and cables different to those specified, with the exception of transducers and cables sold by the appliance and system manufacturer as spare parts, can lead to an increase in emissions or in a decrease of the immunity of the device or system.

<b>Guidance and manufacturer's declaration – electromagnetic Emissions</b>		
The surgical aspirator ASPIMED 2.2 is intended for use in the electromagnetic environment specified below.		
<b>Emissions test</b>	<b>Compliance</b>	<b>Electromagnetic environment - guidance</b>
Irradiated / Conducted emissions CISPR11	Group 1	The surgical aspirator ASPIMED 2.2 only used RF energy only for its internal functioning. Therefore, its RF emissions are very low and are not cause interference in proximity of any Electronic appliances.
Irradiated / Conducted emissions CISPR11	Class [B]	The surgical aspirator ASPIMED 2.2 can be used in all environments, including domestic and those connected directly to the public mains distribution that supplies power to environments used for domestic scopes.
Harmonic emissions IEC/EN 61000-3-2	Class [A]	
Voltage fluctuations / flicker emissions IEC/EN 61000-3-3	Complies	

<b>Guidance and manufacturer's declaration – Immunity Emissions</b>			
The surgical aspirator ASPIMED 2.2 is intended for use in the electromagnetic environment specified below.			
The customers or the user of the surgical aspirator ASPIMED 2.2 should assure that it's used in such an environment.			
<b>Immunity Test</b>	<b>Level indicated by the IEC 60601-1-2</b>	<b>Compliance Level</b>	<b>Electromagnetic environment - guidance</b>
Electrostatic discharge (ESD) IEC/EN 61000-4-2	$\pm 6\text{kV}$ on contact $\pm 8\text{kV}$ in air	The device doesn't change its state	Floors should be wood, concret or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC/EN 61000-4-4	$\pm 2\text{kV}$ power supply lines  $\pm 1\text{kV}$ for input / output lines	The device doesn't change its state	Mains power quality should be that of a typical commercial environment or hospital.
Surge IEC/EN 61000-4-5	$\pm 1\text{kV}$ differential mode	The device doesn't change its state	Mains power quality should be that of a typical commercial environment or hospital.
Loss of voltage, brief voltage interruptions and variations IEC/EN 61000-4-11	$5\%U_T$ ( $>95\%$ dip $U_T$ ) for 0.5 cycle $40\%U_T$ ( $>60\%$ dip $U_T$ ) for 5 cycle $70\%U_T$ ( $>30\%$ dip $U_T$ ) for 25 cycle $<5\%U_T$ ( $>95\%$ dip $U_T$ ) for 5 sec	-	Mains power quality should be that of a typical commercial environment or hospital If the user of the surgical aspirator ASPIMED 2.2 request that the appliance operates continuously, the use of a continuity unit is recommended.
Magnetic field IEC/EN 61000-4-8	3A/m	The device doesn't change its state	The power frequency magnetic field should be measured in the intended installation location to assure that it's sufficiently low.
Nota $U_T$ is the value of the power supply voltage			

**Guidance and manufacturer's declaration – Immunity Emissions**

The surgical aspirator ASPIMED 2.2 is intended for use in the electromagnetic environment specified below.  
The customers or the user of the surgical aspirator ASPIMED 2.2 should assure that it's used in such an environment.

Immunity Test	Level indicated by the IEC 60601-1-2	Compliance Level	Electromagnetic environment - guidance
Conducted Immunity IEC / EN 61000-4-6	3Vrms 150kHz to 80MHz (for non life-supporting devices)	$V_1 = 3 \text{ V rms}$	The portable and mobile RF communication devices, including cables, must not be used closer to the ASPIMED 2.2 device, than the separation distance calculated by the equation applicable to the transmitter frequency. Recommended separation distance $d = [3.5 / V_1] \sqrt{P}$ $d = [12 / E_1] \sqrt{P}$ from 80 MHz to 800MHz $d = [23 / E_1] \sqrt{P}$ from 800 MHz to 2.5 GHz  Where P is the maximum nominal output voltage of the transmitter in Watt (W) depending on the manufacturer of the transmitter and the recommended separation distance in metres (m). The intensity of the field from the fixed RF transmitters, as determined by an electro-magnetic study of the site <sup>9)</sup> , could be lower than the level of conformity of each frequency interval <sup>10)</sup> . It is possible to check for interference in proximity to devices identified by the following symbol: 
Radiated Immunity IEC / EN 61000-4-3	3V/m 80MHz to 2.5GHz (for non life-supporting devices)	$E_1 = 3 \text{ V / m}$	

Note 1: At 80 MHz and 800 MHz the interval with the highest frequency is applied  
Note 2: These guide lines may not be applicable in all situations. The electro-magnetic propagation is influenced by the absorption and by reflection from buildings, objects and people.

- a) The field intensity for fixed transmitters such as the base stations for radiotelephones (mobile and cordless) and terrestrial mobile radio, amateur radio devices, radio AM and FM transmitters and TV transmitters can not be theoretically and accurately foreseen. To establish an electro-magnetic environment generated by fixed RF transmitters, an electro-magnetic study of the site should be considered. If the field intensity measured in the place where the device will be used surpasses the above mentioned applicable level of conformity, the normal functioning of the device should be monitored. If abnormal performance arises, additional measures such as changing the device's direction or positioning may be necessary.
- b) The field intensity on an interval frequency of 150 kHz to 80 MHz should be less than 3 V/m.

**Recommended separation distance between portable and mobile radio-communication devices and the monitor**

The ASPIMED 2.2 surgical aspirator is intended to operate in an electro-magnetic environment where RF irradiated interferences are under control. The client or operator of the ASPIMED 2.2 device can help prevent electro-magnetic interference by keeping a minimum distance between the portable and mobile RF communication devices (transmitters) and the ASPIMED 2.2 device, as recommended below, in relation to the radio-communication maximum output power.

Maximum nominal output power of the Transmitter W	Separation distance from the frequency transmitter (m)		
	150 kHz to 80 MHz $d = [3.5 / V_1] \sqrt{P}$	80 MHz to 800 MHz $d = [12/E_1] \sqrt{P}$	800 MHz to 2.5 GHz $d = [23/E_1] \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters with a maximum nominal output power not shown above, the recommended separation distance in metres (m) can be calculated using the equation applicable to the transmitter frequency, where P is the maximum nominal output power of the transmitter in Watt (W) depending on the transmitter's manufacturer.

Note 1: At 80 MHz and 800 MHz the interval with the highest frequency is applied

Note 2: These guide lines may not be applicable in all situations. The electro-magnetic propagation is influenced by the absorption and by the reflection from buildings, objects and people.

## CLEANING OF ACCESSORIES

To clean the plastic housing of the device wear disposable latex gloves and clean with denaturated alcohol or hypochlorite solutions.

Washing and / or cleaning the autoclavable jar as to be carried out as follows:

- Wear protection gloves and apron (glasses and face mask if necessary) to avoid contact with contaminating substances;
- Disconnect the jar from the device
- Disconnect all tubes from the jar and the protection filter
- Empty and dispose of the content and of the suction catheter according to the laws in force in your country;
- Separate all parts of the cover (overflow valve, o-ring);

After disposing of disposable parts and disassembling the jar wash in running cold water and rinse thoroughly. Then soak in warm water (temperature shall not exceed 60°C).

Wash thoroughly and if necessary use a non-abrasive brush to remove incrustations. Rinse in running warm water and dry all parts with a soft cloth (non-abrasive). The jar and the cover can be autoclaved by placing the parts into the autoclave and running one sterilization stem cycle at 121°C (1 bar relative pressure – 15 min) making sure that the jar is positioned upsidedown.

Mechanical resistance of the jar is guaranteed up to 30 cycles of sterilization and cleaning at the indicated conditions (EN ISO 10079-1). Beyond this limit the physical-mechanical characteristics of the plastic may decrease and replacement of the part is therefore recommended.

After sterilization and cooling at environment temperature of the parts make sure that these are not damaged. Assemble the jar as follows:

- Place the overflow valve into its seat in the cover (under VACUUM connector)
- Insert floating valve keeping the o-ring towards the opening of the cage
- Place the o-ring into its seat around the cover
- After completing assembling operations always make sure that cover seals perfectly to avoid vacuum leakages or liquid exit

The aspiration tubes can be sterilized on autoclave using a sterilization cycle at 120°C (1 bar relative pressure – 15 min). The conical connector can be sterilized on autoclave using a sterilization cycle at 121°C. (1 bar relative pressure – 15 min). The device is ready for a new employment now.

## ACCESSORIES SUPPLIES

DESCRIPTION	CODE
COMPLETE ASPIRATION JAR 1000cc	LRA372
CONICAL FITTING	LRA100
TUBES SET 6mm x 10mm (TRASPARENT SILICON)	LRA120
ASPIRATION PROBE CH20	LRA117
ANTIBACTERIAL AND HYDROPHOBIC FILTER	LRA204

The filter is produced with (PTFE) hydrophobic material to prevent fluids entering the pneumatic circuit.

It should be changed immediately if it becomes wet or if there is any sign of contamination or discolouration.

It should also be changed if the unit is used with a patient whose risk of contamination is unknown.

**Don't use the suction unit without the protection filter.** If the suction unit is used in an emergency or in a patient where the risk of contamination is not known the filter must be changed after each use.

**Suction catheter:** Single-use device to be used on a single patient. Do not wash or re-sterilize after use. Reuse may cause cross-infections. Don't use after lapse of the sell-by date

Check the expiry date on the original packaging of the suction catheter and check the integrity of the sterile packaging. CA-MI declines any liability for injury to the patient correlated to the deterioration of the above-mentioned sterile packaging due to handling of the original packaging by third parties.

**WARNING:** Suction tubes for insertion in the human body purchased separately from the machine should comply with ISO 10993-1 standards on material biocompatibility.

## CLEANING DEVICE

Use a soft dry cloth with not – abrasive and not – solvent detergents.

**BEFORE STARTING CLEANING OPERATIONS, MAKE SURE THAT INTERNAL PARTS OF DEVICE DON'T COME INTO CONTACT WITH LIQUIDS**

## PERIODICAL MAINTENANCE CHECKS

The **ASPIMED 2.2** suction equipment does not need maintenance or lubrication.

It is, however, necessary to inspect the unit before each use. Unpack the instrument and **always check** integrity of plastic parts and feeding cable, they might have been damaged during previous use. Connect the cable to electrical network and turn the switch on.

Close the aspirator outlet with your finger and with suction regulator at maximum check that the vacuum indicators reaches at least -75kPa (-0.75 bar). Rotate the knob from right to left. The vacuum indicator should go down -25kPa (-0.25 bar).

Check that no loud noises are present. A protection fuses (**F 1 x 1.6 A L 250V**) reachable from exterior and situated in the plug protects the instrument. For fuses replacing, always the type and the range.

Type of fault	Cause	Remedy
1. The suction unit doesn't work	Cable is damaged External power source failure	Replace the cable Check the external power source
2. No aspiration	Jar Cap not screwed on correctly	Unscrew the cap, and re-screw it
3. No aspiration	Lid seal not in its seat	Unscrew the cap and insert the seal properly in its seat
4. The Vacuum power on the patient side is either very low or absent	<ul style="list-style-type: none"> <li>• Vacuum regulator set to minimum</li> <li>• Protection filter blocked or damaged</li> <li>• Connection tubes blocked, kinked or disconnected</li> <li>• Shut-off valve blocked or damaged</li> <li>• Pump motor damaged</li> </ul>	<ul style="list-style-type: none"> <li>• Turn the vacuum regulator clockwise and check the value of the vacuum on the gauge</li> <li>• Replace the filter</li> <li>• Replace or reconnect the tubes, check the jar connections</li> <li>• Empty the jar, or disconnect the tube from the jar and unblock the shut-off valve. The unit will only work in the upright position</li> <li>• Refer to authorised service personnel</li> </ul>
5. The float doesn't close	If the cap has been washed, ensure that the float is not partially detached	Fit the float into it's place
6. The float doesn't close	The float it's covered by dirty material	Unscrewed the cap, leave the and put in on autoclave
7. Low suction	Foam inside the jar	Fill the jar to 1/3 full of ordinary water
<b>Faults 1 - 2 - 3 - 4 - 5 - 6 - 7</b>	<b>None of the procedures have achieved the desired results</b>	<b>Contact MORETTI SpA customer service</b>

If the overflow security system fails, continued aspiration will be stopped by the antibacterial filter.

Should this back-up security system also fail, there is a risk of penetration of liquid into the device. Don't attempt to proceed with aspiration under these circumstances.

Return the device to MORETTI SpA for service.

*The Manufacturer will provide upon request electric diagrams, components list, descriptions, setting instructions and any other information that can help the technical assistance staff for product repair.*

**BEFORE EVERY CHECKING OPERATION, IN CASE OF ANOMALIES OR BAD FUNCTIONING, PLEASE CONTACT MORETTI TECHNICAL SERVICE.**

**MORETTI DOES NOT GIVE GUARANTEE IF INSTRUMENT, AFTER THE TECHNICAL SERVICE CHECKING, APPEARS TO BE TAMPERED.**

## RULES FOR RETURNING AND REPAIRING

**COMPLYING WITH THE NEW EUROPEAN RULES, MORETTI INDICATES THE IMPORTANT POINTS TO PROTECT INSTRUMENT AND OPERATORS HYGIENE. THESE RULES MUST BE RESPECTED IN ORDER TO GUARANTEE HYGIENE AND SAFETY TO ALL THE PEOPLE OPERATING WITH THE INSTRUMENT TO OBTAIN QUALITY AND WELL BEING.**

Every returned instrument will be hygienically checked before repairing. If MORETTI finds instrument not suitable for repairing due to clear signs of internal or external contamination, the same will be returned to customer with specification of NOT REPAIRED INSTRUMENT, accompanied by an explanation letter. MORETTI will decide if contamination is due to bad functioning or misuse. If contamination is due to bad functioning, MORETTI will substitute the instrument, only if a SALE RECEIPT and STAMPED GUARANTEE accompany the same.

MORETTI is not responsible for contaminated accessories, they will be substitute at customer's expenses.

For this reason it is **COMPULSORY** to carefully disinfect the external part of the instrument and accessories with a cloth soaked in methylated spirits or hypochlorite-based solutions. Put the instrument and accessories in a bag with indication of disinfecting. We also request to specify the kind of fault, in order to speed up repairing procedures. To this end, please read the instructions carefully in order to avoid damaging the equipment through improper use. Always specify the fault encountered so that MORETTI can establish whether it falls into the category of the faults covered by the guarantee.

### IMPORTANT INFORMATION FOR CORRECT DISPOSAL OF THE PRODUCT IN ACCORDANCE WITH EC DIRECTIVE 2002/96/EC:

At the end of its working life, the product must not be disposed of as urban waste. It must be taken a special local authority differentiated waste collection or to a dealer providing this service.


Disposing of a household appliance separately avoids possible negative consequences for the environment and health deriving from inappropriate disposal and enables the constituent materials to be recovered to obtain significant saving in energy and resources. As a reminder of the need to dispose of separately the product is marked with a crossed-out wheeled dustbin.

## INSTRUCTION FOR USE

- Place the unit on a flat, horizontal surface
- Connect one end of the short silicon tube, with antibacterial filter, to the suction connector on the lid of the jar.
- The other tube already connected to the filter has to be connected to the "VACUUM" jar outlet, where has been fixed the red float (security float). When the 90% of the volume of the jar is reached there is the activation of the security float (the float close the aspiration connector on the jar) to avoid liquid penetration inside the device.


**WARNING:** Ensure that the FLUID SIDE or IN marker on the filter is on the side facing the collection jar lid and fitted into the "VACUUM". A wrong connection causes immediate destruction in case of contact with sucked liquids.

### Filter assembling

- Connect the long silicone tube to the "PATIENT" jar outlet
- Connect the other end of the long silicon tube to the probe plastic connector, then connect the suction probe to it.

- Connect the power cord to the device then connect the plug to the electrical mains supply.
- Push switch on position **I** to start suction.
- Unscrew the lid of the jar and fill the jar 1/3 full or ordinary water (this assists the unit to reach peak vacuum performance and makes clean-up easier) then re-screw the lid on the jar correctly.
- During operation the jar has to be in vertical position to avoid overflow valve to cut off aspiration. Should this happen, switch off the device and disconnect the tube from the jar cover (from "VACUUM" outlet).
- Once finished push switch on **0** position and unplug.
- Remove the accessories and clean.

**NEVER USE THE DEVICE WITHOUT JAR AND / OR PROTECTION FILTER**

**MAKE SURE THAT CHILDREN AND/OR MENTALLY ILL PEOPLE DO NOT USE THE DEVICE WITHOUT ADULT SURVEILLANCE**

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