LEM

LTM500

LTM503

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I Class Medical Device
according to Directive 93/42/EEC and further modifications

CODES

LTM500 2-outlet pump for half leg or half arm garments

LTM503 3-outlet pump for full leg or full arm garments

INTRODUCTION

Thank you for purchasing a LEM compressor therapy by Moretti S.p.A. product range. This user's manual provides some suggestions as how to correctly use the product you have chosen and gives some valuable advice for your safety. Please read through the manual carefully before using the product. Should you have any queries, please contact your retailer for advice and assistance. In case of doubts, we recommend to contact the dealer who can help and advice you properly.

INTENDED USE

The compression therapy is intended for patients who need lymphatic drainage to the limbs.

CAUTION!

The compression therapy is a medical device, do not use without a physician opinion in case of health problems.

- Do not use the product for a purpose not indicated in this manual
- Moretti S.p.A declines all responsibilities for any consequences resulting from an incorrect use of this product and from unauthorized alteration to the frame of the product.
- The manufacturer reserves the right to change the information contained in this document without previous notice.

EC DECLARATION OF CONFORMITY

Moretti S.p.A. hereby declares, under its own exclusive responsibility, that the products manufactured and sold by Moretti S.p.A. and belonging to the kit pump and mattress category are in compliance with the applicable provisions outlined in the Medical Devices Directive 93/42/EEC and further modifications. To this end Moretti S.p.A. guarantees and declares as follows, under its own responsibility:

1. The devices in question satisfy the essential requirements of Annex I of directive 93/42/ EEC as prescribed by Annex VII of the above directive
2. The devices in question ARE NOT MEASUREMENT DEVICES
3. The devices in question ARE NOT INTENDED FOR CLINICAL INVESTIGATION
4. The devices in question are sold in NON STERILE packing
5. The devices in question should be considered as belonging to class I
6. DO NOT use/install the devices in question for any other purpose than that expressly indicated by Moretti S.p.A.
7. Moretti S.p.A. maintains and submits the documentation certifying conformity with Directive 93/42/EEC to the proper authorities and further modifications

STANDARDS

The system has been tested and successfully approved to the following standards:
 EN 60601-1
 EN 60601-1-2

GENERAL WARNINGS

- DO NOT use this product or any available optional equipment without first completely reading and understanding this instruction manual. If you are unable to understand the warnings, cautions or instructions, contact a healthcare professional, dealer or technical personnel before attempting to use this equipment, otherwise injury or damage may occur
- Advise with physician or therapist to determinate the correct adjustment and the correct use of the device
- Keep the packed kit away from heat sources
- SERVICE LIFE- the device use limit is defined by the wear of the parts
- DO NOT allow children to play on or operate the pump

ELECTROMAGNETIC COMPATIBILITY WARNINGS

This equipment has been tested and it complies with medical devices limits belonging to EN 60601-1-2:2007. These limits have been designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guaranteed that the interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference following one or more of the suggestions below:

- Reorient or relocate the receiving device
- Increase the separation between the equipment
- Connect the equipment into an outlet on a circuit different from that to which the other devices are connected
- Consult the manufacturer or field service technician for help

SYMBOLOLOGY

Code	Warning	Temperature limit	"BF" type	Follow instruction manual	Class II isolation	
						IP2X
Production date	Disposal according to EC/2002/96	ON	OFF	Manufacturer	Lot number	The device is protected against foreign objects the size of 12,5mm or more

KIT CONTENT

UNIT CODE	DESCRIPTION
LTM500	Sequential compressor therapy 2-outlet pump ADVANCE 2000
LTM503	Sequential compressor therapy 3-outlet pump ADVANCE 1000
LTM510	3-section full arm *
LTM520	2-section half leg *
LTM515	3-section full leg *
LTM525	2-section half arm *


**Accessories not included. Request according to the needs.*

GENERAL DESCRIPTION

The sequential compression therapy system has been designed for lymphatic drainage. The sections in the garments are made of partially overlap so as to eliminate any gap between one another and reduce the possibility of stagnant fluid. The sections in the garments inflate and deflate sequentially. Compression starts in the distal chamber and moves on to the proximal one, so favoring fluid drainage. The pressure can be adjusted by means of a knob situated on the front of the pump.

-Pump





Note Pictures are only for demonstration and do not refer to a specific product

1. Front panel
2. Mounting brackets
3. Air hose port
4. Power cord

INSTALLATION

1. Place the pump on a flat surface or hang it where you can by means of the hooks
2. Connect the tubes to the pump and to the arm/leg garment. Ensure the air hoses are not kinked or tucked under mattress
3. Plug the power cord into the electrical outlet
4. Turn the switch ON

**WARNING!**

- Make sure the pump unit is suitable for the local power voltage
- Do not position the equipment so that is difficult to disconnect the plug

BEFORE USE

Check the correct assembly, paying particular attention to air hoses connector

- Check for any damage to the garments, which can compromise the operation
- Check the zippers correctly work

USING WARNINGS

Compression therapy is contraindicated for people suffering from the following pathologies:

- Cardiopathy
- Suspected or ascertained thrombosis
- Arteriosclerosis or ischemia
- Dermatitis
- Gangrene
- Local infections

DANGER - to reduce risk of electrocution -

1. Always unplug this product immediately after using
2. Do not use while bathing
3. Do not place or store this product where it can fall or be pulled into a tub or sink
4. Do not place in or drop into water or other liquid
5. Do not reach for a product that has fallen into water. Unplug immediately

WARNING - to reduce risk of burns, electrocution, fine or injury to people -

1. Evaluate patients for entrapment risk according to facility protocol and monitor patient appropriately
2. Close supervision is necessary when this product is used on or near children. Electrical burns or choking accident may occur if a child swallow a small part detached from the device
3. Use this product only for its intended use as described in this manual. Do not use attachment not recommended by the manufacturer
4. Never operate this product if it has a damaged cord or plug, if it is not working properly, or if it has been dropped, damaged, or dropped into water. Return the product to a service centre for examination and repair
5. Keep the cord away from heated surfaces
6. Never block the air openings of this product or place it on a soft surface, such as a bed or couch, where their openings may be blocked. Keep the air openings free of lint, hairs, and other similar particles
7. Never drop or insert any object into any opening or hose
8. Do not modify this equipment without authorization of the manufacturer
9. If there is a possibility of electro-magnetic interference with mobile phones, please increase the distance to 3,3m between devices or turn off the mobile phone

HOW TO USE


1. Turn the pressure-adjusting knob to the minimum
2. Put the garment on and fasten the zipper completely
3. Switch the pump on. The green light indicates that the pump is working. Do not put on or take off the garment while the pump is working
4. Adjust pressure turning the knob on the front of the appliance clockwise to the value indicated by your physician according to the following chart:

Knob	1	2	3	4	5	6	7	8
mmHg	30	38	45	52	59	66	73	80

The chart indicates the millimeters of mercury as units of measure and pressure and the corresponding numbers on the pressure-adjusting knob.

INSTRUCTIONS TO USE 2 FULL LEG GARMENTS WITH 1 PUMP STORAGE

- Fold the garments accurately to prevent creases

**MAINTENANCE**

The LEM devices by MORETTI are checked carefully and supplied with EC mark, once launched on the market. For patient and physician's safety we recommend to let check the device by the manufacturer or by laboratory authorized every 2 years. In case of repair, use only original spare parts and accessories.

CLEANING AND DISINFECTION

Before using the device with people follow this procedure:

- Clean the sleeve with mild soap
- Do not immerse the connection tubes
- Do not immerse or soak pump unit
- After cleaning, dry the mattress without direct sunlight exposure
- Check the connection tubes are not obstructed by dust or else. In this case use a soft swab to clean the tubes

Note Never use acids, alkalis or solvents such as acetone

CONDITIONS OF DISPOSAL

-General conditions of disposal

In case of disposal don't use the inserting container for municipal waste. We recommend to dispose the kit in the appropriate disposal areas for recycling.

**TREATMENT OF BATTERIES (Directive 2002/96/EC)**

At the end of its life, the product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment or at retailers that provide this service. By ensuring these batteries are disposed of correctly, you will help prevent potentially negative consequences for the environment and human health which could otherwise be caused by inappropriate waste handling of the battery. To remark the need to dispose of electrical equipment separately, the products is marked with crossed mobile waste bin

SPARE PARTS/ACCESSORIES

For spare parts and accessories refer to Moretti main catalogue

TROUBLESHOOTING

If your questions can't be answered with below information, please contact your local agent directly. They might require a technician to take care of the problem

PROBLEM	SOLUTION
Pump is not ON	<ul style="list-style-type: none"> • Check if the plug is connected to mains • Check if there is any blown fuse
The garment doesn't inflate	<ul style="list-style-type: none"> • Check if the connection between garment and pump unit is tightly secured • Wait for few minutes before noticing the cycling working of the pump

ELECTROMAGNETIC COMPATIBILITY DECLARATION

Guidance and manufacturer's declaration- Electromagnetic emissions:

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

EMISSION TEST	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT- GUIDANCE
RF emissions CISPR 11	Group 1	The device uses RF Energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR11	Class B	The device is suitable for use in all establishment. Including domestic establishments and those directly connected to the public low-voltage power supply network.
Harmonic emissions IEC61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC61000-3-3	Complies	

Guidance and manufacturer's declaration - Electromagnetic immunity -

The Advance 1000/2000 devices are intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Immunity Test	IEC60601 test level	Compliance	Electromagnetic Environment-Guidance
Electrostatic Discharge (ESD) IEC61000-4-2	± 6kV contact ± 8kV air	± 6kV contact ± 8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/ burst IEC61000-4-4	± 2kV for power supply line ± 1kV for input/out line	± 2kV for power supply line ± 1kV for input/out line	Mains power quality should be that of atypical commercial or hospital environment
Surge IEC61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s)	Mains power quality should be that of atypical commercial or hospital environment.

Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	<5 % UT (>95 % dip in UT) for 0,5 cycle	<5 % UT (>95 % dip in UT) for 0,5 cycle	Mains power quality should be that of atypical commercial or hospital environment. If the user of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
	40 % UT (60 % dip in UT) for 5 cycles	40 % UT (60 % dip in UT) for 5 cycles	
	70 % UT (30 % dip in UT) for 25 cycles	70 % UT (30 % dip in UT) for 25 cycles	
	<5 % UT (>95 % dip in UT) for 5 sec	<5 % UT (>95 % dip in UT) for 5 sec	
Power frequency (50/60Hz) magnetic field IEC61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of atypical location in a typical commercial or hospital environment.
NOTE UT is the a.c. mains voltage prior to the application of the test level			

Guidance and Manufacturer’s Declaration - Electromagnetic Immunity -

These Excel devices are intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Immunity Test	IEC60601 test level	Compliance	Electromagnetic Environment-Guidance
Conducted RF IEC 61000-4-6	3Vrms150 kHz to 80 MHz outside ISM bandsa	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$ 150 kHz to 80 MHz $d = 1.2\sqrt{P}$ 150 kHz to 80 MHz $d = 2.3\sqrt{P}$ 80 MHz to 2.5G MHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). b Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey c, should be less than the compliance level in each frequency ranged. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

- NOTE 1** At 80 MHz and 800 MHz, the higher frequency range applies.
- NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.
- NOTE 3** Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and landmobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device
- NOTE 4** Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and this device:

This device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 KHz to a 80MHz $d=1.2\sqrt{P}$	80MHz to 800MHz $d=1.2\sqrt{P}$	800MHz to 2.5 GHz $d=2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

- NOTE 1** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
- NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people

TECHNICAL FEATURES**-Pump unit**

FEATURES	LTM500/503
Power supply	AC 220-240V~ 50 Hz 0,05A
Power	10 W
Fuse rating	T1A1 250V
Pressure	30-80 mmHg
Dimensions	25x13x9,5 cm
Cycle	3 minutes
Outlets	LTM500 (2 outlets) LTM503 (3 outlets)
Air output	4,5 Liters/min.
Weight	1,6 Kg
Environment temperature	Operation: +10°C / +40°C Storage: -15°C/+50°C Shipping: -15°C/+70°C
Environment humidity	Operation: 10%-90% non- condensing Storage: 10%-90% non- condensing Shipping: 10%-90% non- condensing
Atmospheric pressure	Operation: 700-1013,25 hPa
Classification	Class II, Type BF, IP2X (keep dry) Applied part: air mattress Not suitable for use in the presence of a flammable anesthetic mixture (No AP or APG protection)

-Garments

FEATURES	LTM510	LTM520	LTM515	LTM525
Type	Full arm	Full leg	Half leg	Half arm
Material	Nylon-PU	Nylon-PU	Nylon-PU	Nylon-PU
Sections	3	3	2	2
Dimension	78x22 cm	84x31 cm	50x31 cm	53x20 cm
Circumference	64 cm	72 cm	72 cm	64 cm

WARRANTY

Moretti products are guaranteed from material or manufacturing faults for 2 years from the purchasing date, except possible exclusion or restriction as follows. The warranty shall not be applied in the possible damages caused by improper use, abuse or alteration, and the warranty will not be valid if the instructions for use are not strictly adhered to. The correct intended use is specified in this manual.

Moretti is not responsible for consequent damages, personal injuries or whatever caused by or in relation to wrong installation or improper use.

Moretti warranty does not cover damages resulting from: natural disaster, not authorized maintenance or repairs, faults caused by problems on electricity supply (when necessary), use of spare parts not covered by Moretti, improper use, not authorized alteration, shipment damages (different from original Moretti shipment), or in case of insufficient maintenance as indicated in the manual.

The warranty doesn't cover components subject to wear and tear during the correct use of the device.

REPAIRING

-Warranty repair-

If a Moretti item presents material or manufacturing faults during the warranty period, Moretti will confirm with customer if the fault can be covered from warranty. Moretti, at its unquestionable discretion, can repair or replace the item, by a Moretti dealer or to Moretti headquarters. Labor cost can be charged to Moretti if the repair is covered by warranty. A repair or a replacement doesn't extend the warranty.

-Repair not covered by warranty-

A product out of warranty can be sent after Moretti authorization. The labor and shipping costs for good out of warranty are to be paid by the customer or by the dealer. The repairs are guaranteed for 6 months from the good received.

-Non-defective device-

The customer will be informed if, after the device return and examination, Moretti declares that the device is not faulty. In this case the good will be sent back to customer, the shipping charge will be paid by the customer.

SPARE PARTS

The original Moretti spare parts are guaranteed for 6 months from the receiving date.

EXEMPT CLAUSES

Moretti does not offer any other declaration, explicit or implicit warranty or conditions, including possible declarations, warranties or conditions of merchantability, fitness for a specific purpose, non infringement and non interference, all but what expressly specified in this warranty. Moretti does not guarantee the nonstop and faultless usage. The duration of possible implicit warranties which can be imposed by the law is limited by the warranty period, in the limits of law. Some states or countries don't allow limiting the implicit warranty or the exclusion or the limitation for accidental damages. In such countries, some of those exclusions or limitation may not be applied to the user. The present warranty may be modified without prior notification.

**WARRANTY CERTIFICATE**

Product _____

Purchased on (date) _____

Retailer _____

Address _____ **Town/city** _____

Sold to _____

Address _____ **Town/city** _____

Via Bruxelles, 3 - Meleto 52022 Cavriglia (Arezzo) Tel. +39 055 96 21 11

www.morettispa.com email: info@morettispa.com

MADE IN P.R.C.

Via Bruxelles, 3 - Meleto
52022 Cavriglia (Arezzo)
Tel. +39 055 96 21 11

Fax. +39 055 96 21 200
www.morettispa.com
info@morettispa.com