## GARANZIA

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## WARRANTY

This product has **2-YEAR** warranty from purchase date which is certified by the receipt or invoice. The distributor reserves the right either to repair or replace necessary parts to correct defects in the materials or workmanship without any charge. **The warranty does not apply to damage resulting from failure to follow the operating instructions, accidents, abuse, alteration or disassembly by unauthorized individuals.** In case you need to get the product repaired and you fulfil the conditions of warranty, send it to the nearest retailer together with the purchase receipt.

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**HOT-WATER BOTTLES  
INSTRUCTION MANUAL**

Moretti S.p.A.  
Via Bruxelles 3 - Meleto  
52022 Cavriglia (Arezzo)  
[info@morettispa.com](mailto:info@morettispa.com)  
Telefono +39 055 96 21 11

Fax +39 055 96 21 200  
[www.morettispa.com](http://www.morettispa.com)  
[info@morettispa.com](mailto:info@morettispa.com)

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I class medical device  
according to ECC directive 93/42

## Codes **ST186 Monolaminated bottle** **ST187 Bilaminated bottle**

Thank you for purchasing the hot water bottles **KYARA** by Moretti. We recommend you to contact your physician before the use.

### **Warning!**

Read carefully this instruction manual before using the product.

### **Use**

Hot water bottles kyara is indicated to relieve body pain. Please note that you have to fill it with drinking water.

Verify the surface of the bottle: if you see cracks or damages, don't use it. Clean the internal part of the bottle, if possible. Fill the bottle with hot water (min. 15°C- max. 50°C).

Close the bottle with care using the cap, let the air come out. Verify if the bottle loses water. Place the bottle on the body part using a cloth between the skin and the bottle.

### **Warning**

Do not use the bottle for purposes not indicated in this manual. Moretti SpA is not responsible for purposes not indicated in this manual. We suggest you to use protections when you pour the water into the bottle. Do not put the bottle if you have skin irritations and always use a cloth between the skin and the bottle. To avoid accidental opening pay attention on the bottle if a child use it.

### **Cleaning**

Use lukewarm water and neutral soap to clean the bottle. Dry at the open air, far from sun rays and from heat sources.