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I Class Medical Device according to Directive  
93/42/EEC  
and further modifications

### 1. Codes

*Comby  
mille*

<b>CP210-36</b>	Wheelchair seat 36 cm wheels 60 cm - <b>BE - PL</b>
<b>CP213-36</b>	Wheelchair seat 36 cm wheels 60 cm - <b>BE - PE</b>
<b>CP210-40</b>	Wheelchair seat 40 cm wheels 60 cm - <b>BE - PL</b>
<b>CP215-40</b>	Wheelchair seat 40 cm wheels 60 cm - <b>BE - PL</b>
<b>CP210-43</b>	Wheelchair seat 43 cm wheels 60 cm - <b>BE - PL</b>
<b>CP213-43</b>	Wheelchair seat 43 cm wheels 60 cm - <b>BE - PE</b>
<b>CP215-43</b>	Wheelchair seat 43 cm wheels 30 cm - <b>BE - PL</b>
<b>CP210-46</b>	Wheelchair seat 46 cm wheels 60 cm - <b>BE - PL</b>
<b>CP213-46</b>	Wheelchair seat 46 cm wheels 60 cm - <b>BE - PE</b>
<b>CP215-46</b>	Wheelchair seat 46 cm wheels 30 cm - <b>BE - PL</b>
<b>CP218-46</b>	Wheelchair seat 46 cm wheels 30 cm - <b>BE - PE</b>
<b>CP250-50</b>	Wheelchair heavy duty seat 50cm - <b>BE - PL</b>
<b>CP250-55</b>	Wheelchair heavy duty seat 55cm - <b>BE - PL</b>

*Smarty*

<b>CP600-40</b>	Compact wheelchair 40 cm - <b>RP - BL - PL</b>
<b>CP610-40</b>	Compact wheelchair 40 cm - <b>RP - BC - PL</b>

*light<sup>2</sup>*

<b>CP720-40</b>	Light wheelchair quick release 40cm - <b>RG - BL - PL</b>
<b>CP720-45</b>	Light wheelchair quick release 45cm <b>RG - BL - PL</b>

(**BE** = elevating armrest; **PL** = lateral footrest; **PE** = elevating footrest; **RP** = solid tires; **RG** = pneumatic tires; **BL** = full length armrest; **BC** = desk armrest)

### 2. Introduction

Thank you for purchasing a **ARDEA ONE** wheelchair by Moretti S.p.A. product range. This user's manual provides some suggestions as how to correctly use the product you have chosen and gives some valuable advice for your safety. Please read through the manual carefully before using the product. Should you have any queries, please contact your retailer for advice and assistance. In case of doubts, we recommend to contact the dealer who can help and advice you properly.

**Note:** Check all parts for shipping damage and test before using. In case of damage, do not use the device. Contact your dealer for further instruction.

### 3. Intended use

These wheelchairs are intended for people with movement difficulties.

These wheelchairs can be for self-propulsion (with 60cm casters) or can be wheeled

#### WARNING!

- Do not use the product for a purpose not indicated in this manual
- Moretti S.p.A declines all responsibilities for any consequences resulting from an incorrect use of this product and from unauthorized alteration to the frame of the product
- The manufacturer reserves the right to change the information contained in this document without notice

### 4. EC declaration of conformity

Moretti S.p.A. hereby declares, under its own exclusive responsibility, that the products manufactured and sold by Moretti S.p.A. and belonging to the wheelchair category are compliant with the applicable provisions outlined in the Medical Devices Directive 93/42/EEC and further modifications. To this end Moretti S.p.A. guarantees and declares as follows, under its own responsibility:

1. The devices in question satisfy the essential requirements of Annex I of directive 93/42/EEC as prescribed by Annex VII of the said directive
2. The devices in question ARE NOT MEASUREMENT DEVICES.
3. The devices in question ARE NOT INTENDED FOR CLINICAL INVESTIGATION
4. The devices in question are sold in NON STERILE packing
5. The devices in question should be considered as belonging to class I
6. DO NOT use/install the devices in question for any other purpose than that expressly indicated by Moretti S.p.A.
7. Moretti S.p.A. maintains and submits the documentation certifying conformity with Directive 93/42/EEC and further modifications to the proper authorities


### 5. General warnings

- DO NOT use this product or any available optional equipment without first completely reading and understanding this instruction manual. If you are unable to understand the warnings, cautions or instructions, contact a healthcare professional, dealer or technical personnel before using this equipment otherwise injury or damage may occur.
- Advise with physician or therapist to determinate the correct adjustment and the correct use of the wheelchair
- Keep the packed wheelchair away from heat sources
- SERVICE LIFE- the wheelchair use limit is defined by the wear of the parts
- ALWAYS keep hands and fingers clear of moving parts to avoid injury
- DO NOT allow children to play on or operate the wheelchair
- Close the wheelchair before lifting it.
- DO NOT attempt to lift the wheelchair by any removable parts

## 6. General description


### 6.1 CP2xx parts view

The components can change depending on the model you choose:


(FIG.1)

- A - Frame
- B - Wrenches kit for assembling


(FIG.2)

- C - Lateral footrest


(FIG.3)

- D - Elevating footrest

(FIG.4)

- E - 30cm casters  
(30 cm diameter)

(FIG.5)


- F - 60cm casters  
(60 cm diameter)

(FIG.6)

- E - Casters for narrow way  
(optional CAR280)

### 6.2 CP6xx parts view

- A - Wheelchair
- B - Armrest
- C - Footrest
- D - Wrench kit for assembly and adjustment

### 6.3 CP7xx parts view

## 7. Assembling / disassembling

**Nota:** make reference to safety cautions in this manual

### 7.1 Opening

Get the handles and pull outward, then press the seat until the wheelchair is completely open (FIG.7)


### 7.2 Closing

Lift powerfully the seat and close the wheelchair. Push the handles close (FIG.8)


### 7.3 Install/remove the armrest (CP2xx CP7xx)

- Place the armrest rear side into the rear socket, then hook them into the front socket by overturning them (FIG.9)
- To remove, unlock the armrest by pushing the front release lever, then overturn them and extract by pulling the rear release lever (FIG.10)

### 7.4 Install/remove the armrest (CP6xx)

- To install the armrest insert them into their tubes until the clips is blocked. Before install assure that the sideboard are inward (fig. 11)
- To remove the armrest press the unlock device (A) then lift them upward (fig. 11)

A - Lock device  
B - Lock clip


### 7.5 Install/remove the lateral or elevating footrest

- Assembly the lateral or elevating footrest by installing the hinge plates on the footrest onto the hinge pins on the wheelchair frame. Push the footrest toward inside of the wheelchair until it locks into place (FIG.12)
- Remove the footrest by pushing the footrest release lever inward, rotate outward, lift the footrest off the hinge pins

### 7.6 Adjust the lateral footrest


Adjust the footrest length by loosening the screw as shown in next photo (FIG.13), lower the footrest to the desired length, and securely tighten the screw

(FIG.13)

### 7.7 Adjust the elevating footrest


- Adjust the footrest length by loosening the screw as shown in next photo (FIG.14) (Rif.P), lower the footrest at the desired height, and securely tighten the screw
- Adjust the footrest height by pushing the lever (FIG.14)(Rif.L), and release to lock at the desired height

(FIG.14)

### 7.8 Install/remove the 30cm casters (CP2xx only)

- Place the caster into the hole on the wheelchair frame
- Remove the plastic cap and push the caster until it locks into the place (FIG.15)
- Removing by pushing the button on the centre of the caster, and extract the caster (FIG.15)


(FIG.15)

### 7.9 Adjust the brakes with 30 cm casters (CP2xx only)

- When the 30 cm casters are installed, place the brake on the milling. Lock the device by tightening the Allen wrenches (FIG.16),

**Note** adjust the brake device so that it can make the right pressure to securely lock the caster

### 7.10 Install/remove the 60 cm wheels (CP2xx only)

- Place the caster into the central hole on the wheelchair frame
- Remove the plastic cap and push the caster until it locks into the place (FIG.17)
- Removing by pushing the button on the centre of the caster, and extract the caster (FIG.17)

### 7.11 Adjust the brakes with 60cm casters (CP2xx only)

- When the 60 cm wheels are installed, place the brake on the milling. Lock the device by tightening the Allen wrenches (FIG.18)

**Note** adjust the brake device so that it can make the right pressure to securely lock the caster

(FIG.18)

### 7.12 Adjust the brakes (CP2xx-CP6xx)

To adjust the brakes loose the Allen screws on the brake support, and let the brake slide forward or back until the right pressure on the wheel (FIG.19).

A - Screws on the brake support

(FIG.19)

### 7.13 Install the casters for narrow ways (CAR280 - CP2xx series only)

- Remove the plastic cap on the tip of the frame, place the caster and lock the Allen wrenches (FIG.20)

(FIG.20)

**8. Before use**

- Inspect the frame for any damage. Check the tubing for cracks and rust, in order to guarantee the safe for people ant thing
- Check the correct assembly, in order to guarantee a safe use
- Check the seat is completely open. If not push the seat down from the lateral frame.
- Check the brake correctly works

**9. Using warnings**

- Check the automotive parts usage, in order to guarantee a safe use of the device for people and things
- Assure that the seat is completely open and locked
- Assure that all parts are correctly assembled. Checking carefully the correct assembling of the wheels
- Assure the brakes correctly works

**10. How to use****WARNING!**

- Getting in and out of the wheelchair may be performed independently, provided you have a adequate mobility and upper body strength
- Do not stand on removable partsa

**DO NOT STAND ON THE FOOTRESTS WHEN GETTING IN OR OUT OF THE WHEELCHAIR, OVERTURNING RISK!!!**

**Note:** Make reference to “general warnings” chapter

**10.1 Get in the wheelchair**

1. Position the wheelchair as close as possible to you
2. Engage wheel locks . We suggest to remove the armrest to ease the patient getting in the wheelchair
3. Lean on the frame with arms and slowly go down towards the seat

**10.2 Get out of the wheelchair**

1. Position the wheelchair as close as possible along side the seat to which you are transferring
2. Engage wheel locks. We suggest to remove the armrest to ease the patient getting out of the wheelchair
3. Shift body weight into seat with transfer

**10.3 Stairs and height difference (how to face them)****WARNING!**

The following method must be used only from experienced patients. If you are not experienced, let you help from an assistant!

#### GET ON

Position the wheelchair at the curb to be crossed. Balance the wheelchair on the rear wheels, in order to lift the front wheels, enough to lean them on the obstacle. Bend your trunk forward and vigorously push the handrail. Continue pushing until the rear wheels are climbing over the obstacle.

#### GET OFF

To get off the stair follow the reverse procedure, positioning the rear side of the wheelchair at the curb, keep the handrails vigorously and get off the obstacle slowly; then balance the wheelchair on the rear wheels in order to complete the descent

#### **WARNING!**

We recommend to use this method only in presence of obstacle lower than 10cm!

### 10.4 Use the armrest

#### CP2xx series

The armrest can be adjust for height and depth, and can be tipped over:

1. **Height adjustment** - push the button on the armrest, adjust the height, then release the button to lock the armrest
2. **Depth adjustment** - push the button on the armrest. Slide the armrest. Release the button to lock the armrest
3. **Tip over** - this feature is useful when you have to get in or out of the wheelchair. Push the release level then turnover the armrest

#### CP6xx - CP7xx series

1. **Tip over** - this feature is useful when you have to get in or out of the wheelchair. Push the release level then turnover the armrest

### 10.5 Use the 60cm wheels

The rear 60cm wheels are used when self-propelled. Grasp the handrail and wheel through the desired direction

### 10.6 Use the 30cm casters

The rear 30cm casters are used when wheel by an assistant through the handles

### 10.7 Use the narrow casters (OPTIONAL - CAR280)

The narrow casters are used to cross narrow spaces (elevator, etc.) To use the narrow casters remove the rear 60cm wheels

### 10.8 Carrying

The wheelchair has been designed to be easily carried.

Once the armrest and the footrest are removed, can be closed and easily positioned in the car boot, for example. If this can't be possible, position the wheelchair in front of the car rear seat so that it is stable and cannot slip or slide. Lock the wheelchair with the seatbelt.

## 11. Maintenance

The **ARDEA ONE** devices by Moretti are checked carefully and supplied with EC mark, once launched on the market. For patient and physician's safety we recommend to let check the device by the manufacturer or by laboratory authorized every 2 years. In case of repair use only original spare parts. Before using follow the "**BEFORE USE**" chapter.

## 12. Cleaning and disinfection

### 12.1 Cover and padding

The cover and padding are in Nylon. Clean with warm water and mild soap.

### 12.2 Frame and metallic parts

Clean with warm water and mild soap. Dry with a soft cloth in order to not scratch the painting.

### 12.3 Brakes and wheels

Clean the wheels with warm water and mild soap. Pay attention in drying the wheels, it's dangerous use the wheelchair with wet wheels.

### 12.4 Disinfection

If you need to disinfect the device use a common disinfectant cleaner

**Note** Do not use solvent and abrasive product

## 13. Condition of disposal

In case of disposal don't use the inserting container for municipal waste. We recommend to dispose the wheelchair in the appropriate disposal areas for recycle.

## 14. Spare parts and accessories

For spare parts only advert to MORETTI master catalogue

## 15. Technicals features

### 15.1 Dimensions

MEASUREMENT  
IN CM  
CODES

CP210-36	101					45/55	36		53,5		18,5						
CP210-40	106	89	40	89	71	48	45/55	40		58,5		19					97,5
CP210-43	106						45/55	43	40	61	28	19	20	60			100
CP210-46	106						45/55	46		63,5		19					100
CP213-36	114						45/55	36		53,5		19					97,5
CP213-43	117	89	40	89	71	48	45/55	43	40	61	28	19,5	20	60			100
CP213-46	117						45/55	46		63,5		19,5					100
CP215-40	99						45/55	40		56		17					100
CP215-43	99	73	40	89	71	48	45/55	43	40	58,5	27	17	20	30			100
CP215-46	99						45/55	46		61		17					100
CP218-46	99	73	40	89	71	48	45/55	46	40	61	27	17,5	20	30			100
CP250-50	112,5	94,5	47	91,5	71,5	48	45/55	50	40,5	73	30,5	27	20	60			125
CP250-55	112,5	94,5	47	91,5	71,5	48	45/55	55	40,5	77	30,5	27	20	60			125
CP600-40 CP610-40	105	70	40	93	80	57	40/55	40	35	51	27	13	20	52			100
CP720-40	101	79	40	94	73	50	43/51	40	40	60	30	14,5	20	60			100
CP720-45	101	79	40	94	73	50	43/51	45	40	65	30	14,5	20	60			100

### 15.2 Technical specification

- Steel frame, (aluminium frame CP7xx series only)
- Max weight CP210-36, CP213-36; **97,5 kg** - per CP210-40, CP210-43, CP210-46, CP213-43, CP213-46, CP215-40, CP215-43, CP215-46, CP218-46, CP600-40, CP610-40, CP720-40, CP720-45: **100 kg** per CP250-50, CP250-55: **125 kg**
- PU rear casters, 20cm diameter
- PU rear casters, 30cm diameter
- PU rear wheels, 60cm diameter
- Nylon coating
- Resistance to ignition: the materials used for cover and padding have been treated with methods specified in Directive EN1021-1/2
- Reference regulation: EN12183
- Environmental temperature: -20°C/+ 40°C
- Storage: store in a cool place, away from heat sources, such as sunlight (suggested temperature between +10°C and +40°C)

## 16. Warranty

Moretti products are guaranteed from material or manufacturing faults for 2 years from the purchasing date, except possible exclusion or restriction as follows. The warranty shall not be applied in the possible damages caused by improper use, abuse or alteration, and for the warranty will not be valid if the instructions for use are not strictly adhered to. The correct intended use is specified in this manual. Moretti is not responsible for consequent damages, personal injuries or whatever caused by or in relation to mistaken installation or improper use. Moretti warranty does not cover damages resulting from: natural disaster, not authorized maintenance or repairs, faults caused by problems on electricity supply (when necessary), use of spare parts not covered by Moretti, improper use, alteration not authorized, shipment damages (different from original Moretti shipment), or in case of insufficient maintenance as indicated in the manual. The warranty doesn't cover components subject to wear and tear during the correct use of the device.

### 16.1 Repairing

**Warranty repair:** If a Moretti item presents material or manufacturing faults during the warranty period, Moretti will confirm with customer if the fault can be covered from warranty. Moretti, at its unquestionable discretion, can repair or replace the item, by a Moretti dealer or to Moretti headquarters. Labor cost can be charged to Moretti if the repair is covered by warranty. A repair or a replacement doesn't extend the warranty.

**Repair not covered by warranty:** A product out of warranty can be sent after Moretti authorization. The labor and shipping costs for goods out of warranty are to be paid by the customer or by the dealer. The repairs are guaranteed for 6 months from the goods received..

**Non defective device:** The customer will be informed if, after the device return and examination, Moretti declares that the device is not faulty. In this case the goods will be sent back to customer, the shipping charge will be to be paid by the customer.

### 16.2 Spare parts

The original Moretti spare parts are guaranteed for 6 months from the receiving date.

### 16.3 Exempt clauses

Moretti does not offer any other declaration, explicit or implicit warranty or conditions, including possible declarations, warranties or conditions of merchantability, fitness for a specific purpose, non infringement and non interference, all but what expressly specified in this warranty. Moretti does not guarantee the nonstop and faultless usage. The duration of possible implicit warranties which can be imposed by the law is limited by the warranty period, in the limits of law. Some states or countries don't allow limiting the implicit warranty or the exclusion or the limitation for accidental damages. In such countries, some of those exclusions or limitation may not be applied to the user. The present warranty may be modified without prior notification.

**WARRANTY CERTIFICATE**

Product \_\_\_\_\_

Purchasing date \_\_\_\_\_

Authorized dealer \_\_\_\_\_

Street \_\_\_\_\_ Place \_\_\_\_\_

Sold to \_\_\_\_\_

Street \_\_\_\_\_ Place \_\_\_\_\_

**MORETTI S.P.A.**

Via Bruxelles, 3 - Melegnano 20122 Cavriaglia (Arezzo) Tel. +39 055 96 21 11

**www.morettispa.com** email: **info@morettispa.com****MADE IN P.R.C.**

Distribuito da / Distributed by:

Via Bruxelles, 3 - Melegnano  
52022 Cavriglia (Arezzo)  
Tel. +39 055 96 21 11

Fax. +39 055 96 21 200  
[www.morettispa.com](http://www.morettispa.com)  
[info@morettispa.com](mailto:info@morettispa.com)