EASYRED

- .Istruzioni per l'installazione e l'uso
- .Instructions for installation and use
- .Instruction pour l'installation et l'emploi
- .Instrucciones para la instalación y el uso
- .Informationen für Installation und Gebrauch

MA EM300 01A

Distribuito da/Distributed by:
Moretti S.p.A.
Via Bruxelles 3 - Melegnano
52022 Cavriglia (Arezzo)

Telefono +39 055 96 21 11
Fax +39 055 96 21 200
www.morettispa.com
info@morettispa.com

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Thanks for choosing "MORETTI S.p.A"!

"MORETTI" is committed to continuously improve its products and its service for a full customer satisfaction. Our service centres are at your disposal whatever installation, use or maintenance needs you may have.

Copyright © 1998

The reproduction of images, texts and other parts of this manual without prior written permission by "MORETTI" is prohibited.

This material is protected by the international copyright laws. The information contained in this manual can be modified or updated without notice.

Contents

<i>INTRODUCTION</i>	13
Presentation of the manual	13
Manufacturer	13
Identification	13
Technical introduction	13
Reference regulations	13
Customer Service	13
<i>SAFETY</i>	13
General rules	13
<i>GENERAL DESCRIPTION</i>	14
Technical specifications	14
Name of the main components	14
<i>OPERATION</i>	15
Transport and handling	15
Packing list verification	15
Functionality test	15
<i>OPERATION AND USE</i>	15
Instructions for use	15
Adjustments prior to use	15
Unfolding the chair	15
Folding the chair	15
Adjustment and fastening of the belts	16
Extracting the front transport levers	16
Moving the chair	16
Requirements for operators and their training	17
<i>CLEANING AND DISINFECTION</i>	17
Sanitization products	17
Sanitization intervals	17
<i>MAINTENANCE</i>	17
Ordinary maintenance	17
Requirements for operators and their training	17
<i>CLEANING AND DISINFECTION</i>	17
Sanitization products	17
Sanitization intervals	17
<i>MAINTENANCE</i>	17
Ordinary maintenance	17
Storage	17
<i>WARRANTY CONDITIONS</i>	17
Limitations of warranty	17
Limitations of liability	17
Authorisation to return	17
Setting the chair aside	18
Disposal	18
<i>EEC DECLARATION OF CONFORMITY FOR MEDICAL DEVICES</i>	19

INTRODUCTION

Presentation of the manual

This manual is a part of the article to which it is referred.

Please read the instructions carefully as they provide important information on SAFETY DURING USE AND MAINTENANCE PROCEDURES.

Descriptions and images contained in this publication are not binding; "MORETTI" has the right to modify any information for enhancement purposes without any obligation to update this manual.

Graphics and images contained in this manual are only examples and may consequently be different from real situations.

Manufacturer

The device described in this manual is produced by:


Identification

The product can be identified from the side label which contains the following data:

Name and address of the manufacturer
Product Type
Serial Number
Batch Number
Year of manufacturing
Reference regulations

Technical introduction

The EM 300 folding chair, made of aluminium, has been realised taking into consideration the needs of patients and first aid operators.

It is an essential device for carrying patients to difficult locations which the stretcher cannot reach, such as stairs, corridors and lifts.

The chair is light and easy to handle but solid; it is equipped with four levers to improve transport in any situation and two fixed wheels for movement on smooth grounds.

An automatic safety locking allows to open and close the chair in order to obtain extremely reduced overall dimensions when folded.

The chair is covered with an orange-coloured PVC sheet, water-proof and easy to clean, equipped with a pocket at the back of the seat.

The chair has two fast-release belts to immobilise the patient.

Reference regulations

The folding chair has been designed and built by "MORETTI" according to the safety rules established by the EEC/93/42 directive dated 14.06.93 (class I) concerning medical equipment.

All the materials used have dimensions suitable to work properly within the limits reported in the technical specifications; they are fire-proof, they do not release any toxic or polluting substances and they are protected to obtain high resistance to wear and corrosion.

Customer Service

The customer and product assistance services are important parts of "MORETTI" company structure.

"MORETTI" is the only authorised agent for the management, assistance service and repair of the "MORETTI" products. If you need to order spare parts, take note of the Product Type, the Serial Number and the Part Number of the folding chair. This will allow "MORETTI" to give you faster and better quality assistance.

The Customer service is at your disposal for further information on use, maintenance and service of these products.

SAFETY

General rules

The chairs are designed to allow the movement of a patient:
through the wheels and the levers of the chair;
in an ambulance.

WARNING! The vehicle must be equipped with a tested locking system in order to prevent accidental displacement of the stretcher during handling

Operators must read this manual carefully, follow the instructions and familiarise with the correct operating and maintenance procedures of the chair.

The chair can be operated by two average height persons in good physical conditions. In particular operators must possess sufficient strength in their back, arms and legs to lift, pull and push in addition to good physical balance, reflexes and muscular co-ordination. Unusual circumstances or excessive weight of the patient can require additional personnel to ensure comfort and safety to the patient.

Always tell the patient when you adjust the chair.

Never leave the chair unattended if a patient is sitting on it.

Always use safety belts.

Improper use or maintenance can be harmful to people or objects. Always use and maintain the chair as prescribed by this manual. Use only "MORETTI" spare parts and service.

Keep this manual for reference and to support personnel training. Transfer it together with the product if you sell or give the product to new users.

WARNING! IN CASE OF BLOOD-TRANSMITTED DISEASES : To reduce the risk of exposure during use, follow the maintenance instructions contained in this manual and the regulations on personnel safety established by the person in charge of the Medical Emergency Service.

GENERAL DESCRIPTION

Technical specifications

WIDTH	56 cm
HEIGHT (CHAIR OPEN)	95 cm
HEIGHT (CHAIR CLOSED)	97 cm
WEIGHT	10 kg
MAXIMUM LOAD	170 kg

- 1) Rear transport levers
- 2) Handles for transportation on wheels
- 3) Belts for immobilisation of the patient
- 4) Arm-rest
- 5) Main frame
- 6) Unlocking pawls
- 7) Front transport levers
- 8) Wheels
- 9) Inclination hub
- 10) Nylon seat/back
- 11) Seat handle

OPERATION

WARNING! Do not exceed the specified workload; this can cause damage to the patient and/or to the chair.

Name of the main components

Transport and handling

WARNING! Lifting and handling manoeuvres must be carried out by personnel specialised and trained for this purpose.

Transport can be carried out by air, sea, or land.

The weight of the chair is reported in the technical specifications.

The chair can be moved manually.

Packing list verification

Remove the content of the package and

check it against the packing list.

If all the components are present and undamaged, the package can be disposed or recycled. It must be kept in appropriate location, inaccessible to children or animals.

If the chair has been damaged during shipment, keep the package and call the shipping agent within 48 hours from delivery.

Functionality test

Before using the folding chair:

- verify if the chair opens and closes correctly referring to the chapter "Operating mode" of the present manual;
- if the check result is positive the chair is ready for regular use, if it is negative call "MORETTI" customer service.

WARNING! The above mentioned check must be repeated at regular intervals in order to verify the efficiency of the product as specified in the section "Maintenance".

- The chair is delivered ready for use and does not require assembling.

OPERATION AND USE

Instructions for use

- The folding chair is exclusively for professional use.
- Use the chair only as indicated in this manual.
- Before putting the chair in regular service read the instructions carefully and practise the operating procedures.
- Use the folding chair only when in perfect operating conditions.

Adjustments prior to use

The EM 300 folding chair DOES NOT require any adjustment prior to use.

Unfolding the chair

Bend and hold the EM/300 chair by handles (2) and (11) as indicated by Figure A of the side sequence.

Unfold with a movement similar to opening a book as indicated in Figure B until reaching the position indicated in Figure C.

Set the chair upright so that the front of the frame lays on the ground and apply a pressure downwards on the seat (11) handle.

The chair is now in the position indicated by Figure D; lift and rotate the pawls simultaneously (6) in order to let the teeth correspond to the slots.

If the teeth did not reach the bottom of the pawls, seize the chair at point (2) and (11) and carry out the adjustment by slightly forcing the opening or closing so that that the pins of the pawls reach the bottom of their seats.

WARNING! Always make sure that the teeth of the pawls have completely entered their seats before moving the chair.

Folding the chair

Extract the pawls (6) and rotate them.

Bend the chair as illustrated by Figure C, hold the handles (2) and (11) as in Figure B then fold the chair as illustrated by Figure A.

Adjustment and fastening of the belts

WARNING! ALWAYS tell the patient

A

B

C

D

when carrying out adjustments

The EM 300 chair has two quick-release

to immobilise the patient during transport. One belt is located on the back of the seat, to block the chest zone and the second one is located on the seat itself to block the legs.

The quick-release belt is composed of two parts: a fixed part (female) and an adjustable part (male) equipped with two hooks (Figure G).

Release the belts as illustrated in Figure E-F-G.


The fast-release procedure must be carried out by applying a pressure on the hooks as illustrated in Figure F and by pulling outwards the two parts.

Release the belts by letting the exceeding part slide on the hook.

Join the two parts until you hear the typical fastening sound (click).

Pull the free part of the belt until complete and correct immobilisation of the patient.

Extracting the front transport levers


WARNING! ALWAYS tell the patient when carrying out adjustments

To carry out this operation you must:

- lift the pawl (6) with one hand;
- extract the front transport lever (7) with the other hand;
- release the pawl and move the transport lever slightly back or forth until the locking device (6) snaps shut (click).

Figure I and L show all the steps.

Moving the chair

WARNING! The chair must be operated by two average height persons in good physical conditions. In particular operators must possess sufficient strength in their back, arms and legs to lift, pull and push in addition to good physical balance, reflexes and muscular co-ordination.

WARNING! When moving the chair while the patient is sitting on it, always insert the locking pawls and fasten the safety belts.

move the chair on smooth surfaces:

- lock the wheel by putting your foot on the inclination hub (9);
- bend the chair backwards for forty degrees while maintaining the fulcrum on the hub and lever on the handles (2);
- push or pull the chair by holding the handles (2);
- at the end of the way re-lock the wheels by acting on the hub (9) and set upright slowly, without sharp movements;

To move the chair in uncomfortable locations (stairs), just do as follows:

- the first operator must lift the bar of the rear transport levers (1);
- the second operator must extract and lock the front transport levers (7) as indicated in the previous section;
- both must lift the chair and the patient slowly but simultaneously;
- once the obstacle has been overcome, put the chair on the ground and close the front transport levers.

WARNING! Before moving the chair verify that the locking pawls are correctly inserted.

Requirements for operators and their training

The operators:

- must have a practical knowledge of moving patients in emergency;
- must be physically capable of assisting a patient;
- must read this manual and become familiar with the procedures described;
- must be familiar with the folding chair before using it in real situations with a real patient;
- must follow a training programme approved by the person in charge of the training of the Medical Emergency Service;
- must have their training verified;
- must do appropriate exercises and be careful to avoid the risk of damaging their spinal cord.

WARNING! Keep a record of the training

CLEANING AND DISINFECTION

Sanitization products

WARNING! Sanitization agents are corrosive.

Scrupulously observe the instructions provided by the producer of the sanitization agent during use. Whenever possible, consult the producer for guarantees regarding the degree of corrosiveness of the solutions utilised.

It is extremely important to follow the specifications regarding the concentration, temperature and reaction times.

The modification of these characteristics may damage the article.

Sanitization intervals

The sanitization intervals to be adopted are defined by the user on the basis of requirements while bearing in mind the information provided in this manual and the instructions listed on the sanitization products employed.

MAINTENANCE

Ordinary maintenance

Check the folding chair once a month or more frequently if the product is used in heavy-duty conditions.

The inspection must include the search for anything which could compromise the integrity and correct operation of the folding chair, e.g., bent or broken parts.

During the inspection it is necessary to carry out the same tests described in “Check before operation”.

WARNING! If there is any sign of damage, put the chair out of service until the necessary repairs or replacements have been carried out.

WARNING! Keep a record of the maintenance activities carried out

every 30 days:

- clean the metal parts with water and soap, rinse and dry completely;
- clean the rubber parts with water and soap, rinse and dry completely;
- carefully clean the wheels and verify their efficiency;
- lubricate all the joints and hinges with oil (the type used for lubricating motors);
- lubricate all the moving parts with white lithium-base lubricant;
- remove the excess lubricant before operation;
- carry out a complete operation test.

WARNING! IN CASE OF BLOOD-TRANSMITTED DISEASES : To reduce the risk of exposure during use, follow the maintenance instructions contained in this manual and the regulations on personnel safety established by the person in charge of the Medical Emergency Service.

Spare parts list

To order a spare part it is necessary to provide the Product Type and the Serial Number of the folding chair in addition to the part number defined in the following Table.

Storage

Long term storage must be performed according to the following conditions:

- the chair must be packed;
- the storage area must be dry and not exposed to direct sunlight;
- do not stack more than six items.

Setting the chair aside

In case the chair must be set aside for a long period it is necessary to protect the parts which could be damaged by dust deposits.

Grease the parts that could be damaged by oxidation or excessive drying .

Disposal

If you decide not to use the chair any longer, it is necessary to render it idle. We also recommend to make those parts which can be a source of damage or danger harmless.

Evaluate the class of the object according to the disposal degrees.

Dispose together with iron material and send it to the appropriate waste disposal centres.

If you consider the chair as a special waste, disassemble it and divide it into similar parts, then dispose it according to the legislation in force.

WARRANTY CONDITIONS

Limitations of warranty

The present certificate must be kept until its expire date. It must be shown together with the invoice, receipt, shipping documents and the proof of servicing which contains the name of every technical service deemed necessary by the dealer who sold the product. Warranty will be invalid to users who are unable to show proper documentation.

All "MORETTI" products are guaranteed against material or manufacturing defects for two year from the date of purchase, except for in the following circumstances: non metallic parts or products (surface treatments, decalcomanias, labels, tissues, vinyl parts or foam, etc.) are guaranteed for 120 days. Wear or improper use of the product, are not covered by the guarantee. Repairs and services are guaranteed for 120 days or until the end of the above-mentioned period, depending on which is longest. These conditions are dependant on the product being used according to the prescribed instructions contained in this manual. The warranty period starts 30 days after the packing list verification (delivery). Shipment fees are not covered by the warranty. "MORETTI" cannot be held liable for damages which occur during shipment or which occur from improper use of the product.

Limitations of liability

"MORETTI" is only liable, under the above mentioned conditions, for replacing or repairing a product or a defective part or, at its own choice, for the total reimbursement of the purchase price. Shipment fees for estimated defective parts or products to and from

"MORETTI" are not covered. Outside the above mentioned warranty conditions "MORETTI" does not recognise other types of warranty. "MORETTI" makes no warranty of any kind with regard to this material, including, but not limited to, the implied warranties of merchantability and fitness for a particular purpose. "MORETTI" shall not be liable for incidental or consequential damages in connection with the furnishing, performance or use of its product. Repair or replacement of a defective product or refund of the purchase price are exclusive and all-inclusive. In no case the liability of "MORETTI" will exceed the price of the defective product. Each legal action for failure to comply with the above mentioned warranty conditions must begin within one year from the day the defect has been noted. The court of Bolzano (Italy) will be competent for any dispute.

Authorisation to return

No product will be accepted without previous agreement with "MORETTI". The items accepted for repair or replacements will be evaluated according to previously defined criteria, including age and conditions of the product and conformity to the original performance. If it is possible to repair the product, a written cost estimate will be submitted to the customer for approval before carrying out the repair.

EEC DECLARATION OF CONFORMITY FOR MEDICAL DEVICES

The undersigned MORETTI S.p.A Company, with headquarters in Via, Bruxelles 3 - Meleto, 52022 Cavriglia (Italy), responsible for planning, manufacturing, packaging and labelling the technical aids: "Folding chair EM 300", such as those disciplined by the Council's 93/42/EEC Directive, dated June 14, 1993, concerning class 1 medical devices subject to the procedure established in Art. 11 caption 5, of the Directive itself,

GUARANTEES AND DECLARES

that the above-mentioned aids, put in service and on the market starting from the date of this document, integrally comply with the above-mentioned Directive, which is applicable to the aids.

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Copyright © 1998 "MORETTI"