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DYNA SERIES
INSTRUCTIONS MANUAL

MORETTI S.P.A.
Via Bruxelles 3 - Melegnano
52022 Cavriglia (Arezzo)
Tel. +39 055 96 21 11

Fax +39 055 96 21 200
www.morettispa.com
info@morettispa.com

INDICE

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I Class Medical Device
according to Directive 93/42/EEC
and further modifications

CODES

LTM610 Kit mattress and pump Dyna Best 5000
LTM615 Kit mattress and pump Dyna Flo 8000

INTRODUCTION

Thank you for purchasing LEVITAS kit mattress and pump Dyna by Moretti S.p.A. product range. This user's manual provides some suggestions as how to correctly use the product you have chosen and gives some valuable advice for your safety. Please read through the manual carefully before using the product. Should you have any queries, please contact your retailer for advice and assistance. In case of doubts, we recommend to contact the dealer who can help and advice you properly.

INTENDED USE

CAUTION!

- Do not use the product for a purpose not indicated in this manual
- Only qualified personnel or trained for the treatment and prevention of bedsores, can operate with the device
- Moretti S.p.A declines all responsibilities for any consequences resulting from an incorrect use of this product and from unauthorized alteration to the frame of the product
- The manufacturer reserves the right to change the information contained in this document without previous notice

The kit pump and mattress have been designed for the bedsores treatment and prevention, which may occur during long hospital and/or house confinement.

STANDARDS

The system has been tested and successfully approved to the following standards:

- EN 60601-1
- EN 60601-1-2
- EN55011 Class B
- IEC61000-3-2
- IEC61000-3-3

GENERAL WARNINGS

- DO NOT use this product or any available optional equipment without first completely reading and understanding this instruction manual. If you are unable to understand the warnings, cautions or instructions, contact a healthcare professional, dealer or technical personnel before attempting to use this equipment, otherwise injury or damage may occur
- Advise with physician or therapist to determinate the correct adjustment and the correct use of the device
- Keep the packed kit away from heat sources
- SERVICE LIFE- the device use limit is defined by the wear of the parts
- DO NOT allow children to play on or operate the pump

ELECTROMAGNETIC COMPATIBILITY WARNINGS

This equipment has been tested and it complies with medical devices limits belonging to EN 60601-1-2:2007. These limits have been designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guaranteed that the interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference following one or more of the suggestions below:

- Reorient or relocate the receiving device
- Increase the separation between the equipment
- Connect the equipment into an outlet on a circuit different from that to which the other devices are connected
- Consult the manufacturer or field service technician for help

SYMBOLS

 Code	 Warning	 Temperature limit		 "BF" type	 Follow instruction manual	 Class II isolation
 Production date	 Disposal according to EC/2002/96	 ON	 OFF	 Manufacturer	 Lot number	IP21  The device is protected against foreign objects the size of 12,5mm or more. Protected from vertical water fall.

KIT CONTENT

The Dyna series is composed by two kind of kit which can satisfy all the requirement about treatment and prevention of complications caused by long hospital and/or house confinement, in particular bedsores.

LTM610 Kit Dyna Best 5000- interchangeable elements mattress (12,7cm h) and pump with pressure regulation


LTM608 Kit Dyna Flo 8000- interchangeable elements mattress (20,3 cm h) and pump with pressure regulation

UNIT CODE	DESCRIPTION	LTM610	LTM615
LTM966	Pump with pressure adjustment	●	
LTM969	Pump with pressure adjustment		●
LTM930	Mattress with interchangeable cells (12,7cm h) including cover	●	
LTM935	Mattress with interchangeable cells (20,3cm h) including cover		●
—————	Instruction manual	●	●

GENERAL DESCRIPTION


-Mattress and pump

N.B. Pictures are only for demonstration and do not refer to a specific product


1. CPR valve
2. Pump unit
3. Mattress system
4. Quick connector


-Front view

1. Quick connector
2. Alarm reset
3. Front panel
4. Power cord
5. Fuses
6. Mounting brackets
7. Air filter

-Front panel

1. Pressure adjust knob
2. Normal pressure indicator
3. Low pressure indicator
4. On/Off button
5. Static function

INSTALLATION


Unpack the box and then inspect for any damage which may have occurred during shipment. If there are any damages, please contact your dealer immediately.

1. Place the mattress on a common mattress, with the air hose connectors oriented to the footboard

WARNING!

- Lean the Dyna Best 5000 mattress on a 5cm mattress
- Lean the Dyna Flo 8000 mattress on the bed

2. Hang the pump onto the bed rail (footboard side), or place it on a flat surface

3. Connect air hose connectors from air mattress to the pump unit. When a “click” sound is felt or heard, the connection is completed and secured

WARNING!

- Check and ensure the air hoses are not kinked or tucked under mattress

4. Plug the power cord into electrical outlet

WARNING!

- Make sure the pump unit is suitable for the local power voltage
- Do not position the equipment so that it is difficult to disconnect the device

5. Turn the main power to ON position

BEFORE USE

Check the correct assembly, paying particular attention to air hoses connector

USING WARNINGS**DANGER -To reduce risk of electrocution-**

1. Always unplug this product immediately after using
2. Do not use while bathing
3. Do not place or store this product where it can fall or be pulled into a tub or sink
4. Do not place in or drop into water or other liquid
5. Do not reach for a product that has fallen into water. Unplug immediately

WARNING -To reduce risk of burns, electrocution, fine or injury to people-

1. Evaluate patients for entrapment risk according to facility protocol and monitor patient appropriately
2. This system is not intended to be used with patients who have spinal cord injury
3. Close supervision is necessary when this product is used on or near children. Electrical burns or choking accident may occur if a child swallow a small part detached from the device
4. Use this product only for its intended use as described in this manual. Do not use attachment not recommended by the manufacturer
5. Never operate this product if it has a damaged cord or plug, if it is not working properly, or if it has been dropped, damaged, or dropped into water. Return the product to a service centre for examination and repair
6. Keep the cord away from heated surfaces
7. Never block the air openings of this product or place it on a soft surface, such as a bed or couch, where their openings may be blocked. Keep the air openings free of lint, hairs, and other similar particles
8. Never drop or insert any object into any opening or hose
9. Do not modify this equipment without authorization of the manufacturer
10. Mattress cover has passed skin sensitization and skin irritation test. However, if you suspect that you may have had or are having an allergic reaction, please consult a physician immediately
11. If there is a possibility of electro-magnetic interference with mobile phones, please increase the distance to 3,3m between devices or turn off the mobile phone

HOW TO USE**-Turn on the pump**

1. Turn the main power switch to ON position
2. The low pressure indicator will light up when the mattress is not yet full inflated
3. For Dyna Best series: the low pressure indicator (yellow LED) will light up when the mattress is not yet fully inflated

For Dyna Flo series: the low pressure (yellow LED) will flash twice during system check up. There is no light during the setup phase until the pressure has reached the factory preset value

4. For Dyna Best series: when the appropriate pressure is reached the low pressure indicator (yellow LED) will go off and normal pressure indicator (green LED) will light up automatically
- For Dyna Flo series: once the pressure is reached, the normal pressure indicator (green LED) will light automatically

NOTE When the mattress pressure is verified, in case it is low the yellow LED will light up automatically

WARNING!

- Every time the mattress is set up for use, it is recommended that the pressure first to be set to the max
- Ensure the mattress is duly inflated before laying the patient down
Always keep the pump turned on during all the period of confinement

NOTE- ensure the mattress is duly inflated before laying the patient down.
Always keep the pump turned on during all the period of confinement

-CPR Valve

When there is an emergency to perform CPR on the patient, quickly pull the CPR valves to release air from mattress. The CPR valves are located at the head-end, right-hand-side of the mattress. The quick connector on the pump unit can be disconnected for even faster deflation.

-Adjust the pressure of the mattress

Users can adjust pressure of air mattress to a desired softness by adjusting the pressure knob. Please consult your health care professional for a suitable setting. Patient's can lie on the mattress after Normal Pressure indicator lights up, this indicates the pressure is reached to the desired softness.

WARNING!

- Check to see if the suitable pressure is selected by sliding one hand between the deflated air cells and the patient to feel's patient's buttock. Users should be able to feel the minimum contact

-Low Pressure and Alarm Function

When an abnormal low pressure occurs, the Low Pressure indicator (yellow LED) will light up. Check the CPR valve is correctly close

WARNING!

- If the pressure level is consistently low, check for any leakage (tubes or connecting hoses). If necessary, replace any damaged tubes or hoses. Or contact local qualified dealer for repair. If the pump unit is equipped without alarm buzzer, the Low Pressure indicator will light up until the low-pressure fault condition is resolved. If the pump unit is equipped with alarm buzzer, the alarm will buzz and the Low Pressure indicator light up when air pressure is below normal. To mute the alarm, simply press the Alarm Reset button on the side. However, the low pressure indicator remains until the low-pressure fault condition is solved. To re-activate the audible alarm function, please press Alarm Reset button again on the side to turn on the function

-Static Mode

If the pump unit is equipped with STATIC Low Air Loss only function, press the orange button to suspend the alternating mode of mattress. By pressing this button again, the mattress will be back to alternating therapy mode.

-Fuse Replacement

1. Disconnect the plug from mains power when a blown fuse is suspected
2. Remove the cover of the fuse holder by means of a small screwdriver
3. Insert a new fuse of the correct rating in, and replace the cover of the fuse holder back
The fuse should be rated as T1AL/250V type

STORAGE

- Fold the mattress accurately to prevent creases
- Roll from the head end towards the foot end with the valve open

MAINTENANCE

The LEVITAS devices by MORETTI are checked carefully and supplied with EC mark, once launched on the market. For patient and physician's safety we recommend to let check the device by the manufacturer or by laboratory authorized every year. In case of repair, use only original spare parts and accessories.

CLEANING AND DISINFECTION**WARNING!**

- Unplug the pump before starting the cleaning procedure

Before using the device with people follow this procedure:

- Clean the mattress with mild soap
- Do not immerse or soak pump unit
- After cleaning, dry the mattress without direct sunlight exposure.
- Check the connection tubes are not obstructed by dust or else. In this case use a soft swab to clean the tubes

NOTE Never use acids, alkalis or solvents such as acetone

CONDITIONS OF DISPOSAL

-General conditions of disposal

In case of disposal don't use the inserting container for municipal waste. We recommend to dispose the kit in the appropriate disposal areas for recycling.

TREATMENT OF BATTERIES (Directive 2002/96/EC)

At the end of its life, the product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment or at retailers that provide this service. By ensuring these batteries are disposed of correctly, you will help prevent potentially negative consequences for the environment and human health which could otherwise be caused by inappropriate waste handling of the battery. To remark the need to dispose of electrical equipment separately, the products is marked with crossed mobile waste bin.

SPARE PARTS/ACCESSORIES

For spare parts and accessories refer to main catalogue

TROUBLESHOOTING

If your questions can't be answered with above information, please contact your local agent directly. They might require a technician to take care of the problem

PROBLEM	SOLUTION
Power is not ON	<ul style="list-style-type: none"> • Check if the plug is connected to mains • Check for a blown fuse
Alarm is on (Audible & visual)	<ul style="list-style-type: none"> • Check if the power is suddenly shut down • Check if the CPR is sealed • Check if the connection between air tube connector to pump unit is tightly secured • Check if all tubing connections along mattress are secured
Patient is bottoming out	<ul style="list-style-type: none"> • Pressure setting might be inadequate for the patient, adjust comfort range 1 to 2 levels higher and wait for a few minutes for best comfort
Mattress form is loose	<ul style="list-style-type: none"> • Check if all the snap buttons or straps of mattress are all securely fastened • Check if the mattress is fixed to the bed frame by straps
The air hose connection doesn't connect to the pump	<ul style="list-style-type: none"> • Gently press the connector before connect to the pump • Ensure the connector is duly connected to the pump (a "click" sound should be heard)

ELECTROMAGNETIC COMPATIBILITY DECLARATION

Guidance and manufacturer's declaration- Electromagnetic emissions-

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Emission test	Compliance	Electromagnetic environment- guidance
Rf emissions CISPR 11	Group 1	The device uses RF Energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR11	Class B	The device is suitable for use in all establishment. Including domestic establishments and those directly connected to the public low-voltage power supply network.
Harmonic emissions IEC61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC61000-3-3	Complies	

Guidance and manufacturer's declaration- Electromagnetic immunity

The Supra devices are intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Immunity Test	IEC60601 test level	Compliance	Electromagnetic Environment-Guidance
Electrostatic Discharge (ESD) IEC61000-4-2	± 6kV contact ± 8kV air	± 6kV contact ± 8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/ burst IEC61000-4-4	± 2kV for power supply line ± 1kV for input/out line	± 2kV for power supply line ± 1kV for input/out line	Mains power quality should be that of atypical commercial or hospital environment
Surge IEC61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s)	Mains power quality should be that of atypical commercial or hospital environment.

<p>Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11</p>	<p><5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 sec</p>	<p><5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 sec</p>	<p>Mains power quality should be that of atypical commercial or hospital environment. If the user of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.</p>
<p>Power frequency (50/60Hz) magnetic field IEC61000-4-8</p>	<p>3 A/m</p>	<p>3 A/m</p>	<p>Power frequency magnetic fields should be at levels characteristic of atypical location in a typical commercial or hospital environment.</p>
<p>NOTE UT is the a.c. mains voltage prior to the application of the test level</p>			

Guidance and Manufacturer's Declaration - Electromagnetic Immunity:

This Supra devices are intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Immunity Test	IEC60601 test level	Compliance	Electromagnetic Environment-Guidance
<p>Conducted RF IEC 61000-4-6</p>	<p>3Vrms150 kHz to 80 MHz outside ISM bandsa</p>	<p>3 Vrms</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$ 150 kHz a 80 MHz $d = 1.2\sqrt{P}$ 150 kHz a 80 MHz $d = 2.3\sqrt{P}$ 80 MHz a 2.5G MHzWhere P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).b Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey c, should be less than the compliance level in each frequency ranged. Interference may occur in the vicinity of equipment marked with the following symbol: </p>
<p>Radiated RF IEC 61000-4-3</p>	<p>3 V/m 80 MHz to 2.5 GHz</p>	<p>3 V/m</p>	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

NOTE 3 Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and landmobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device

NOTE 4 Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and this device:

This device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as recommended below, according to the maximum output power of the communications equipment. For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz a 80 MHz $d = 1.2\sqrt{P}$	150 kHz a 80 MHz $d = 1.2\sqrt{P}$	80 MHz a 2.5G Hz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

TECHNICAL FEATURES
-Pump unit

FEATURES	LTM610	LTM615
Power supply	AC 220-240V~ 50 HZ 0,17A	AC 220-240V~50 HZ 0,17A
Power	38 W	38 W
Fuse rating	T1AL 250V	T1AL 250V
Pressure	30-60 mmHg	10-40 mmHg
Dimensions	30 x 10,5 x 22,5 cm	30 X 10,5 X 22,5 CM
Cycle	8 minuti	10 minuti
Weight	3,4 Kg	3,4 Kg
Environment temperature	Operation: +10°C / +40°C Storage: -15°C/+50°C Shipping: -15°C/+70°C	Operation: +10°C / +40°C Storage: -15°C/+50°C Shipping: -15°C/+70°C
Environment humidity	Operation: 10%-90% non- condensing Storage: 10%-90% non- condensing Shipping: 10%-90% non- condensing	Operation: 10%-90% non- condensing Storage: 10%-90% non- condensing Shipping: 10%-90% non- condensing
Atmospheric pressure	Operation: 700-1013,25 hPa	Operation: 700-1013,25 hPa
Classification	Class II, Type BF, IP21 Applied part: air mattress Not suitable for use in the presence of a flammable anesthetic mixture (No AP or APG protection)	Class II, Type BF, IP21 Applied part: air mattress Not suitable for use in the presence of a flammable anesthetic mixture (No AP or APG protection)

-Mattress

FEATURES	MATTRESS FOR LTM610	MATTRESS FOR LTM615
Dimensions (without flaps)	200 x 90 x 12,7(h) cm	200 x 90 x 20,3(h) cm
Material	NYLON + PU	NYLON + PU
Elements number	17	20
Elements with micro-openings	9	9
Max weight	140 Kg	160 Kg
Weight	5,2 Kg	10 Kg
Pressure	30-60 mmHg	10-40 mmHg

WARRANTY

Moretti products are guaranteed from material or manufacturing faults for 2 years from the purchasing date, except possible exclusion or restriction as follows. The warranty shall not be applied in the possible damages caused by improper use, abuse or alteration, and the warranty will not be valid if the instructions for use are not strictly adhered to. The correct intended use is specified in this manual.

Moretti is not responsible for consequent damages, personal injuries or whatever caused by or in relation to wrong installation or improper use.

Moretti warranty does not cover damages resulting from: natural disaster, not authorized maintenance or repairs, faults caused by problems on electricity supply (when necessary), use of spare parts not covered by Moretti, improper use, not authorized alteration, shipment damages (different from original Moretti shipment), or in case of insufficient maintenance as indicated in the manual.

The warranty doesn't cover components subject to wear and tear during the correct use of the device.

REPAIRING

-Warranty repair-

If a Moretti item presents material or manufacturing faults during the warranty period, Moretti will confirm with customer if the fault can be covered from warranty. Moretti, at its unquestionable discretion, can repair or replace the item, by a Moretti dealer or to Moretti headquarters. Labor cost can be charged to Moretti if the repair is covered by warranty. A repair or a replacement doesn't extend the warranty.

-Repair not covered by warranty-

A product out of warranty can be sent after Moretti authorization. The labor and shipping costs for good out of warranty are to be paid by the customer or by the dealer. The repairs are guaranteed for 6 months from the good received.

-Non-defective device-

The customer will be informed if, after the device return and examination, Moretti declares that the device is not faulty. In this case the good will be sent back to customer, the shipping charge will be paid by the customer.

SPARE PARTS

The original Moretti spare parts are guaranteed for 6 months from the receiving date.

EXEMPT CLAUSES

Moretti does not offer any other declaration, explicit or implicit warranty or conditions, including possible declarations, warranties or conditions of merchantability, fitness for a specific purpose, non infringement and non interference, all but what expressly specified in this warranty. Moretti does not guarantee the nonstop and faultless usage. The duration of possible implicit warranties which can be imposed by the law is limited by the warranty period, in the limits of law. Some states or countries don't allow limiting the implicit warranty or the exclusion or the limitation for accidental damages. In such countries, some of those exclusions or limitation may not be applied to the user. The present warranty may be modified without prior notification.

**WARRANTY CERTIFICATE**

Product _____

Purchased on (date) _____

Retailer _____

Address _____ **Town/city** _____

Sold to _____

Address _____ **Town/city** _____

**APEX MEDICAL CORP.**

EC REP

APEX MEDICAL S.L. Maximo Aguirre 18 Bis, 8a planta 48011 Bilbao, Vizcaya. Spain*Distribuito da/Distributed by* **MORETTI S.P.A**

Via Bruxelles, 3 - Meleto 52022 Cavriglia (Arezzo) Tel. +39 055 96 21 11

www.morettispa.com email: **info@morettispa.com****MADE IN P.R.C.**

Distribuito da / Distributed by:

Via Bruxelles, 3 - Melegnano
52022 Cavriglia (Arezzo)
Tel. +39 055 96 21 11

Fax. +39 055 96 21 200
www.morettispa.com
info@morettispa.com