**KYARA**

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 Physio-Med Services 7-23 Glossop Brook Business Park, Surrey Street, Glossop,  
Derbyshire, SK13 7AJ, England.

Distribuito da / Distributed by:  
MORETTI S.P.A.

Via Bruxelles, 3 - Meleto 52022 Cavriglia (Arezzo) Tel. +39 055 96 21 11  
[www.morettispa.com](http://www.morettispa.com) email: [info@morettispa.com](mailto:info@morettispa.com)

MADE IN TAIWAN

# INDEX

1. INTRODUCTION
  - 1.1 Pain explanation
  - 1.2 TENS explanation
  - 1.3 How TENS works
  - 1.4 Common pathology treated with TENS
  - 1.5 Electrodes placement
2. CAUTIONS AND WARNINGS
3. GENERAL DESCRIPTION
4. DEVICE VIEW
5. TECHNICAL FEATURES
6. SPARE PARTS
7. ACCESSORIES
8. GRAPHIC SYMBOLS
9. PARAMETERS CONTROL
  - 9.1 Pulse duration
  - 9.2 Pulse frequency
  - 9.3 Treatment Mode
10. ATTACHMENT OF ELECTRODE LEAD WIRES
11. LEAD MAINTENANCE
12. ELECTRODE OPTIONS
13. ELECTRODE PLACEMENT
  - 1.1 Contiguous placement
14. SKIN CARE ADVICE
15. APPLICATION OF RE-USABLE SELF ADHESIVE ELECTRODES
  - 1.1 Application
  - 15.2 Removal
  - 15.3 Care and storage
  - 1.4 Important warnings
16. ADJUSTING THE CONTROL
  - 1.4 Slide cover
  - 1.5 Display Led
  - 1.6 On/Off Switch and Intensity Control
  - 1.7 Lead connector
  - 1.8 Mode control
  - 1.9 Pulse frequency control
  - 1.10 Pulse width control
  - 1.11 TENS with Timer (LTM711) or without Timer (LTM710)

- 17. **BATTERY INFORMATIONS**
  - 1.4 Check/ Replacement Battery
  - 1.5 Precautions
  - 1.6 Rechargeable battery (do not included)
  - 1.7 How to charge battery
- 18. **MAINTENANCE, TRANSPORTATION AND STORAGE OF TENS DEVICE**
- 19. **SAFETY-TECHNICAL CONTROL**
- 20. **MALFUNCTION**
- 21. **ELECTRIC DEVICE AND BATTERY DISPOSAL**
- 22. **CONFORMITY TO SAFETY STANDARDS**
- 23. **WARRANTY**

**II Class Medical Device**  
according to ECC directive 93/42

**CODES:**

LTM710                                      two channels TENS without timer

LTM711                                      two channels TENS with timer

## **CHAPTER 1: INTRODUCTION**

Thank you for purchasing TENS KYARA by Moretti . Their design and their quality guarantee maximum comfort, liability and safety. This instruction manual contains some useful tips for a correct use of the device you have bought as well as some important advises for your safety. In case of doubts, please contact your retailer who will be able to help you correctly.

### **1.1 Pain explanation**

Pain is a warning system and the body's method of telling us that something is wrong. Pain is important; without it abnormal conditions may go undetected, causing damage or injury to vital parts of our bodies. Even though pain is a necessary warning signal of trauma or malfunction in the body, nature may have gone too far in its design. Aside from its value in diagnosis, long-lasting persistent pain serves no useful purpose. Pain does not begin until a coded message travels to the brain where it is decoded, analyzed, and then reacted to. The pain message travels from the injured area along the small nerves leading to the spinal cord. Here the message is switched to different nerves that travel up the spinal cord to the brain. The pain message is then interpreted, referred back and the pain is felt.

### **1.2 Tens explanation**

Transcutaneous Electrical Nerve Stimulation is a non-invasive, drugfree method of controlling pain. TENS uses tiny electrical impulses sent through the skin to nerves to modify your pain perception. TENS does not cure any physiological problem; it only helps control the pain. TENS does not work for everyone; however,

in most patients it is effective in reducing or eliminating the pain, allowing for a return to normal activity.

### 1.3 How Tens works

There is nothing “magic” about Transcutaneous Electrical Nerve Stimulation (TENS). TENS is intended to be used to relieve pain. The TENS unit sends comfortable impulses through the skin that stimulate the nerve (or nerves) in the treatment area. In many cases, this stimulation will greatly reduce or eliminate the pain sensation the patient feels. Pain relief varies by individual patient, mode selected for therapy, and the type of pain. In many patients, the reduction or elimination of pain lasts longer than the actual period of stimulation (sometimes as much as three to four times longer). In others, pain is only modified while stimulation actually occurs. You may discuss this with your physician or therapist.

### 1.4 Common pathology treated with Tens

The TENS has been used successfully for several pathology, in which:

#### **CHRONIC PAIN**

Cervical- amputation- phantom limb- headache- lumbar ache- lumbago- legs ache- sciatica-arthritis

#### **ACUTE PAIN**

Post surgery- muscles and articulation- tendinitis- fracture- tennis elbow

### 1.1 Electrodes placement

Different electrodes placement should be test to achieve optimal results. The effective points should be found in the area treatment or near of it. The electrodes distance must be from 2” to 20”

## CHAPTER 2: CAUTIONS AND WARNINGS

- Read operation manual before use of TENS.
- We emphasize that patient with an implanted electronic device (for example, a pacemaker) should not undergo TENS treatment without first consulting a doctor. The same applies to patients with any metallic implants.
- If TENS therapy becomes ineffective or unpleasant, stimulation should be discontinued until its use is reevaluated by the physician or therapist.
- Avoid adjusting controls while operating machinery or vehicles.
- Turn the T.E.N.S. off before applying or removing electrodes.
- Stimulate T.E.N.S. devices have no AP/APG protection.

### **Do not use it in the presence of explosive atmosphere and flammable mixture.**


- Caution should be used in applying TENS to patients suspected of having heart disease. Further clinical data is needed to show there are no adverse results.
- Electrical stimulation safety has not been established during pregnancy. Do not use TENS during pregnancy.
- Do not place electrodes on the front of the throat as spasm of the Laryngeal and Pharyngeal muscle may occur. Do not stimulate over the carotid nerve, particularly with patients with known sinus reflex sensitivity.
- Care should be taken so that when operating potentially dangerous machinery the stimulator controls are not changed abruptly.
- Cases of skin irritation at the electrode site have been reported. Stimulation should be stopped and electrodes removed until the cause of the irritation can be determined.
- Electrodes should not be placed over the eyes, in the mouth, or internally.
- Keep this device out of the reach of children.

## CHAPTER 3: GENERAL DESCRIPTION

The TENS is a battery operated pulse generator that sends electrical impulses electrodes to the body and reach the nerves causing pain. The device is provided with two controllable output channels, each independent of each other. An electrode pair can be connected to each output channel.

The electronics of theTENS create electrical impulses whose Intensity, duration, number per second and modulation may be altered with the controls or switches. Dial controls are very easy to use and the slide cover prevents accidental changes in the setting.

## CHAPTER 4: DEVICE VIEW

## CHAPTER 5: TECHNICALS FEATURES

1. Channel: two isolated between them
2. Wave form: asymmetrical bi-phases-square pulse
3. Intensity: adjustable from 14mA rms on a charge of 500ohm
4. Voltage: from 0 to 7V<sub>rm</sub> adjustable
5. Power supply: 1 battery 9V (6F 22 type)
6. Size: 95x65x23.4 cm
7. Weight: 115gr with battery
8. Pulse frequency: adjustable from 2 to 150Hz
9. Pulse width :30 - 260  $\mu$ s adjustable
10. Burst: 2 Burst/second. Frequency 100Hz
11. Modulation: pulse width vary automatically according to a regular cycle during 6.5 sec.
  - a) Pulse width decrease linearly during 0.5 seconds, from the selected rate to a smaller rate of 60%
  - b) The minor width has a duration of 2 seconds, after which the rate comes back to the original
  - c) The cycle is repeated
12. Functioning: Temperature from 0° to 40°C Relative humidity from 30° to 75° Atmospheric pressure 700hPa
13. Tolerance untill +/-20%
14. TIMER (only LTM711) 15-30 minutes or continuous. At the end the device will turn it off automatically

## CHAPTER 6: SPARE PARTS

Only use original MORETTI accessories.

Check Moretti catalogue to find spare parts and accessories.

## CHAPTER 7: ACCESSORIES

Each set TENS are completed with standard accessories and standard label as given below

DESCRIPTION	Q.TY
Adhesive Electrodes	4 pieces
Electrodes Leads	2 pieces
V Battery, type 6F22	1 piece
Instruction Manual	1 piece
Carrying Case	1 piece

### LABEL

The label attached to the back of device contains important message about this device model, serial number, supply voltage, the name of manufacturer, CE number and classification. Please do not remove.

## CHAPTER 8: GRAPHIC SYMBOL

Warning

Degree of Electrical Protection BF

Do not insert the plug into AC power supply socket.

Consult Instructions for use

Direct Current (DC power source)

Manufacturer

Serial Number

## CHAPTER 9: PARAMETER CONTROLS

### 9.1 Pulse duration

Wider pulse duration settings will deliver stronger stimulation for any given intensity setting. As mentioned in the Controls section, by using a combination of intensity and pulse duration, it is felt that various pulse widths are capable of stimulating different groups of nerve fibres.

The wider pulse duration is needed to recruit motor fibres, whereas the narrow pulse duration is used on the more sensory fibres. The choice of which pulse duration to use is partially dependent upon the Treatment Mode and Protocol selected (refer to the appropriate section).

PDF

### 9.2 Pulse frequency

The Pulse Frequency (hertz or pulses per second) chosen depends greatly upon the type of electrode placement given to the patient. When using contiguous and dermatome electrode placements (i.e. stimulating directly through the area of pain or localized enervation), a quick pulse rate (setting greater than 80Hz on the Pulse Rate Control) is desired. The patient should not perceive individual pulses but rather have the sensation of steady continuous stimulation.

When using point treatments, it has been suggested that slow pulses be utilized (less than 10Hz). With this setting the patient should be able to slightly perceive individual pulses.

When using multiple electrode placement strategies, such as combinations of point and contiguous electrode placements, the quicker pulse rates are suggested.

Despite above recommendations, these individual patients may require slight variations of the above settings, according to the nature of their condition.

### 9.3 Treatment mode

Normal or Conventional TENS offers the practitioners complete control over all the various treatment parameters of the instrument.

***Burst Mode*** is analogous to the Low Rate TENS technique except the low frequency individual pulses are replaced by individual “bursts” of 7-10 individual pulses. It is thus a combination of Conventional TENS and Low Rate TENS. In Burst Mode, the treatment frequency is fixed by the instrument and is not adjustable with the Frequency Rate control.


***Modulated Mode*** attempts to prevent nerve accommodation by continuously cycling the treatment intensity. When using Modulated Mode increase the intensity only when the unit is at the maximum intensity of the modulation cycle. If the intensity is increased during a low intensity period of the cycle, the patient may turn up the control very slowly, so that they may feel the intensity any higher.

## **DURATION**

The relief should start after the determination of intensity. In some cases to join the relief could be necessary even 30 minutes. TENS are used for long period, with a minimum of 20-30 min. and in post-surgery cases even for 36 hour. The relief could disappear in 30 min after stimulation cessation.

## **CHAPTER 10: ATTACHMENT OF ELECTRODE LEAD WIRES**

The wires provided with the system insert into the jack sockets located on top of the device. Holding the insulated portion of the connector, push the plug end of the wire into one of the jacks (see drawing); one or two sets of wires may be used. After connecting the wires to the stimulator, attach each wire to an electrode. Use care when you plug and unplug the wires. Jerking the wire instead of holding the insulated connector body may cause wire breakage.

## CAUTION!

Do not insert the plug of the patient lead wire into the AC power supply socket.

## CHAPTER 11: LEAD WIRE MAINTENANCE

Clean the wires by wiping with a damp cloth. Coating them lightly with talcum powder will reduce tangling and prolong life.

## CHAPTER 12: ELECTRODE OPTIONS

Your clinician will decide which type of electrode is best for your condition. Follow application procedures outlined in electrode packing, to maintain stimulation and prevent skin irritation. Use of CE marked Adhesive Electrodes is recommended.

## CHAPTER 13: ELECTRODE PLACEMENT

The electrodes are disposable and should be routinely replaced when they start to lose their adhesive nature. If you are unsure of your electrode adhesive properties, order replacement electrodes. Replacement electrodes should be re-ordered through or on the

advice of your physician to ensure proper quality. Follow application procedures outlined in electrode packing, to maintain optimal stimulation and to prevent skin irritation.

### 13.1 Contiguous placement

This is the most common placement technique. It involves placing the electrodes alongside the area of localized pain site, in such a way as to direct the flow of current through or around the area of pain.

In a single channel application, this would involve placing each pad on either side of the pain site if the pain is localized on a limb and deep within the tissue. Pad placement on the posterior and anterior aspects of the affected limb will allow the current to flow completely through the limb and thus through the endogenous pain site.

With a two channels application, the clinician may either direct the current flow to cross through the pain site or, in what is called the “bracket” method allowing the current flow on either side of the painful area, generally through the nerve branches that feed into the pain site.

## CHAPTER 14: SKIN CARE ADVICE

To avoid skin irritation, especially if you have sensitive skin, follow these suggestions:

1. Wash the area of skin where you will be placing the electrodes, using mild soap and water before applying electrodes, and after taking them off. Be sure to rinse soap off thoroughly and dry skin well.
2. Excess hair may be clipped with scissors; do not shave stimulation area.
3. Wipe the area with the skin preparation your clinician has recommended. Let this dry. Apply electrodes as directed.
4. Many skin problems arise from the “pulling stress” from adhesive patches that are excessively stretched across the skin during application. To prevent this, apply electrodes from centre outward; avoid stretching over the skin.

5. To minimize “pulling stress”, tape extra lengths of lead wires to the skin in a loop to prevent tugging on electrodes.
6. When removing electrodes, always remove by pulling in the direction of hair growth.
7. It may be helpful to rub skin lotion on electrode placement area when not wearing electrodes.
8. Never apply electrodes over irritated or broken skin.


## CHAPTER 15: APPLICATION OF RE-USABLE SELF ADHESIVE ELECTRODES

### 15.1 Application

1. Clean and dry the skin at the prescribed area thoroughly with soap and water prior to application of electrodes.
2. Insert the lead wire into the pin connector on the pre-wired electrodes.
3. Remove the electrodes from the protective liner and apply the electrodes firmly to the treatment site.

### 15.2 Removal

1. Lift at the edge of electrodes and peel; do not pull on the lead wires because it may damage the electrodes.
2. Place the electrodes on the liner and remove the lead wire by twisting and pulling at the same time.

### 15.3 Care and Storage

1. Between uses, store the electrodes in the resealed bag in a cool dry place.
2. It may be helpful to improve repeated application by spreading a few drops of cold water over the adhesive and turn the surface up to air dry. Over Saturation with water will reduce the adhesive properties.


### 15.4 Important Warnings

1. Do not apply to broken skin.
2. The electrodes should be discarded when they are no longer adhering.
3. The electrodes are intended for single patient use only.
4. If irritation occurs, discontinue use and consult your clinician.
5. Read the instruction for use of self-adhesive electrodes before application.

## CHAPTER 16: ADJUSTING THE CONTROLS

### 16.1 Slide Cover:

A slide-on panel cover covers the controls for Pulse Width, Pulse Rate, Mode Selector and Modulation Selector. Your medical professional may wish to set these controls for you and request that you leave the cover in place.

## 16.2. Display Led

Each of the leds illuminates whenever the electronics of the device create a current impulse. Due to the capacity of the human eye, the illumination of the lamp can only be recognized up to a frequency of approximately 30 Hz.

At higher frequencies, the lamp will appear to be constantly illuminated.

## 16.3. On/Off Switch and Intensity Control

If both controls are in the off-position (white markings on the housing), the device is switched off. By turning the controls clockwise, the appropriate channel is switched on and the impulse display led will illuminate and begin to pulse according to the frequency set. The current strength of the impulses transmitted to the electrodes increases further when the control is turned clockwise. To reduce the current strength or switch the device off, turn the controls counter clockwise to the required setting or off-position.

## 16.4. Lead Connector

Connection of the electrodes is made with two lead connector. The device must be switched off before connecting the cables. Both intensity controls must be at the Off position. Electrodes must be pressed firmly on the skin.

## 16.5. Mode Control

Expose the controls by sliding front cover down from top of unit. This switch has 3 positions: **B** for **Burst** stimulation, **N** for **Constant** stimulation, and **M** for **modulation** stimulation. Push the Mode Selector until engaged in position desired.

## 16.6. Pulse Rate Control

This dial determines how many electrical impulses are applied through the skin each second. By turning these controls, the number of current impulses per second (Hz) for both channels can be continually adjusted. Unless otherwise instructed, turn the pulse rate control to the 70-120 Hz range.

## 16.7.Pulse Width Control

This dial adjusts the length of time each electrical signal is applied through the skin, which controls the strength and sensation of the stimulation. If no instructions regarding the pulse width are given in therapy, set the control to the suggested 70-120  $\mu$ s setting.

## 16.8 TENS with TIMER (LTM711) or without TIMER (LTM710)

TENS are available in two versions: with timer (LTM711) and without timer (LTM710). Timer function permit to choose the treatment duration between 15 and 30 min. or continuous.


At the end the device will turn it off automatically.

## CHAPTER 17: BATTERY INFORMATIONS

### 17.1 Check/Replace the Battery

Over time, in order to ensure the functional safety of TENS, changing the battery is necessary.

1. Make sure that both intensity controls are switched to off position.
2. Slide the battery compartment cover and remove.
3. Remove the battery from the compartment.
4. Insert the battery into the compartment. Note the polarity indicated on the battery and in the compartment.
5. Replace the battery compartment cover and slide to close

### 17.2 Precautions

1. Remove battery if equipment is not likely to be used for some time.
2. Please recycle the used battery in accordance with domestic regulation.

3. Do not throw the used battery into fire.  
If you use rechargeable batteries, please follow the instructions.

### 17.3 Rechargeable batteries (not included)

Prior to the use of a new unit, the rechargeable battery should be charged according to the battery manufacturer's instructions. Before using the battery charger, read all instructions and cautionary markings on the battery and in this instruction manual. After being stored for 60 days or more, the batteries may lose their charge. After long periods of storage, batteries should be charged prior to use.

### 17.4 Battery charging

- (1) Plug the charger into any working 110 or 220/240v mains electrical outlet. The use of any attachment not supplied with the charger may result in the risk of fire, electric shock, or injury to persons.
- (2) Follow the battery manufacturer's instructions for charging time.
- (3) After the battery manufacturer's recommended charging time has been completed, unplug the charger and remove the battery.
- (4) Batteries should always be stored in a fully charged state.

**To ensure optimum battery performance, follow these guidelines:**

- (a) Although overcharging the batteries for up to 24 hours will not damage them, repeated overcharging may decrease useful battery life.
- (b) Always store batteries in their charged condition. After a battery has been discharged, recharge it as soon as possible. If the battery is stored more than 60 days, it may need to be recharged.
- (c) Do not short the terminals of the battery. This will cause the battery to get hot and can cause permanent damage. Avoid storing the batteries in your pocket or purse where the terminals may accidentally come into contact with coins, keys or any metal objects.

### **(d) WARNINGS:**

1. Do not attempt to charge any other types of batteries in your charger, other than the nickel-cadmium rechargeable batteries. Other types of batteries may leak or burst.
2. Do not incinerate the rechargeable battery as it may explode!

## **CHAPTER 18: MAINTENANCE, CARRING AND STORAGE OF TENS DEVICE**

1. Non-flammable cleaning solution is suitable for cleaning the device. Note: Do not smoke or work with open lights (for example, candles, etc.) when working with flammable liquids.
2. Stains and spots can be removed with a cleaning agent.
3. Do not submerge the device in liquids or expose it to large amounts of water.
4. Return the device to the carrying box with sponge foam to ensure that the unit is well-protected before transportation.
5. If the device is not to be used for a long period of time, remove the batteries from the battery compartment (acid may leak from used batteries and damage the device). Put the device and accessories in carrying box and keep it in cool dry place.
6. The packed EMS device should be stored and transported under the temperature range of  $-20^{\circ}\text{C} \sim +60^{\circ}\text{C}$ , relative humidity 20% ~ 95%, Atmosphere pressure 500hPa ~ 1060hPa.

## CHAPTER 19: SAFETY-TECHNICAL CONTROLS

For safety reasons, check your TENS each week based on the following checklist.

1. Check the device for external damage.
  - deformation of the housing.
  - damaged or defective output sockets.
2. Check the device for defective operating elements.
  - legibility of inscriptions and labels.
  - make sure the inscriptions and labels are not distorted.
3. Check Led
  - led must be illuminated when switched on.
4. Check the usability of accessories.
  - patient cable undamaged.
  - electrodes undamaged.

**Please consult your distributor if there are any problems with device and accessories.**

## CHAPTER 20: MALFUNCTIONS

Should any malfunctions occur while using the Stimat TENS, check

- whether the switch/control is set to the appropriate form of therapy. Adjust the control correctly.
- whether the cable is correctly connected to the device. The cables should be inserted completely into the sockets.
- whether the impulse display led is illuminated. If necessary, insert a new battery.
- for possible damage to the cable. Change the cable if any damage is detected.

**\* If there is any other problem, please return the device to your distributor. Do not try to repair a defective device.**

## CHAPTER 21: ELECTRIC DEVICE AND BATTERY DISPOSAL

### TREATMENT OF BATTERIES (Directive 2002/96/Ec)

At the end of its life, the product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment or at retailers that provide this service. By ensuring these batteries are disposed of correctly, you will help prevent potentially negative consequences for the environment and human health which could otherwise be caused by inappropriate waste handling of the battery. To remark the need to dispose of electrical equipment separately, the products is marked with crossed mobile waste bin

### DISPOSAL OF WASTE BATTERIES (Directive 2006/66/Ec)

This symbol on the battery or on the packaging indicates that the battery provided with this product shall not be treated as household waste. By ensuring these batteries are disposed of correctly, you will help prevent potentially negative consequences for the environment and human health which could otherwise be caused by inappropriate waste handling of the battery. The recycling of the materials will help to conserve natural resources. At the end of their life hand the batteries over to the applicable collection points for the recycling of waste batteries. For more detailed information about recycling of this product or battery, please contact your local Civic Office, your household waste disposal service or the shop where you purchased the product.

## CHAPTER 22: CONFORMITY TO SAFETY STANDARDS

The Stimote N601 TENS devices are in compliance with the EN 60601-1-2:2001 and EN 60601-1:1990+A1:1993+A2:1995+A13:1996 safety standards.

## WARRANTY

The product has 2-YEAR warranty from purchase date which is certified by the receipt or invoice. The distributor reserves the right either to repair or replace necessary parts to correct defects in the materials or workmanship without any charge. The warranty does not apply to damage resulting from failure to follow the operating instructions, accidents, abuse, alteration or disassembly by unauthorized individuals. In case you need to get the product repaired and you fulfil the conditions of warranty, send it to the nearest retailer together with the purchase receipt.

## WARRANTY CERTIFICATE

Product \_\_\_\_\_

Purchasing date \_\_\_\_\_

Authorized dealer \_\_\_\_\_

Street      Place

Sold to

Street      Place

 Everyway Medical Instruments Co. Ltd 3Fl., No.5, Lane 155, Sec. 3, Bei-Sheng Rd.,  
Shen Keng Hsiang, Taipei Hsien 222 , Taiwan

 Physio-Med Services 7-23 Glossop Brook Business Park, Surrey Street, Glossop,  
Derbyshire, SK13 7AJ, England.

Distribuito da / Distributed by:  
MORETTI S.P.A.

Via Bruxelles, 3 - Meleto 52022 Cavriglia (Arezzo) Tel. +39 055 96 21 11  
[www.morettispa.com](http://www.morettispa.com) email: [info@morettispa.com](mailto:info@morettispa.com)

MADE IN TAIWAN