# GIMA

PROFESSIONAL MEDICAL PRODUCTS

Gima S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) Italy  
Italia: tel. 199 400 401 - fax 199 400 403  
Export: tel. +39 02 953854209/221/225 - fax +39 02 95380056  
gima@gimaitaly.com - export@gimaitaly.com  
[www.gimaitaly.com](http://www.gimaitaly.com)

## PODOSCOPIO - PODOSCOPE - PODOSCOPE PODOSCOPIO - PODOSKOP

Manuale d'uso - User Manual - Notice d'utilisation  
Manual del usuario - Bedienungs-und Wartungsanleitung

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## **INDEX**

### **Chapter 1 GENERAL DESCRIPTION**

1.1	General Identification data.....	page 17
1.2	General characteristics.....	page 18
1.3	Technical features and voltage .....	page 18
1.4	Classification and reference norms .....	page 19
1.5	Electrical features and voltage.....	page 20
1.6	Components identification.....	page 20
	Electrical schema.....	page 21
1.7	General safety warning .....	page 21

### **Chapter 2 INSTALLATION**

2.1	Transport and stocking .....	page 22
2.2	Conditions for installation .....	page 22
2.3	Electric connection .....	page 25

### **Chapter 3 ISTRUCTIONS FOR USE**

3.1	Use of the podoscope .....	page 25
-----	----------------------------	---------

### **Chapter 4 MAINTENANCE**

4.1.	Ordinary maintenance.....	page 26
------	---------------------------	---------

# Chapter 1

## GENERAL DESCRIPTION

### 1.1 General identification data

This instruction handbook is considered as an essential document of podoscope 27354. Therefore they must not be separated in case of cessation to third parties, in accordance with EEC Directive 93/42. This handbook gives the operator all the relevant information to carry out the following operations correctly: installation, usage and maintenance. Instructions for a correct installation are contained in the par. 2.2 Conditions for the installation.

Follow the instructions for use

Read carefully this manual for the installation, the usage and the maintenance of the podoscope.

**Commercial trademark:** Gima

**Typology:** Podoscope

**Model:** 27354

**Manufacturer:** Gima s.p.a. - Via Marconi, 1 - 20060 Gessate (MI) Italy

### Identification and labels, symbols explication

Danger of breaking with loads over 135 kg

 Via Marconi, 1 - 20060 Gessate (MI) Italy		MADE IN ITALY		
PODOSCOPIO PODOSCOPE <b>REF</b> 27354	 	SERIAL NUMBER    2010 <input type="checkbox"/> 2011 <input type="checkbox"/> <div style="background-color: #cccccc; width: 150px; height: 20px; margin: 5px 0;"></div> 2012 <input type="checkbox"/> 2013 <input type="checkbox"/>		
MAINS VOLTAGE 220 V	FREQUENCY 50/60 HZ	CURRENT 0,2 A	 2 x T 0,5A - 5x20mm	

On each device you will find a similar plate with the manufacturer's identification.

Device with applied part Type B

Follow the instructions for use

**Attention:** danger

## 1.2 Description and general characteristics

The particular system of light diffusion of the Podoscope 27354 enables you to study the patient's foot plant and to mark the parts with the heaviest pressure.


Its structure is very simple: the methacrylate structure is an excellent light conductor and it is also endowed with an inner light source: in this way you can discover variations of colour intensity on the foot plant, corresponding to higher or lower plant pressure.

The base is of ABS.

A separable net cable with plug, a jack, inner wiring conductors and a fluorescent lamp on its feeder (reactor and socket) are its electric components.

## 1.3 Technical features and voltage

Material	Transparent methacrylate
Tension	220 V single-phase
Frequency	50/60 Hz
Current	0,162 A
Fluorescent lamp (Dulux G23 Osram)	9 W
Nr. of included plugs	n. 1
	2 Fuses:
	5x20 mm: 2 x T 0,5 A - 250 V
Conditions of use	temperature: from - 10°C up to 40°C
	umidity: from da 20% up to 85%
Stocking conditions:	temperature: from 0°C up to 40°C
	umidity: from 10% up to 90%
Height	22 cm (8.67 in)
Total weight (with package)	Kg 9,8 (21,61 lb)
Maximum charge	Kg 135 (297 lb)

## 1.4 Classification and reference norms

Reference harmonized norms:

- CEI EN 60601-1 "Safety of electromedical devices"
- CEI EN 60601-1-2 "Collateral harmonized norm for electromedical devices"
- Electromagnetic compatibility."

Classification according to safety norm EN 60601-1:

- portable device
- Class I
- type B
- ordinary protection against humidity
- do never use with inflammable anaesthetics
- voltage: 220 V +/- 10%, 50/60 Hz, with safety ground-wire.

Moreover, the particular features of this device grant adequate protection and isolation against direct or indirect electric shocks and against electric or mechanic surcharge.

The thermic power implied in this podoscope avoid dangerous overheating for materials and components. Electric isolation and well-balanced distances between the surfaces avoid the formation of electric arches.

In addition to this, the podoscope has a solid mechanical structure and an excellent stability on the ground. It has no sharp edges, therefore there exists no danger of mechanical nature both for operators and patients. Each part of the inner circuits are assembled and fixed in such a solid way so as to stand up to the often adverse transport conditions.

As to indirect shocks, this reliable device has been supplied with a safety wire grounding: therefore it can be classified in class 1. As to current leakage it can be classified in class B.

No essential performance is defined for this device.

## 1.5 Electrical features and voltage

The VOSSLOH SCHWABE feeder for fluorescent or discharge lamp functions with single-phase current (220V 50/60 Hz); the manufacturer VOSSLOH SCHWABE GmbH grants the conformity of its product both as to safety prerequisites (reference norms EN 60920 and EN 60922) and to electromagnetic compatibility (reference norms EN 61000-3-2 and EN 55015).

The external cable – H05VV-F – 3 mt (9.84 ft) long, (3x1x3m (9,84 ft)) bears the trademark IMQ, has a three-pole plug and it is connected to the podoscope with a double-fuse “basin” jack.

The internal connectors of the external cable have been harmonised in conformity with norm CEI 20-20 (type H05VV-F, 1 square mm section).


## 1.6 Components identification

Perspective general view

- 1 - Transparent methacrylate surface
- 2 - Basement in ABS
- 3 - Mirror
- 4 - Jack
- 5 - Lamp switch
- 6 - Bayonet fuse holder
- 7 - Identification plate

## Electrical schema

### 1.7 General safety warning

As far as safety is concerned, this podoscope has undergone all the necessary checks and inspections as foreseen.

In order to respect the safety prerequisites during the whole life of your device, we suggest that you always take the following precautions:

- A) Always unplug the podoscope during maintenance operations
- B) Do not move the podoscope while using it
- C) Each time before the usage, check the functioning of the cable and of its electric parts.
- D) Operators must not intervene or take any initiatives on the podoscope which are not foreseen in this handbook.

**Attention:** do not charge the podoscope with loads over 135 kg (297 lb) in order to avoid any risk of breaking

**Attention:** in order to avoid any risk of electric shock, this device has to be connected exclusively to a power supply network with an efficient grounding system.

**Attention:** not complying with warnings will relieve GIMA company from any kind of responsibility.

## Chapter 2 INSTALLATION

### 2.1 Trasport and stocking

Before shipping the podoscope 27354 is covered by a film and packed in a cardboard with protective panels inside.

- The box contains the assembled podoscope and its cable (which must be obviously connected to the apparatus before using it); the caoutchou is not included
- In case of a longer storage, please store the whole package clean and dry.

### 2.2 Conditions for installation


- Unpack the podoscope completely
- Disposal of packing elements (plastic foams, polystyrene, polyethylene, cardboard) must follow your local waste disposal guidelines
- Verify that the podoscope has not been damaged during transport. In case, please contact the manufacturer immediately

**IMPORTANT** Please take notice of the fact that light stripes and/or microscopic luminescent dots on the transparent surface of the podoscope are due to the particular nature of its material and do not impair the equipment functionality.

### Right positioning and safety installation of the podoscope

Position the podoscope so as to leave enough room around it to be able to safely operate and to connect and disconnect easily the device to the power supply network.

We suggest to attach a handle on the wall as indicated in the scheme below to make easier for patients to go on and down the podoscope.

**Electromagnetic compatibility (EMC)**

The podoscope needs special precautions concerning electromagnetic compatibility (EMC) and must be installed and operated according to the EMC information contained in this manual.

In order to avoid possible risks of electromagnetic interferences, do not use radio frequency portable or mobile devices near the podoscope. Generally the podoscope should not be used near or together other equipments; in the case this would not be possible, it's necessary to observe its operativity to verify the normal functioning.

GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC EMISSIONS		
The equipment <b>27354</b> is intended for use in the electromagnetic environment specified below. The customer or the user of the equipment <b>27354</b> should assure that it is used in such an environment		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group1	The 27354 is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Compliant	
Voltage fluctuations/ flicker emissions  IEC 61000-3-3	Compliant	

**GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC IMMUNITY**

The equipment **27354** is intended for use in the electromagnetic environment specified below.  
The customer or the user of the equipment **27354** should assure that it is used in such an environment

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±2 KV contact ±4 KV contact ±6 KV contact ±2 KV air ±4 KV air ±8 KV air	±2 KV contact ±4 KV contact ±6 KV contact ±2 KV air ±4 KV air ±8 KV air	Floor should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at less 30 %
Electrical fast transient/burst IEC 61000-4-5	±2 KV for power supply lines	±2 KV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment
Surge IEC 61000-4-5	±1 KV differential mode ±2 KV common mode	±1 KV differential mode ±2 KV common mode	Mains power quality should be that of a typical commercial or hospital environment
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% $U_T$ (>95% dip in $U_T$ ) for 0,5 cycle  40% $U_T$ (60% dip in $U_T$ ) for cycle  70% $U_T$ (30% dip in $U_T$ ) for cycle  <5% $U_T$ (>95% dip in $U_T$ ) for 5 s	<5% $U_T$ (>95% dip in $U_T$ ) for 0,5 cycle  40% $U_T$ (60% dip in $U_T$ ) for cycle  70% $U_T$ (30% dip in $U_T$ ) for cycle  <5% $U_T$ (>95% dip in $U_T$ ) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of 27354 requires continued operation during power mains interruptions, it is recommended that the 27354 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 HZ) magnetic field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic field should be at levels characteristic of a typical location in a typical commercial or hospital environment.

## 2.3 Electric connection

Connect your podoscope verifying that the indications on the plate correspond to the local electric voltage, then plug it into the socket which must be supplied with a surcharge safety device.

The use of accessories, transducers and cables, except the ones sold by the producer as spare parts, can increase the emissions and decrease the EM (electromagnetic) device immunity.

# Chapter 3

## INSTRUCTIONS FOR USE

### 3.1 Use of the podoscope

Plug the podoscope in to light it up, then switch on P to position I.

Let the patient step on the transparent basis of the podoscope and make sure that the feet are right in the centre of it. Then proceed with a sight examination on the mirror.

If have not been supplied with caoutchouc the examination can be performed only if the patient is barefoot, otherwise lay the caoutchouc on the transparent surface. This material is also useful to check the validity of arch support treatments.

At the end of the examination, press P again to put it back to the 0 position and turn it off.

**IMPORTANT:** please help elderly patients step on the podoscope and remain steady

## Chapter 4

# MAINTENANCE

### 4.1 Ordinary maintenance

- Clean the podoscope at least once a week, always after having unplugged it. Use clean clothes, slightly moistened with water and detergent for glasses. **WARNING:** never use alcohol or solvents for this purpose.
- The frequency of the checks on internal wiring may vary according to external conditions (i.e. presence of dust) and repeated usage.
- As a general rule, we recommend a sight inspection every six months or at least within one year after the purchase of the podoscope, to ascertain whether there is no wear and tear in the cable or in the electric connections.
- If you have to replace the lamp, check that the latter is of the same kind and class as the original.  
Disconnect the cable from the electric socket; put onto the work table a suitable protection (fabric or anything else) where you will lie the methacrylate side of the podoscope. Now the device foot will be upwards. With a cross-head screwdriver unscrew the four screws at the top of the rectangular plastic, covering the fluorescent lamp housing.  
Raise the rectangular plastic supporting the reactor connected to the fluorescent lamp, and put it down by turning it towards the podoscope inner basis. Before going on, observe carefully how the lamp which has to be replaced, is assembled. With small nippers cut the fastening strip holding the lamp to the power supply reactor's connector. Take out the lamp from the power supply reactor's connector. Prepare the new lamp by threading the new fastening strip through it. Take the new lamp near the reactor' socket and let the thinnest strip cable thread through the reactor's connector.  
Press the lamp into the reactor' socket; close by locking the new holding strip which will hold the lamp to the reactor's connector; than snip off the excess strip. Place again the assembled components into the podoscope lodgement and tighten down the 4 screws at the top of the rectangle. Now the podoscope is ready to work again.
- If you have to replace one of the fuses, open the fuse holder placed in the position indicated in the picture on page 21. Then get out the blown fuse with a suitable instrument and replace it with a T 0,5 A – 5x20 250V fuse.
- In the case of troubles or needs to replace spare parts, always contact the producer.

**Warning:** all maintenance should be done only after disconnecting the device from the electrical mains.

**Disposal:** *The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment. For further information on recycling points contact the local authorities, the local recycling center or the shop where the product was purchased. If the equipment is not disposed of correctly, fines or penalties may be applied in accordance with the national legislation and regulations.*

## **GIMA WARRANTY CONDITIONS**

Congratulations for purchasing a GIMA product. This product meets high qualitative standards both as regards the material and the production. The warranty is valid for 12 months from the date of supply of GIMA. During the period of validity of the warranty, GIMA will repair and/or replace free of charge all the defected parts due to production reasons. Labor costs and personnel traveling expenses and packaging not included. All components subject to wear are not included in the warranty. The repair or replacement performed during the warranty period shall not extend the warranty. The warranty is void in the following cases: repairs performed by unauthorized personnel or with non-original spare parts, defects caused by negligence or incorrect use. GIMA cannot be held responsible for malfunctioning on electronic devices or software due to outside agents such as: voltage changes, electro-magnetic elds, radio interferences, etc. The warranty is void if the above regulations are not observed and if the serial code (if available) has been removed, cancelled or changed. The defected products must be returned only to the dealer the product was purchased from. Products sent to GIMA will be rejected.

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**WARNING**  
**NEVER CLEAN THE PODOSCOPE**  
**WITH ALCOHOL**

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*Results to comply with the following EU directive:*

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**93/42/EEC DIRECTIVE "Medical Devices" of 14 June 1993 and following changes and integrations  
(2007/47/EC DIRECTIVE)**

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**EN 60601-1-2 Class. CEI 62-50 (3rd edition - 2007) "Electromedical devices part 1: General prescriptions for fundamental safety and essential performances – Collateral norm: electromagnetic compatibility".**

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**GIMA S.p.A.**