

**SUCTION FILTER CAT 28229**

**BACTERIAL/VIRAL PERFORMANCE**

**CASING: TRANSLUCENT STYROLUX**

**CONNECTIONS: 8mm HOSE BARBS**

**HYDROPHOBIC PTFE MEDIA DISC**

**99.9999% BACTERIAL AND VIRAL EFFICIENCY  
ON NELSON LABORATORIES USA**

**EFFECTIVE FILTRATION AREA = 24.6 cm<sup>2</sup>**

**OTHER CONNECTORS AVAILABLE**

**MAXIMUM OPERATING TEMPERATURE = 121 DEGREES C**

**MAXIMUM OPERATING PRESSURE = 20 psi**

**FILTERS DOWN TO .027 MICRON**

**FOR STERILE APPLICATIONS ETO  
(ETHYLENE OXIDE) IS STANDARD**



**GIMA Spa - Via Monza, 102 - 20060 Gessate (MI) - Italia**

**ITALIA: Tel. 199 400 401 (8 linee r.a.) - Fax 199 400 403**

**E-mail: gima@gimaitaly.com - www.gimaitaly.com**

**INTERNATIONAL: Tel. ++39 02 953854209 - Fax ++39 02 95380056**

**E-mail: export@gimaitaly.com - www.gimaitaly.com**

