



Numero Verde
800 172547

Da cellulare
030 7777942

SENSORE DA DITO - NELLCOR - 0.9 m - tipo Y - neonatale

Codice di riferimento: 35255

Marca: GIMA

Prezzo: 91.5€

CARATTERISTICHE:

OEM - Marchio: Nellcor N: 6 Tipo di sensore: Tipo Y neonatale Lunghezza cavo: 0.9 m Cavo d'estensione: (cod. 35304)
*cavo d'estensione disponibile. Tutti i prodotti sono forniti con manuale in GB, FR, IT, ES.

CARATTERISTICHE TECNICHE:

Compatibilità Artema: Diascope Bionet: BMS Bionics: BPM-200 Datascope: (Mindray): Accutorr Plus, Expert, Passport, Accutorr V, DPM Central Station, DPM1, 2, 3, 4, 5, 6, 7, VS-800, PM-50 Drager: Cicero, Narkomed, Vitalert 1000, Infinity Vista GE Medical: Corometrics 118, 120, 129, 188, 510, 556, Critikon 9710, 9720, Plus, Compact, Dinamap Select, MPS, Oxyshuttle, Pro 300/400, Eagle Series, Dash Series, Hellige Eagle, 1000, 3000, Eagle SMK, SMU, SMV, SOLAR Series, TRAM Modular x50SL Series Invivo Research: MI2 (Nellcor SpO2), Millennia 4500 Plus 2 Invivo/MDE Escort 100 OPT11A Escort 300A, Escort II OPT30, Prism (Nellcor SpO2) Ivy Biomedical: 405A, 4050, 405P, 405T, Vital-Guard 450 CN Medtronic Physio Control: Lifepak 12 Mennen: Horizon XL (Nellcor SpO2), MR 1300 Mercury 1200; Nellcor: N100 (req. preamp cable), N180, N200 (req. preamp cable), N-20, NPB-40, NPB-75, NPB-290, NPB- 295, N-390, N-395, N-3000, N6000 Nihon Kohden: BSM-1102 Life Scope EC, BSM2304 iPro (Life Scope i), BSM-4102, 4112 Procyon LT (Life Scope P), BSM-4104, 4114 Procyon (Life Scope P), BSM-9510 Life Scope M (Nellcor SpO2), OGS2000 PocketCare (Nellcor SpO2) Pace Tech: Vitalmax 800 Series (cara negra Nellcor SpO2) Philips 78352C, 78354C, 78834C, MI020A, M1025B, M1205A, M1350B, M1 722A/B, M1732A/B, M2475B, M2601A, M3000A

MARCA E MODELLO:

GIMA 35255

CONTENUTO DELLA CONFEZIONE:

1 pz.

